

thermo**scientific**



2018/2019 Microbiology Products Catalog

Clinical • Food & Beverage • Pharmaceutical • Environmental • Veterinary

ThermoFisher
SCIENTIFIC

HOW TO ORDER

Contact details

Thermo Fisher Specialty Diagnostics Ltd
Wade Road
Basingstoke
Hampshire
RG24 8PW
UK

Tel: +44 (0) 1256 841144

Fax: +44 (0) 1256 334994

Email: oxoid.orders@thermofisher.com

Web: www.thermofisher.com/microbiology

UK Office Hours:

Mon-Fri 08.00 - 17.00

Sat & Sun closed

Technical Support helpline

Our web site is intended to make it easier for you to ask questions or raise issues with us. Simply go to www.thermofisher.com/microbiology, and choose the most appropriate option in "Contact Us".

Technical enquiries can also be made by telephone, by fax or online:

Tel: +44 (0) 1256 694238

Fax: +44 (0) 1256 329728

Email: microbiology.techsupport.uk@thermofisher.com

Antimicrobial Susceptibility Testing	6	Diagnostics	61
AST MIC Strips	6	Biochemical Identification	61
AST Accessories	7	Immunological Tests	63
AST Media	7	Blood and Serum Tests	65
Antimicrobial Susceptibility Testing Discs	8	Toxin Detection Kits	66
Antifungal Susceptibility Testing Discs	10	Agglutination Tests	66
AST Microbroth Dilution Plates	11	Enzyme Immunoassays	69
AST Microbroth Dilution Instrumentation and Accessories	12	ELISA Accessories	71
AST Microbroth Dilution Media	12	Food Allergen ELISA Assays	71
Incubators	13	Immunofluorescence Assays	71
Atmosphere Generation	14	Lateral Flow Assays	72
AGS	14	Reagents and Stains	73
Blood Culture	15	Laboratory Supplies	74
Manual Blood Culture Systems	15	Molecular Products	75
Automated Blood Culture System	15	Sample Prep for Molecular and Non-Molecular Workflows	75
Culture System Veterinary	16	Food Molecular Detection and Quantification	76
Miscellaneous	17	SureTect System	78
Transportation Swabs Aluminium Caps	17	Hygiena Bax System	79
Culture Media	18	Microbial Characterization System	80
Dehydrated Culture Media	18	PathoProof System	81
Chromogenic Media	21	Quality Control Organisms	82
Pre-Weigh Dehydrated Culture Media	40	Quality Control Organism, qualitative	82
Supplementary Reagents	40	Quality Control Organism Sets, qualitative	86
Antibiotic Single Supplements	41	Quality Control Organism, quantitative	88
Blood Products	41	Quality Control Organism Sets, quantitative	88
Ready Poured Plates	42	Water Quality Control Organism and Sets	88
Disposable Analytical Filters	45	WaterBugs	88
Filter Membranes and Accessories	45		
Reusable Filter Holders and Funnels	46		
Bottled Media	51		
Bagged Media	54		
Dip Slides	56		
Biochemical Reagents	56		
Collection and Transportation Products	57		
Laboratory Preparations and Biological Extracts	57		
Veggiatones	58		
Process Simulation	59		
Dehydrated Media	59		
BioProcess Containers	59		
Environmental Monitoring	60		
Air Sampler	60		
Contact Plate Applicator	60		
Environmental Monitoring Swabs	60		
Prepared Media	60		
Dehydrated Media	60		

Power up your Salmonella Screening Program – rapidly detect three targets in one simple test

Thermo Scientific RapidFinder Salmonella Multiplex PCR Assay

Make faster decisions with simultaneous screening for Salmonella species, S. Typhimurium and S. Enteritidis in poultry, pork and environmental samples following a straightforward workflow with few handling steps.

The ThermoScientific™ RapidFinder™ Salmonella Multiplex workflow has been independently validated and delivers results for all three targets in as few as 16 hours for most sample types.



More testing flexibility

Extended validations for Thermo Scientific SureTect PCR Assays

With extended AOAC-RI Performance Tested Methods certification and NF VALIDATION granted by AFNOR certification, Thermo Scientific™ SureTect™ PCR Pathogen Detection Assays, as well as all our other food testing PCR Kits, can now be run on the Applied Biosystems™ 7500 Fast 96-well Real-Time PCR Instrument.



Greater laboratory efficiency

Expanded range of convenience format dehydrated culture media

New format culture media offer greater efficiency and convenience in your laboratory. Thermo Scientific™ FitBag™ and QuickBag™ Bagged Enrichment Media remove the need for weighing of powdered culture media and heat sterilization prior to use, therefore simplifying the media preparation process and reducing overall media preparation time from several hours to just a few minutes.

Thermo Scientific™ Pre-Weigh Dehydrated Culture Media are available in a choice of 1L and 8L 'rip and tip' pouches that contain the same high quality media as our standard DCM range.



Accurate and compliant coliform detection in water

Thermo Scientific Chromogenic Coliform Agar (ISO)

Thermo Scientific™ Chromogenic Coliform Agar (ISO) for water testing is compliant with the ISO9308-1:2014 standard and offers easy to read results in samples with low level bacterial contamination because of the large colony size produced and vibrant coloration of the target organisms.



Expand treatment options with new antimicrobials for Gram negative infections

Thermo Scientific Sensititre Gram Negative DKMGN and EURGNCOL MIC Plates

Test with confidence and meet EUCAST recommendations for antimicrobial susceptibility testing via broth microdilution with our new Thermo Scientific™ Sensititre™ Gram Negative MIC Plates, including colistin and the latest antimicrobial combinations, ceftazidime/avibactam and ceftolozane/tazobactam, for reliable MDRO confirmation testing.



The latest antimicrobials for the most critical cases

Ceftazidime/Avibactam now available for AST with Thermo Scientific™ Oxoid™ Antimicrobial Susceptibility Disks and Thermo Scientific Sensititre Custom MIC Plates


Our new Thermo Scientific AST products give microbiology laboratories the ability to perform susceptibility testing using ceftazidime/avibactam, the latest antimicrobial combination for treating multidrug-resistant Gram negative infections. This antimicrobial combination provides a new therapy option for clinicians struggling to treat complicated bacterial infections that are increasingly resistant to antibiotics.



One plate, two critical tests

Group B Streptococci (GBS) and Gram positive organism screening on one plate

Combined in a bi-plate format, Thermo Scientific Oxoid Columbia CNA Agar and Thermo Scientific™ Brilliance™ GBS Agar enable simple screening of GBS and Gram positive organisms on one plate. The combination of media allows laboratories to get a full clinical picture for GBS and Gram positive organisms from just one sample.





Accurate MIC values with M.I.C.E. Strips

Antimicrobial gradient strips for establishing Minimum Inhibitory Concentration (MIC) values. Thermo Scientific™ Oxoid™ M.I.C.Evaluator™ Strips provide a gradient of stabilized antimicrobial covering 15 doubling dilutions that conform to international standards (CLSI and ISO).

Easy to read MIC values

- Scale conforms to international standards
- Large font simplifies reading

Simple to perform

- Easy to manipulate
- Strip handle is presented when sachet is peeled apart

Convenient

- Single strip packaging – just open and use
- 10 and 50 pack sizes available
- All antimicrobials store at 2–8°C

Excellent results

- International trials have reported equivalent performance to standard methods.



M.I.C.Evaluator Strips (M.I.C.E.™)

for the determination of the minimum inhibitory concentration (MIC) of a test organism to an antimicrobial

Amikacin (AK)		
10	256–0.015µg/mL	MA0135D
50	256–0.015µg/mL	MA0135F
Amoxicillin (AML)		
10	256–0.015µg/mL	MA0109D
50	256–0.015µg/mL	MA0109F

Amoxicillin/clavulanic acid (AMC)		
10	256–0.015µg/mL	MA0107D
50	256–0.015µg/mL	MA0107F
Ampicillin (AMP)		
10	256–0.015µg/mL	MA0110D
50	256–0.015µg/mL	MA0110F
Cefotaxime (CTX)		
10	32–0.002µg/mL	MA0111D
50	32–0.002µg/mL	MA0111F
Cefotaxime (CTX)		
10	256–0.015µg/mL	MA0112D
50	256–0.015µg/mL	MA0112F
Ceftaroline (CPT)		
10	32–0.002µg/ml	MA0127D
50	32–0.002µg/ml	MA0127F
Ceftazidime† (CAZ)		
10	256–0.015µg/mL	MA0120D
50	256–0.015µg/mL	MA0120F
Ceftriaxone (CRO)		
10	32–0.002µg/mL	MA0122D
50	32–0.002µg/mL	MA0122F
Ciprofloxacin (CIP)		
10	32–0.002µg/mL	MA0104D
50	32–0.002µg/mL	MA0104F
Clindamycin (DA)		
10	256–0.015µg/mL	MA0119D
50	256–0.015µg/mL	MA0119F
Daptomycin (DAP)		
10	256 - 0.015µg/mL	MA0125D
50	256 - 0.015µg/mL	MA0125F
Erythromycin (E)		
10	256–0.015µg/mL	MA0108D
50	256–0.015µg/mL	MA0108F
Gentamicin (CN)		
10	256–0.015µg/mL	MA0116D
50	256–0.015µg/mL	MA0116F
Gentamicin (CN)		
10	1024–0.06µg/mL	MA0117D
50	1024–0.06µg/mL	MA0117F
Imipenem (IPM)		
10	32–0.002µg/mL	MA0115D
50	32–0.002µg/mL	MA0115F
Levofloxacin (LEV)		
10	32–0.002µg/mL	MA0113D
50	32–0.002µg/mL	MA0113F
Linezolid (LZD)		
10	256–0.015µg/mL	MA0106D
50	256–0.015µg/mL	MA0106F
Meropenem (MEM)		
10	32–0.002µg/mL	MA0121D
50	32–0.002µg/mL	MA0121F
Metronidazole (MTZ)		
10	256–0.015µg/mL	MA0103D
50	256–0.015µg/mL	MA0103F
Oxacillin (OX)		
10	256–0.015µg/mL	MA0114D
50	256–0.015µg/mL	MA0114F
Penicillin G (P)		
10	32–0.002µg/mL	MA0100D
50	32–0.002µg/mL	MA0100F

Penicillin G (P)		
10	256–0.015µg/mL	MA0101D
50	256–0.015µg/mL	MA0101F
Teicoplanin (TEC)		
10	256–0.015µg/mL	MA0118D
50	256–0.015µg/mL	MA0118F
Tetracycline (TE)		
10	256–0.015µg/mL	MA0105D
50	256–0.015µg/mL	MA0105F
Tigecycline (TGC)		
10	256–0.015µg/mL	MA0124D
50	256–0.015µg/mL	MA0124F
Vancomycin (VA)		
10	256–0.015µg/mL	MA0102D
50	256–0.015µg/mL	MA0102F

Thermo Scientific™ Oxoid™ AST Media Diagnostic Sensitivity Test Agar (DSTA)

500g	for 12.5 litres medium	CM0261B
2.5kg	for 62.5 litres medium	CM0261R
5kg	for 125.0 litres medium	CM0261T

Diagnostic Sensitivity Test Agar	
10x90mm	PO0130A ⁺
Diagnostic Sensitivity Test (D.S.T.) with chocolate Horse Blood	
10x90mm	PB0133A

Haemophilus Test Medium (HTM)

HTM Base		
500g	for 11.6 litres medium	CM0898B
HTM Supplement Store at -20°C to 0°C		
10 vials	for 500mL medium	SR0158E
Haemophilus Test Medium		
10x90mm	PO0228A ⁺	
10x90mm	PO5138A ⁺	

Iso-Sensitest Agar

a semidefined medium for antimicrobial susceptibility testing

500g	for 15.9 litres medium	CM0471B
2.5kg	for 79.5 litres medium	CM0471R
5kg	for 159.0 litres medium	CM0471T

Iso-Sensitest Agar	
10x90mm	PO0144A
Iso-Sensitest Agar (Deep Fill)	
10x90mm	PO0779A
Iso-Sensitest Agar with Chocolate Horse Blood	
10x90mm	PB0147A
Iso-Sensitest Agar with Horse Blood	
10x90mm	PB0146A
Iso-Sensitest Agar with Horse Blood and 20mg/L NAD	
10x90mm	PB0378A
Iso-Sensitest Agar with Lysed Horse Blood	
10x90mm	PB0145A
Iso-Sensitest Agar with Sheep Blood	
10x90mm	PB5003A

Iso-Sensitest Broth

a semidefined medium for antimicrobial susceptibility testing

500g	for 21.4 litres medium	CM0473B
------	------------------------	---------

Mueller-Hinton Agar

conforming to international standards for antimicrobial susceptibility testing media and tested to CLSI M6-A2

500g	for 13.2 litres medium	CM0337B
2.5kg	for 66.0 litres medium	CM0337R
5kg	for 132.0 litres medium	CM0337T

Mueller-Hinton Agar	
10x90mm	PO5007A ⁺
10x90mm	PO0152A [^]
10x140mm	PO1191S ⁺
Mueller-Hinton Agar with 2% Salt	
10x90mm	PO0799A [^]
10x90mm	PO5139A ⁺
Mueller-Hinton Agar with 5% Sheep Blood	
10x90mm	PB0431A [^]
10x90mm	PB5007A ⁺
Mueller-Hinton Agar with Horse Blood	
10x90mm	PB1229A [^]
10x90mm	PB5303A ⁺

Mueller-Hinton Broth

conforming to international standards for antimicrobial susceptibility testing media

500g	for 23.8 litres medium	CM0405B
------	------------------------	---------

Sensitest Agar

500g	for 15.6 litres medium	CM0409B
2.5kg	for 78.0 litres medium	CM0409R

Wilkins-Chalgren Anaerobe Agar

recommended for antimicrobial susceptibility testing of anaerobes

500g	for 11.6 litres medium	CM0619B
------	------------------------	---------

Anaerobe Blood Agar (Wilkins Chalgren) with Neomycin	
10x90mm	PB0112A
Anaerobe Blood Agar (Wilkins Chalgren) with Nalidixic Acid and Tween™	
10x90mm	PB0113A

Wilkins-Chalgren Anaerobe Broth

recommended for antimicrobial susceptibility testing of anaerobes

500g	for 15.2 litres medium	CM0643B
------	------------------------	---------

AST Accessories

Disc Dispenser (90mm)		
each	for 6 cartridges	ST6090
Disc Dispenser (90mm)		
each	for 8 cartridges	ST8090
Disc Dispenser (100mm)		
each	for 6 cartridges	ST6100
Disc Dispenser (100mm)		
each	for 8 cartridges	ST8100
Ejector for using cartridges individually		
6		HP0053A
Desiccant Sachets		
2		HP0063A
Replacement Desiccant (Dispenser)		
each		ST9100
Cartridge Container (complete with desiccant)		
each		HP0055A
Turbidometer for use in the UK and Europe		
each		AU0107A

Antimicrobial Susceptibility Testing

Antimicrobial Susceptibility Testing Discs

Thermo Scientific™ Oxoid™ Antimicrobial Susceptibility Testing Discs

Each pack contains 5 cartridges. Each cartridge contains 50 discs. Discs in routine use should be stored at 2 – 8°C. Longer term storage should be at -20°C. Allow discs to reach room temperature before opening.

When ordering please quote the 'symbol', 'level' and order code required.

Antibiotic	Symbol	Level	Code	CLSI	EUCAST	OTHER
Amikacin	AK	30µg	CT0107B	X	X	X
Amoxicillin	AML	2µg	CT0060B			X
Amoxicillin	AML	10µg	CT0161B		X	X
Amoxicillin	AML	25µg	CT0061B			X
Amoxicillin/ clavulanic acid	AMC	3µg	CT0538B		X	X
Amoxicillin/ clavulanic acid	AMC	30µg	CT0223B	X	X	X
Ampicillin	AMP	2µg	CT0002B		X	X
Ampicillin	AMP	10µg	CT0003B	X	X	X
Ampicillin	AMP	25µg	CT0004B			X
Ampicillin/sulbactam 1:1	SAM	20µg	CT0520B	X	X	
Ampicillin/sulbactam 2:1	SAM	30µg	CT1653B			X
Apramycin	APR	15µg	CT0545B	X		
Azithromycin	AZM	15µg	CT0906B	X		X
Aztreonam	ATM	30µg	CT0264B	X	X	X
Bacitracin	B	10 units	CT0005B			
Cefaclor	CEC	30µg	CT0149B	X	X	X
Cefadroxil	CFR	30µg	CT0453B		X	
Cefamandole	MA	30µg	CT0108B	X		X
Cefepime	FEP	30µg	CT0771B	X	X	X
Cefixime	CFM	5µg	CT0653B	X	X	X
Cefoperazone	CFP	30µg	CT0193B			X
Cefoperazone	CFP	75µg	CT0249B	X		
Cefoperazone/ sulbactam 2:1	SCF	105µg	CT1727B			
Cefotaxime	CTX	5µg	CT0407B		X	X
Cefotaxime	CTX	30µg	CT0166B	X		X
Cefotetan	CTT	30µg	CT0665B	X		X
Cefovecin	CVN	30µg	CT1929B	X		
Cefoxitin	FOX	10µg	DD0026B			X
Cefoxitin	FOX	30µg	CT0119B	X	X	X
Cefpirome	CPO	30µg	CT1412B			X
Cefpodoxime	CPD	10µg	CT1612B	X	X	X
Cefprozil	CPR	30µg	CT1647B	X		

KEY

CLSI:	Clinical and Laboratory Standards Institute
EUCAST:	European Committee on Antimicrobial Susceptibility Testing
OTHER:	BSAC: British Society for Antimicrobial Chemotherapy
	SFM: Société Française de Microbiologie
	DIN: Deutsches Institut für Normung

Antibiotic	Symbol	Level	Code	CLSI	EUCAST	OTHER
Cefsulodin	CFS	30µg	CT0263B			X
Ceftaroline	CPT	5µg	CT1942B		X	
Ceftaroline	CPT	30µg	CT1941B	X		
Ceftazidime	CAZ	10µg	CT1629B		X	
Ceftazidime	CAZ	30µg	CT0412B	X		X
Ceftazidime/ Avibactam (10:4)	CZA	14µg	CT1948B		X	
Ceftazidime/ Avibactam (30:20)	CZA	50µg	CT1947B	X		
Ceftibuten	CFT	30µg	CT1662B	X	X	
Ceftiofur	EFT	30µg	CT1751B	X		
Ceftizoxime	ZOX	30µg	CT0477B	X		X
Ceftriaxone	CRO	5µg	CT1743B			X
Ceftriaxone	CRO	30µg	CT0417B	X	X	X
Cefuroxime	CXM	5µg	CT0406B			X
Cefuroxime	CXM	30µg	CT0127B	X	X	X
Cephalexin	CL	30µg	CT0007B		X	X
Cephalothin	KF	30µg	CT0010B	X		X
Cephazolin	KZ	30µg	CT0011B	X		
Cephradine	CE	30µg	CT0063B			X
Chloramphenicol	C	10µg	CT0012B			X
Chloramphenicol	C	30µg	CT0013B	X	X	
Chloramphenicol	C	50µg	CT0014B			
Ciprofloxacin	CIP	1µg	CT0623B			X
Ciprofloxacin	CIP	5µg	CT0425B	X	X	X
Ciprofloxacin	CIP	10µg	CT1615B			
Clarithromycin	CLR	2µg	CT1599B			X
Clarithromycin	CLR	5µg	CT1623B			X
Clarithromycin	CLR	15µg	CT0693B	X		
Clindamycin	DA	2µg	CT0064B	X	X	X
Clindamycin	DA	10µg	CT0015B			
Cloxacillin	OB	5µg	CT0016B			
Colistin	CT	10µg	CT0017B	X		
Colistin	CT	25µg	CT0065B			X
Colistin	CT	50µg	CT0664B			
Compound sulphonamides	S3	300µg	CT0059B			
Doripenem	DOR	10µg	CT1880B	X	X	X
Doxycycline	DO	30µg	CT0018B	X		X

X = the antibiotic concentration is listed in the applicable AST method guidance
Not all products available in all territories. Please contact your local office for details.

Antibiotic	Symbol	Level	Code	CLSI	EUCAST	OTHER
Enrofloxacin	ENR	5µg	CT0639B	X		X
Ertapenem	ETP	10µg	CT1761B	X	X	X
Erythromycin	E	5µg	CT0066B			X
Erythromycin	E	10µg	CT0019B			
Erythromycin	E	15µg	CT0020B	X	X	
Erythromycin	E	30µg	CT0021B			
Florfenicol	FFC	30µg	CT1754B	X		
Fluconazole	FCA	25µg	CT1806B	X		
Flumequine	UB	30µg	CT0666B			X
Fosfomicin	FOS	50µg	CT0183B			X
Fosfomicin/ trometamol	FOT	200µg	CT0758B			X
Framycetin	FY	100µg	CT0071B			
Fusidic acid	FD	5µg	CT0493B			
Fusidic acid	FD	10µg	CT0023B	X	X	X
Fusidic acid	FD	50µg	CT1617B			
Gentamicin	CN	10µg	CT0024B	X	X	X
Gentamicin	CN	30µg	CT0072B		X	
Gentamicin	CN	120µg	CT0794B	X		
Gentamicin	CN	200µg	CT0695B			X
Imipenem	IPM	10µg	CT0455B	X	X	X
Kanamycin	K	5µg	CT0025B			
Kanamycin	K	30µg	CT0026B	X		
Levofloxacin	LEV	1µg	CT1586B			X
Levofloxacin	LEV	5µg	CT1587B	X	X	X
Lincomycin	MY	2µg	CT0027B			
Lincomycin	MY	10µg	CT0123B			
Lincomycin	MY	15µg	CT0028B			X
Lincomycin/ neomycin	LN	75µg	CT1757B			
Lincomycin/ spectinomycin	LS	109µg	CT1758B			
Linezolid	LZD	10µg	CT1649B		X	X
Linezolid	LZD	30µg	CT1650B	X		
Lomefloxacin	LOM	10µg	CT1661B	X		
Mecillinam	MEL	10µg	CT0096B	X	X	X
Mecillinam	MEL	25µg	CT0091B			
Meropenem	MEM	10µg	CT0774B	X	X	X
Metronidazole	MTZ	5µg	CT0067B			X
Metronidazole	MTZ	50µg	CT0466B			
Mezlocillin	MEZ	30µg	CT0174B			X
Mezlocillin	MEZ	75µg	CT0192B	X		X
Minocycline	MH	30µg	CT0030B	X	X	X

Antibiotic	Symbol	Level	Code	CLSI	EUCAST	OTHER
Moxifloxacin	MXF	1µg	CT1683B			X
Moxifloxacin	MXF	5µg	CT1633B	X	X	X
Mupirocin	MUP	5µg	CT0522B			X
Mupirocin	MUP	20µg	CT1826B			X
Mupirocin	MUP	200µg	CT0523B		X	
Nafcillin	NF	1µg	CT1895B	X		
Nalidixic acid	NA	30µg	CT0031B	X	X	X
Neomycin	N	10µg	CT0032B			X
Neomycin	N	30µg	CT0033B	X		
Netilmicin	NET	10µg	CT0424B		X	X
Netilmicin	NET	30µg	CT0225B	X		
Nitrofurantoin	F	50µg	CT0069B			
Nitrofurantoin	F	100µg	CT0034B		X	
Nitrofurantoin	F	200µg	CT0035B			X
Nitrofurantoin	F	300µg	CT0036B	X		
Norfloxacin	NOR	2µg	CT0687B			X
Norfloxacin	NOR	5µg	CT0668B			
Norfloxacin	NOR	10µg	CT0434B	X	X	
Novobiocin	NV	5µg	CT0037B			
Novobiocin	NV	30µg	CT0038B			
Nystatin	NS	100 units	CT0073B			
Ofloxacin	OFX	5µg	CT0446B	X	X	X
Oleandomycin	OL	15µg	CT0039B			
Oxacillin	OX	1µg	CT0159B	X	X	X
Oxacillin	OX	5µg	CT0040B			X
Oxolinic acid	OA	2µg	CT0181B			
Oxytetracycline	OT	30µg	CT0041B			
Pefloxacin	PEF	5µg	CT0661B		X	
Penicillin G	P	1 unit	CT0152B		X	X
Penicillin G	P	1.5 units	CT0042B			
Penicillin G	P	2 units	CT0088B			
Penicillin G	P	5 units	CT0124B			
Penicillin G	P	10 units	CT0043B	X		
Penicillin V	PV	10µg	CT1900B			
Penicillin/novobiocin	PNV	40	CT1755B	X		
Pipemidic acid	PIP	20µg	CT0180B			
Piperacillin	PRL	30µg	CT1619B		X	
Piperacillin	PRL	75µg	CT0261B			X
Piperacillin	PRL	100µg	CT0199B	X		
Piperacillin/ tazobactam	TZP	36µg	CT1616B		X	
Piperacillin/ tazobactam	TZP	40µg	CT1628B			
Piperacillin/ tazobactam	TZP	85µg	CT0720B			X

Antimicrobial Susceptibility Testing

Antimicrobial and Antifungal Susceptibility Testing Discs

Antibiotic	Symbol	Level	Code	CLSI	EUCAST	OTHER
Piperacillin/tazobactam	TZP	110µg	CT0725B	X		
Pirlimycin	PIR	2µg	CT1668B	X		
Polymyxin B	PB	300 units	CT0044B	X		
Quinupristin/dalfopristin	QD	15µg	CT1644B	X	X	X
Rifampicin	RD	2µg	CT0078B			X
Rifampicin	RD	5µg	CT0207B	X	X	X
Rifampicin	RD	30µg	CT0104B			X
Spectinomycin	SH	10µg	CT0046B			
Spectinomycin	SH	25µg	CT0411B			X
Spectinomycin	SH	100µg	CT0823B	X		
Spiramycin	SP	100µg	CT0232B			X
Streptomycin	S	10µg	CT0047B	X		X
Streptomycin	S	25µg	CT0048B			
Streptomycin	S	300µg	CT1897B			X
Sulbactam/ampicillin 1:1	SAM	20µg	CT0520B	X	X	
Sulbactam/ampicillin 1:2	SAM	30µg	CT1653B			X
Sulfisoxazole	see Sulphafurazole					
Sulphafurazole (Sulfisoxazole)	SF	300µg	CT0075B	X		
Sulphamethoxazole	RL	25µg	CT0051B			
Sulphamethoxazole	RL	100µg	CT0074B			
Sulphamethoxazole/trimethoprim 19:1	SXT	25µg	CT0052B	X	X	
Sulphonamides compound	S3	300µg	CT0059B			
Teicoplanin	TEC	30µg	CT0647B	X	X	X
Telithromycin	TEL	15µg	CT1714B	X	X	X
Temocillin	TEM	30µg	CT1952B			
Tetracycline	TE	10µg	CT0053B			X
Tetracycline	TE	30µg	CT0054B	X	X	
Ticarcillin	TIC	75µg	CT0167B	X	X	X
Tigecycline	TGC	15µg	CT1841B	X	X	X
Tilmicosin	TIL	15µg	CT1756B	X		
Tobramycin	TOB	10µg	CT0056B	X	X	X
Tobramycin	TOB	30µg	CT1618B			
Trimethoprim	W	1.25µg	CT0057B			
Trimethoprim	W	2.5µg	CT0070B			X
Trimethoprim	W	5µg	CT0076B	X	X	X
Trimethoprim/sulphamethoxazole 1:19	SXT	25µg	CT0052B	X	X	
Vancomycin	VA	5µg	CT0188B		X	X
Vancomycin	VA	30µg	CT0058B	X		
Voriconazole	VOR	1µg	CT1807B	X		

Antifungal Susceptibility Testing Discs

Antibiotic	Symbol	Level	Code	CLSI	EUCAST	OTHER
Fluconazole	FCA	25µg	CT1806B	X		
Nystatin	NS	100 units	CT0073B			
Voriconazole	VOR	1µg	CT1807B	X		

Sensititre ID/AST System

Meeting the unique testing demands of all laboratory sizes and volumes, the Thermo Scientific™ Sensititre™ ID/AST System includes a full range of standard and custom MIC plates that allow you to eliminate offline testing while meeting FDA, CLSI, and EUCAST breakpoint requirements.

Combined with modular equipment options, the Sensititre System allows you to create the ultimate tests program tailored to your laboratory requirements.

More antimicrobials, more testing options for clinical and veterinary microbiology laboratories. The Sensititre System includes over 240 antimicrobials, including 40+ veterinary-specific antimicrobials, for susceptibility testing with true minimum inhibitory concentration (MIC) results, and superior quality and reproducibility.

As the system of choice for global surveillance programs including National Antimicrobial Resistance Monitoring System (NARMS) coordinated via FDA-CVM, USDA and the CDC, the Sensititre System provides standardized results for confident susceptibility testing.

Sensititre Plates

Anaerobe AST Plate
Box of 10 RUO ANO2B
Gram-Negative AST Plate
Box of 10 IVD GN2F*
Gram-Negative AST Plate
Box of 10 IVD GN3F*
Gram-Negative AST Plate
Box of 10 IVD GN4F*
Gram-Negative AST Plate with Colistin, CZA, C/T - EUCAST
Box of 10 IVD DKMGN
Gram-Negative AST Plate with CZA, C/T
Box of 10 IVD GN6F
Gram-Negative AST Plate with Colistin
Box of 10 RUO GNX2F*
Gram-Negative AST Plate with Colistin
Box of 10 RUO GNX3F
Gram-Negative Non-Fermenter AST Plate
Box of 10 IVD NF*
ESBL AST Plate
Box of 10 IVD ESB1F*
Gram-Negative Identification Plate
Box of 10 IVD GNID
Gram-Positive Identification Plate
Box of 10 IVD GPID
Gram-Positive AST Plate
Box of 10 IVD GPN3F*
Gram-Positive AST Plate
Box of 10 IVD GPALL1F*
Staph AST Plate with Telavancin - EUCAST
Box of 10 IVD EUSTAPF
Fastidious AST Plate
Box of 10 IVD HPB1*
Strep AST Plate
Box of 10 IVD STP6F*
YeastONE MIC Plate with Micafungin
Box of 10 IVD YO10
Myco TB MIC Plate
Box of 10 IVD MYCOTBI
Myco Non-TB MIC Plate - Slow-Growing
Box of 10 IVD SLOMYCOI
Myco Non-TB MIC Plate - Rapid-Growing
Box of 10 IVD RAPMYCOI
Vet Avian AST Plate
Box of 10 Vet AVIAN1F*
Vet Bovine AST Plate
Box of 10 Vet BOPO6F*
Vet Urine AST Plate
Box of 10 Vet CMV1BURF*
Vet Companion Gram-Positive AST Plate with Pradofloxacin
Box of 10 Vet COMPGP1F*

Vet Companion Gram-Negative AST Plate with Pradofloxacin
Box of 10 Vet COMPGN1F*
Vet Companion AST Plate
Box of 10 Vet COMPAN1F*
Vet Companion AST Plate
Box of 10 Vet COMPAN2F*
Vet Equine AST Plate
Box of 10 Vet EQUIN1F*
Vet Campy AST Plate
Box of 10 Vet CAMPY*
Vet Mastitis AST Plate
Box of 10 Vet CMV1AMAF*
NARMS Gram-Positive MIC Plate
Box of 10 Surveillance CMV3AGPF*
NARMS Gram-Negative MIC Plate
Box of 10 Surveillance CMV4AGNF
Vet Topicals AST Plate
Box of 10 Vet JOEYE2*
Salmonella MIC Plate
Box of 10 Surveillance EUMVS2
Gram-Positive MIC Plate
Box of 10 Surveillance EUST
Campy MIC Plate
Box of 10 Surveillance EUCAMP
Salmonella/E. coli MIC Plate
Box of 10 Surveillance EUVSEC
ESBL MIC Plate
Box of 10 Surveillance EUVSEC2
Staph MIC Plate
Box of 10 Surveillance EUVENC
Campy MIC Plate
Box of 10 Surveillance EUCAMP2



* Products are not stock items. Please contact your local representative for details of availability.

Not all products available in all territories. Please contact your local office for details.

Sensititre Equipment

Thermo Scientific™ Sensititre™ ARIS™ 2X ID/AST System

- Fully automated, bench-top incubating and reading system reduces workload and facilitates efficient workflow
- Accommodates up to 192 tests on a single instrument
- Heated carousel individually incubates plates, ensuring optimal growth conditions

Thermo Scientific™ Sensititre™ OptiRead™ Automated Fluorometric Plate Reading System

- Fast, accurate, automated plate reading for maximum consistency
- Quickly transfer results for processing, interpretation and reporting, improving efficiency and productivity
- Compact, light-weight design maximizes bench space

Thermo Scientific™ Sensititre™ Vizion™ Digital MIC Viewing System

- Combines automation with visual results, allowing you to consolidate offline testing on a single instrument
- Easy-to-read digital images reduce reading errors
- Store digital images of plate results for training purposes and complete traceability

Thermo Scientific™ Sensititre™ AIM™ Automated Inoculation Delivery System

- Quickly and accurately doses 96-well plates, eliminating skipped wells and costly repeat tests
- Icon-driven touch screen facilitates quick dosing selection for simple, intuitive operation


Nephelometer		
Each	Automated turbidity measurement	V3011
AIM		
Each	Automated plate inoculation	V3020
Manual Viewbox		
Each	Manual plate reading	V4007
Vizion System		
Each	Digital MIC Viewing System	V2021
OptiRead		
Each	Automated plate reading	V3030
ARIS 2X		
Each	Automated plate reading and incubation	V3090
SWIN Database and Computer System		
Each	Automated result interpretation	SW4000
SWIN Software Epidemiology Module		
Each		SW120

Sensititre Accessories

AIM Doseheads		
Box of 100	For use with Sensititre AIM System	E3010
Mineral Oil 250mL		
Each		MINOIL
0.5 McFarland Standard 5mL		
Each		E1041
0.9 McFarland Standard 5mL		
Each		E1039

Sensititre Media

YeastOne Broth		
Box of 10	11mL	Y3462
Mueller-Hinton Broth with TES		
Box of 100	11mL	T3462
Mueller-Hinton Broth with TES		
Box of 100	5mL	T3462-05
Mueller-Hinton Broth with TES		
Box of 10	11mL	T3462-10
Mueller-Hinton Broth with Lysed Horse Blood		
Box of 10	11mL	CP112-10
Mueller-Hinton Broth with Lysed Horse Blood for Autoread		
Box of 10	11mL	CP114-10
Demineralized Water		
Box of 10	5mL	T3339-10
Demineralized Water		
Box of 100	5mL	T3339
HTM Broth		
Box of 10	11mL	T3470
Brucella Broth for Anaerobes		
Box of 10	11mL	T3450
Sterile Water with Glass Beads		
Box of 10	5mL	T4393
Demineralized Water with Glass Beads		
Box of 10	11mL	T8006
Saline Tween with Glass Beads		
Box of 10	5mL	T3491
Middlebrook 7H9 with OADC		
Box of 10	10mL	T3441
Veterinary Fastidious Medium		
Box of 10	11mL	T3460



Thermo Scientific™ Heratherm™
Microbiological Table Top Incubators

Compact

IMC18 Heratherm Compact Incubator; Mechanical Convection;	
0.65 cu.ft. (18 L); 17 to 40 C; 120-230VAC; 50/60Hz - EU standard, UK, Australia cable/plug	HT50125882

Gravity Convection

IGS60 Heratherm General Protocol Incubator	
2.6 cu.ft. (75 L); ambient +5 C to 75 C; 230VAC; 50/60Hz	HT51028130
IGS100 Heratherm General Protocol Incubator	
4.0 cu.ft. (117 L); ambient +5 C to 75 C; 230VAC; 50/60Hz	HT51028131
IGS180 Heratherm General Protocol Incubator	
6.85 cu.ft. (194 L); ambient +5 C to 75 C; 230VAC; 50/60Hz	HT51028132

Dual Convection

IMH60 Heratherm Advanced Protocol Incubator	
2.3 cu.ft. (66 L); ambient +5 C to 105 C; 230VAC; 50/60Hz	HT51028133
IMH100 Heratherm Advanced Protocol Incubator	
3.67 cu.ft. (104 L); ambient +5 C to 105 C; 230VAC; 50/60Hz	HT51028134
IMH180 Heratherm Advanced Protocol Incubator	
6.3 cu.ft. (178 L); ambient +5 C to 105 C; 230VAC; 50/60Hz	HT51028135
IMH60-S Heratherm Advanced Protocol Security Incubator	
2.3 cu.ft. (66 L); ambient +5 C to 105 C; 230VAC; 50/60Hz	HT51028136
IMH100-S Heratherm Advanced Protocol Security Incubator	
3.67 cu.ft. (104 L); ambient +5 C to 105 C; 230VAC; 50/60Hz	HT51028137
IMH180-S Heratherm Advanced Protocol Security Incubator	
6.3 cu.ft. (178 L); ambient +5 C to 105 C; 230VAC; 50/60Hz	HT51028138
IMH60 SS Heratherm Advanced Protocol Incubator	
2.3 cu.ft. (66 L); ambient +5 C to 105 C; 230VAC; 50/60Hz; stainless steel exterior	HT51028717
IMH100 SS Heratherm Advanced Protocol Incubator	
3.67 cu.ft. (104 L); ambient +5 C to 105 C; 230VAC; 50/60Hz; stainless steel exterior	HT51028718
IMH180 SS Heratherm Advanced Protocol Incubator	
6.3 cu.ft. (178 L); ambient +5 C to 105 C; 230VAC; 50/60Hz; stainless steel exterior	HT51028719
IMH60-S SS Heratherm Advanced Protocol Security Incubator	
3 cu.ft. (66 L); ambient +5 C to 105 C; 230VAC; 50/60Hz; stainless steel exterior	HT51028541
IMH100-S SS Heratherm Advanced Protocol Security Incubator	
3.67 cu.ft. (104 L); ambient +5 C to 105 C; 230VAC; 50/60Hz; stainless steel exterior	HT51028542
IMH180-S SS Heratherm Advanced Protocol Security Incubator	
6.3 cu.ft. (178 L); ambient +5 C to 105 C; 230VAC; 50/60Hz; stainless steel exterior	HT51028543



Table Top Incubators Accessories

Shelves

Stainless steel perforated shelf for IMC18	
For IMC18 / IMC18-CHINA models	HT50125605
Stainless steel perforated shelf for IGS60, including 2 shelf supports	
For General Protocol 60 liters incubator	HT50127770
Stainless steel perforated shelf for IGS100, including 2 shelf supports	
For General Protocol 100 liters incubator	HT50127771
Stainless steel perforated shelf for IGS180, including 2 shelf supports	
For General Protocol 180 liters incubator	HT50127772
Stainless steel perforated shelf for IMH60 / IMH60-S, including 2 shelf supports	
For Advanced Protocol / Advanced Protocol Security 60 liters incubators	HT50127773
Stainless steel perforated shelf for IMH100 / IMH100-S, including 2 shelf supports	
For Advanced Protocol / Advanced Protocol Security 100 liters incubators	HT50127774
Wire mesh shelf for OMS60/ OMH60 / OMH60-S, including 2 shelf supports	
For all 60 liters incubator models	HT50127764
Wire mesh shelf for OMS100/ OMH100 / OMH100-S, including 2 shelf supports	
For all 100 liters incubator models	HT50127765
Wire mesh shelf for OMS180/ OMH180 / OMH180-S, including 2 shelf supports	
For all 180 liters incubator models	HT50127766
Shelf with holders for Petri dishes for all 60 l incubators, including 2 shelf supports	
For all 60 liters incubator models; 90mm diameter; stainless steel	HT50128816
Shelf with holders for Petri dishes for all 100 l incubators, including 2 shelf supports	
For all 100 liters incubator models; 90mm diameter; stainless steel	HT50128818
Shelf with holders for Petri dishes for all 180 l incubators, including 2 shelf supports	
For all 180 liters incubator models; 90mm diameter; stainless steel	HT50128819
Shelf with holders for Petri dishes for all 60 l incubators, including 2 shelf supports	
For all 60 liters incubator models; 50mm diameter; stainless steel	HT50128793
Shelf with holders for Petri dishes for all 100 l incubators, including 2 shelf supports	
For all 100 liters incubator models; 50mm diameter; stainless steel	HT50128794
Shelf with holders for Petri dishes for all 180 l incubators, including 2 shelf supports	
For all 180 liters incubator models; 50mm diameter; stainless steel	HT50128815
Stainless steel drip tray for all 60 l incubators; including 2 shelf supports	
For all 60 liters incubator models	HT50128683
Stainless steel drip tray for all 100 l incubators, including 2 shelf supports	
For all 100 liters incubator models	HT50128791
Stainless steel drip tray for all 180 l incubators; including 2 shelf supports	
For all 180 liters incubator models	HT50128792

Sample Sensor

Sample sensor for IMH60-S / IMH100-S / IMH180-S	
For all Heratherm Advanced Protocol Security incubators	HT50127768

Support Stands

Support stand with coasters for Heratherm 60L / 2 cu.ft. Models	
For all Heratherm 60 liter models	HT50127741
Support stand with coasters for Heratherm 100L / 3.5 cu.ft. Models	
For all Heratherm 100 liter models	HT50127742
Support stand with coasters for Heratherm 180L / 6.4 cu.ft. Models	
For all Heratherm 180 liter models	HT50127743

AGS

Safe and speedy generation of optimal growth conditions for anaerobes, microaerophiles and CO₂-dependent organisms

Oxoid AnaeroGen System

Thermo Scientific™ Oxoid™ AnaeroGen™ is an anaerobic atmosphere generation system

AnaeroGen 2.5 litre for use in 2.5 litre jar	
10 sachets	AN0025A
AnaeroGen 3.5 litre for use in 3.5 litre jar	
10 sachets	AN0035A
AnaeroGen Compact	
10 sachets and 10 pouches	AN0010C
AnaeroGen W-Zip Compact	
10 sachets and 10 W-Zip pouches	AN0010W
AnaeroGen Compact for use in plastic pouches and/or W-Zip pouches	
10 sachets	AN0020D

Oxoid CampyGen System

Thermo Scientific™ Oxoid™ CampyGen™ produces suitable gaseous atmosphere for the growth of microaerophilic organisms

CampyGen 2.5 litre for use in 2.5 litre jar	
10 sachets	CN0025A
CampyGen 3.5 litre for use in 3.5 litre jar	
10 sachets	CN0035A
CampyGen Compact for use in plastic pouches and/or W-Zip pouches	
20 sachets	CN0020C

Oxoid CO₂ Gen System

produces suitable conditions for organisms which require or grow better in a CO₂ enriched atmosphere

CO₂ Gen 2.5 litre for use in 2.5 litre jar	
10 sachets	CD0025A
CO₂ Gen Compact for use in plastic pouches and/or W-Zip pouches	
20 sachets	CD0020C

Compact Accessories

W-Zip Seal Pouches (integral seal)	
20 pouches	AG0060C
Plastic Pouches	
20 pouches	AG0020C
Sealing Clips for Plastic Pouches	
5 clips	AN0005C



Oxoid AnaeroJar 2.5l Jar

Thermo Scientific™ Oxoid™ AnaeroJar™ with a 2.5 litre capacity, and supplied complete, holds up to 12 plates

1 jar	AG0025A
-------	---------

Rectangular AnaeroBox™ Container

Thermo Scientific™ Rectangular AnaeroBox™ Container with a 2.5 litre capacity, and supplied complete, holds up to 12 90mm plates or 21 55mm. (100x130x190mm)

1 jar	AB0025A
-------	---------

AnaeroJar Accessories

AnaeroJar Base	
1	AG0026A
AnaeroJar Lid	
1	AG0027A
AnaeroJar Handle	
1	AG0028A
AnaeroJar Plate Carrier	
1	AG0029A
AnaeroJar 'O' Ring	
5	AG0030A
AnaeroJar Clips	
2	AG0031A

Anaerobic Indicator (resazurin)

100 sachets	BR0055B
-------------	---------

Oxoid Anaerobic 3.5 l Jar

Thermo Scientific™ Oxoid™ Anaerobic Jar with a 3.5 litre capacity holds up to 15 plates. Note – for connection to a vacuum pump and gas cylinder, Schrader valve chucks and clips are also required.

1 jar	HP0011A
-------	---------

Oxoid Modified Anaerobic Jar

Thermo Scientific™ Oxoid™ Anaerobic Jar with single Schrader valve and without pressure gauge.

1 jar	HP0031A
-------	---------

Anaerobic Jar Accessories

Anaerobic Catalyst	
5	BR0042A
Clamp and Screw	
1	HP0025A
Lid, complete with clamp and screw	
1	HP0014A
'O' Ring	
5	HP0018A
Plate Carrier – Stainless Steel	
1	HP0026A
Polycarbonate Jar (base only)	
1	HP0013A
Pressure Gauge	
1	HP0017A
Pressure Release Valve	
1	HP0016A
Schrader Valve Chuck and Clip	
2	HP0020A
Schrader Valve Cores	
10	HP0022A
Schrader Valve Kit	
2	HP0015A

Not all products available in all territories. Please contact your local office for details.

Manual Blood Culture Systems

Isolator™ Microbial Tubes

A high sensitive, direct-draw blood culture system using lysis-centrifugation/direct plating technology

Isolator Microbial 10 Tubes for testing adult samples	
1x25	BC0507C
Isolator Microbial 1.5 Tubes for testing paediatric samples (no centrifugation step)	
1x25	BC0505C

Accessories for Isolator Microbial 10 Tubes only

Isostat Press and Rack	
each	BC0502C
Isostat Consumable Pack contains 100 each of vent caps, supernatant and concentrate pipettes	
1x100	BC0509C
Isostat Centrifuge Containment Adaptor	
each	BC0501C

Signal Blood Culture System

The Thermo Scientific™ Signal™ One Bottle System is formulated to grow aerobic, anaerobic and microaerophilic organisms. Each pack comprises 20 Signal bottles and 20 Signal growth indicator devices.

Signal Blood Culture System	
1x20	BC0100M
Trays white disposable shaker/incubator trays to hold 10 Signal systems during shaking and incubation	
1x5	BC0104A



Automated Blood Culture System

Automated Microbial Detection System

Ultimate versatility. Superior organism recovery. Thermo Scientific™ VersaTREK™ Automated Microbial Detection System is the only system capable of detecting any gas produced or consumed by organisms. The result is faster recovery and detection of a much wider range of common and fastidious organisms. Speed and performance in the lab leads to reduced length of stay and therapy costs, and ultimately, better patient care.

Accessories for VersaTREK Instrument

VersaTREK Drawer-Left (Must order drawers 2/time; 1 left & 1 right)	
Each	6030-30
VersaTREK Drawer-Right (Must order drawers 2/time; 1 left & 1 right)	
Each	6031-30
Myco Drawer - Right (must order drawers 2/time; 1 left & 1 right)	
Each	6029-30
Myco Drawer - Left (must order drawers 2/time; 1 left & 1 right)	
Each	6028-30
Myco Adaptors, 25/box	
Box	6040-30
Disabled Location Caps, Red	
Each	6044-30
Bracket, 528 Seismic Restrainer Kit	
Each	6320-30
VersaTREK Operators Manual (CD)	
Each	6400-30
Stylus (3/pack)	
Pack	6420-30
Stylus Holder	
Each	6422-30
Myco Carrying Rack	
Each	6423-30
Bottle Carrying Rack	
Each	6424-30
REDOX Transporter (10/box)	
Box	6428-30
External Speaker w/Swivel Bracket	
	6826-30
External Speaker Adapter Cable	
	6825-30
Extension Cable, 20' Mons for Remote Alarm	
Each	6431-30
Internal Speaker	
Each	6432-30
Thermometer	
Each	6440-30
Absorbent Sheet	
Each	6444-30
VersaTREK Label, Drawer Numbers 1-22	
Each	6452-30

Blood Culture

Automated Blood Culture System

VersaTREK Starter Kit		
Each	Includes (5) 6044-30, (1) 6420-30, (1) 6424-30, (3) 6440-30, (3) 6444-30, (1) 6452-30 & (6) 6705-30	6499-31
VersaTREK Myco Starter Kit		
Each	Includes (5) 6044-30, (1) 6420-30, (1) 6452-30, (1) 6440-30, (1) 6423-30, (3) 6444-30, (1) 6172-30, (1) 6460-30 & (1) 6422-30	6497-31
Thermometer Mounts		
Each		6705-30
Conditioned UPS, 1440VA		
Each		6856-31



Computer Package for VersaTREK Instruments

VTI Computer Package Kit		
Purchase	Includes: Computer with VersaTREK software, ethernet switch and cable	6128-31

Consumables for VersaTREK Instruments

VersaTREK Myco Growth Supplement (Myco GS) 50ml (250 tests/box)		
Box		7112-42
Myco PVNA 25ml, (250tests/box)		
Box of 5		7113-42
Myco Antibiotic Supplement (Myco AS) 25ml, (250 tests/box)		
Box of 5		7114-42
Myco Susceptibility Kit; 2x25ml Rifampin; 2x25ml Isoniazid; 3x25ml Ethambutol. 80-100 tests per kit		
Kit		7115-60
Streptomycin Kit; 20 tests per kit		
Kit		7120-30
Myco PZA, 50 tests per kit		
Kit		7116-70
Connector for Aerobic, Anaerobic, Myco Bottles		
Box of 50		7150-44
Thermo Scientific™ REDOX™ 1 (Aerobic) 80 ml Bottles		
Box of 50	for 0,1 - 10 ml	7102-44
REDOX 2 (Anaerobic) 80 ml		
Box of 50	for 0,1 - 10 ml	7103-44
REDOX 1 EZ Draw™ (Aerobic) 40 ml		
Box of 50	for 0,1 - 5 ml	7106-44
REDOX 2 EZ Draw (Anaerobic) 40 ml		
Box of 50	for 0,1 - 5 ml	7107-44
REDOX 2 M EZ Draw 40 ml		
Box of 50	for 0,1 - 5 ml	7110-44

Culture System Veterinary

Accessories for Thermo Scientific™ para-JEM™ Applications

Rack, para-JEM, 24 pos.	
Each	6423-30
VersaTREK EZ-ID Labels. 100 ea.	
Each	6450-35
para-JEM Starter Kit	
Each	6498-31
Printer, Barcode Thermal Transfer	
Each	SW3000
Labels, para-JEM Barcode, 1310/RL	
Each	SW1101

Consumables for para-JEM Applications

para-JEM BLUE	
Each	00-030
para-JEM Kit (250 tests/kit)	
Kit	7130-42
para-JEM Broth	
Box of 50	7131-42
para-JEM GS	
Box of 5	7132-42
para-JEM AS	
Box of 5	7133-42
para-JEM EYS	
Box of 5	7134-42
para-JEM Kit w/BLUE (250 tests/kit)	
Kit	7137-42

Transportation Swabs, Aluminium Caps

Aluminium Test Tube Caps

Aluminium Test Tube Caps (12mm)	
1x12	TA0001A
Aluminium Test Tube Caps (16mm)	
1x12	TA0002A
Aluminium Test Tube Caps (19mm)	
1x12	TA0003A
Aluminium Test Tube Caps (25mm)	
1x12	TA0004A

Transport Swabs

Amies Agar Gel	
10x50	TS0001A
Amies Agar Gel with Charcoal	
10x50	TS0002A
Amies Liquid	
10x50	TS0003A
Amies Agar Gel Mini Tip: Firm Applicator Wire	
10x50	TS0012A

Thermo Scientific™ Oxoid™ Dehydrated Culture Media

a wide range of general purpose, enrichment, selective and differential media, diluents and supplements. NB: information provided here is intended only as an aid to purchasing. When making up media, please, refer to the relevant instructions for use. Supplementary Reagent products should be stored at 2–8°C unless otherwise stated.

Aeromonas Medium (Ryan)

a selective, differential medium for the isolation of *Aeromonas hydrophila*

Aeromonas Medium Base (Ryan)		
500g	for 8.5 litres medium	CM0833B
Ampicillin Supplement		
10 vials	for 500mL medium	SR0136E

AFPA

a selective identification medium for the detection of *Aspergillus flavus* and *A. parasiticus*

AFPA Base		
500g	for 11.0 litres medium	CM0731B
Chloramphenicol Supplement		
10 vials	for 500mL medium	SR0078E
10 vials	for 2.0 litres medium	SR0078H

Alkaline Peptone Water

a broth medium for the enrichment of *Vibrio* spp. from food, water and clinical samples

500g	for 16.7 litres medium	CM1028B
------	------------------------	---------

Alkaline Saline Peptone Water (ISO)

a high pH, high salt medium for the selective enrichment of *Vibrio* spp. from food samples; complies with ISO/TS 21872

500g	for 12.5 litres medium	CM1117B
------	------------------------	---------

ALOA® – see Chromogenic Listeria Agar (ISO)

Amies Transport Medium

for the preservation of organisms in transit

250g	for 12.5 litres medium	CM0425J
------	------------------------	---------

Ampicillin Supplement

10 vials	25mg/vial	SR0136E
----------	-----------	---------

★ Anaerobe Basal Agar

for the growth of fastidious and other anaerobes; can be made selective using supplements

500g	for 10.9 litres medium	CM0972B
Neomycin Supplement 150mg/vial		
10 vials	for 2.0 litres medium	SR0163H

Anaerobe Basal Broth

a medium for general growth of anaerobes

500g	for 14.1 litres medium	CM0957B
------	------------------------	---------

Andrade's Peptone Water – see Peptone Water (Andrade)

Antibiotic Medium No. 1 (Seed agar)

a medium recommended for the seed layer in preparation of plates for the microbiological assay of antibiotics

500g	for 18.5 litres medium	CM0327B
------	------------------------	---------

Antibiotic Medium No. 3 (Assay broth)

a medium for use in the serial dilution assay of penicillin and other antibiotics

500g	for 28.6 litres medium	CM0287B
------	------------------------	---------

Arcobacter Broth

for the selective growth of *Arcobacter* spp.

Arcobacter Broth Base		
500g	for 20.8 litres medium	CM0965B
CAT Supplement		
10 vials	for 500mL medium	SR0174E
10 vials	for 2.0 litres medium	SR0174H
OR		
CCDA Selective Supplement		
10 vials	for 500mL medium	SR0155E
10 vials	for 2.0 litres medium	SR0155H
10 vials	for 9.0 litres medium	SR0155N

Aspergillus Medium – see AFPA

Azide Blood Agar Base

a selective medium for the detection and isolation of streptococci and staphylococci from faeces, sewage and other specimens

500g	for 16.7 litres medium	CM0259B
------	------------------------	---------

Azide Dextrose Broth (Rothe Broth)

for the detection of enterococci in water, sewage and foods

500g	for 14.0 litres medium	CM0868B
------	------------------------	---------

Bacillus cereus Selective Agar (PEMBA)

a selective and diagnostic medium for the isolation and enumeration of *Bacillus cereus*. This medium can be made more selective by the addition of cycloheximide solution.

Bacillus cereus Selective Agar Base		
500g	for 12.2 litres medium	CM0617B
2.5kg	for 61.0 litres medium	CM0617R
Polymyxin B Supplement		
10 vials	for 500mL medium	SR0099E
Egg Yolk Emulsion		
100mL		SR0047C
Cycloheximide 0.1% Solution		
100mL		SR0222C

Bacillus cereus MYP Agar – see MYP Agar

Baird-Parker Egg Yolk Agar (ISO)

a selective medium for the isolation and enumeration of coagulase-positive staphylococci from food, with formulation that complies with ISO 6888-1:1999

Baird-Parker Agar (ISO) Base for use with Egg Yolk Tellurite Emulsion. Do not use with RPF supplement		
500g	for 7.9 litres medium	CM1127B
5kg	for 79 litres medium	CM1127T
Egg Yolk Tellurite Emulsion		
100mL		SR0054C
OR		
Egg Yolk Emulsion		
100mL		SR0047C
Potassium Tellurite 3.5%		
10x2mL		SR0030J

Baird-Parker Egg Yolk Tellurite Agar

selective medium for the isolation and enumeration of coagulase-positive staphylococci from food

Baird-Parker Agar Base for use with Egg Yolk Tellurite Emulsion. Do not use with RPF supplement		
500g	for 7.9 litres medium	CM0275B
2.5kg	for 39.5 litres medium	CM0275R
5kg	for 79 litres medium	CM0275T
Egg Yolk Tellurite Emulsion		
100mL		SR0054C
OR		
Egg Yolk Emulsion		
100mL		SR0047C
Potassium Tellurite 3.5%		
10x2mL		SR0030J

Baird-Parker RPF Agar

selective medium for the isolation and enumeration of coagulase-positive staphylococci from food

Baird-Parker Agar Base (RPF) for use with RPF supplement. Do not use with Egg Yolk Tellurite Emulsion.		
500g	for 7.1 litres medium	CM0961B
RPF Supplement only to be used with Baird-Parker Agar Base (RPF). Store below 0°C.		
10 vials	for 100mL medium	SR0122A

BCYE – see Legionella BCYE Agar

BiGGY Agar (Nickerson Medium)

for the isolation and presumptive identification of *Candida* spp.

500g	for 11.9 litres medium	CM0589B
------	------------------------	---------

Bile Aesculin Agar

a differential medium for the isolation and presumptive identification of enterococci

500g	for 11.2 litres medium	CM0888B
------	------------------------	---------

Bismuth Sulphite Agar (Modified Wilson & Blair Medium)

a medium for the isolation of *Salmonella* Typhi and other salmonellae. Particularly useful for the isolation of lactose-fermenting salmonellae. **N.B.** maximum shelf life of 12 months

500g	for 25.0 litres medium	CM0201B
------	------------------------	---------

* Blaser-Wang Selective Medium

a *Campylobacter* medium which can inhibit the growth of *Candida albicans*; can be supplemented with Campylobacter Growth Supplement

Columbia Blood Agar Base		
500g	for 12.8 litres medium	CM0331B
OR		
Blood Agar Base No. 2		
500g	for 12.5 litres medium	CM0271B
Campylobacter Selective Supplement (Blaser-Wang)		
10 vials	for 500mL medium	SR0098E
Campylobacter Growth Supplement		
10 vials	for 500mL medium	SR0232E

BLEB – see Buffered Listeria Enrichment Broth

* Blood Agar Base

a nonselective, general purpose medium which may be enriched with blood or serum

500g	for 12.5 litres medium	CM0055B
5kg	for 125.0 litres medium	CM0055T

* Blood Agar Base No. 2

an improved Blood Agar Base possessing enhanced nutritional properties suitable for the cultivation of fastidious organisms

500g	for 12.5 litres medium	CM0271B
2.5kg	for 62.5 litres medium	CM0271R
5kg	for 125.0 litres medium	CM0271T

Blood Agar Base Sheep Blood – see Sheep Blood Agar Base

Blood Agar Base with Tryptose – see Tryptose Blood Agar Base

BMPA - see Legionella BMPA-α Agar

* Bolton Selective Enrichment Broth

for selective enrichment of *Campylobacter* spp. in food samples; can be supplemented with Campylobacter Growth Supplement

Bolton Broth Base		
500g	for 18.1 litres medium	CM0983B
Bolton Broth Selective Supplement		
10 vials	for 500mL medium	SR0183E
OR		
Modified Bolton Broth Selective Supplement alternative supplement containing amphotericin B as a replacement for cycloheximide		
10 vials	for 500mL medium	SR0208E
Campylobacter Growth Supplement		
10 vials	for 500mL medium	SR0232E

* Bordetella Selective Medium

for the selective isolation of *Bordetella pertussis* and *B. parapertussis*

Charcoal Agar Base		
500g	for 9.8 litres medium	CM0119B
Bordetella Selective Supplement		
10 vials	for 500mL medium	SR0082E

Brain Heart Infusion Agar

highly nutritious medium, made using porcine brain, for the growth of fastidious organisms

500g	for 10.6 litres medium	CM1136B
------	------------------------	---------

Brain Heart Infusion Broth

highly nutritious medium, made using porcine brain, for the growth of fastidious organisms, recommended for blood culture work; also see Heart Infusion Broth

500g	for 13.5 litres medium	CM1135B
2.5kg	for 67.5 litres medium	CM1135R

Brewer's Thioglycollate Medium – see Thioglycollate Medium (Brewer)

Brilliance range of chromogenic agars see pages Chromogenic Media

Thermo Scientific™ *Brilliance*™ Chromogenic Dehydrated Culture Media

Brilliance *Bacillus cereus* Agar

a chromogenic medium for the isolation and differentiation of *Bacillus cereus* from food samples

<i>Brilliance</i> <i>Bacillus cereus</i> Agar Base		
500g	for 12.2 litres medium	CM1036B
<i>Brilliance</i> <i>Bacillus cereus</i> Selective Supplement		
10 vials	for 500mL medium	SR0230E

Brilliance CampyCount Agar - see page Chromogenic Media

Brilliance *Candida* Agar

a chromogenic medium for the isolation and differentiation of clinical *Candida* spp., including *C. albicans*

<i>Brilliance</i> <i>Candida</i> Agar Base		
500g	for 16.0 litres medium	CM1002B
<i>Brilliance</i> <i>Candida</i> Agar Selective Supplement		
10 vials	for 500mL medium	SR0231E

Brilliance *E. coli*/coliform Agar

for food and environmental samples

500g	for 9.0 litres medium	CM0956B
------	-----------------------	---------

Brilliance *E. coli*/coliform Selective Agar

for food, water and environmental samples

500g	for 17.8 litres medium	CM1046B
2.5kg	for 89.0 litres medium	CM1046R

Brilliance *Enterobacter sakazakii* Agar (DFI)

a chromogenic medium for the isolation and differentiation of *Cronobacter* spp. from food and dairy samples, according to the formulation by Druggan, Forsythe and Iverson

500g	for 11.6 litres medium	CM1055B
------	------------------------	---------

Brilliance ESBL Agar - see page Chromogenic Media

Brilliance *Listeria* Agar

a chromogenic medium for the selective growth and differentiation of *Listeria monocytogenes* and *Listeria* spp., and the isolation medium used in the **Listeria Precis™** method

<i>Brilliance</i> <i>Listeria</i> Agar Base		
500g	for 7.4 litres medium	CM1080B
5kg	for 74 litres medium	CM1080T
<i>Brilliance</i> <i>Listeria</i> Selective Supplement		
10 vials	for 500mL medium	SR0227E
<i>Brilliance</i> <i>Listeria</i> Differential Supplement		
10 vials	for 500mL medium	SR0228E
1 Vial	for 2.5L of medium	SR0228C

Brilliance MRSA 2 Agar - see page Chromogenic Media

Brilliance *Salmonella* Agar

a novel chromogenic medium incorporating Inhibigen technology, for the presumptive identification of *Salmonella* spp. from clinical and food samples, and the isolation medium used in the **Salmonella Precis** method

<i>Brilliance</i> <i>Salmonella</i> Agar Base		
500g	for 9.3 litres medium	CM1092B
5kg	for 93 litres medium	CM1092T
<i>Salmonella</i> Selective Supplement		
10 vials	for 500mL medium	SR0194E

Brilliance Staph 24 Agar - see page Chromogenic Media

Brilliance UTI Agar

a chromogenic medium to aid differentiation of organisms found in urinary tract infections (UTI)

400g	for 8.9 litres medium	CM0949C
------	-----------------------	---------

Brilliance VRE Agar - see page 19

Brilliance UTI Clarity Agar

an improved, clear, chromogenic UTI medium for the presumptive identification and differentiation of all the main microorganisms that cause UTIs

500g	for 13.5 litres medium	CM1106B
------	------------------------	---------

Brilliant Green Agar (Kauffmann Medium)

a selective medium for the isolation of salmonellae other than *Salmonella* Typhi. This can be made more selective by the addition of sulphamandelate.

500g	for 10.0 litres medium	CM0263B
------	------------------------	---------

Brilliant Green Agar (Modified) (Edel-Kampelmacher Medium)

a selective and diagnostic agar for salmonellae (other than *Salmonella* Typhi) from food and feeds. This can be made more selective by the addition of sulphamandelate.

500g	for 9.6 litres medium	CM0329B
2.5kg	for 48.0 litres medium	CM0329R
5kg	for 96.0 litres medium	CM0329T

Brilliant Green Bile 2% Broth

a medium for the detection or confirmation of members of the coli-aerogenes group; recommended for the 44°C confirmatory test for *E. coli*


500g	for 12.5 litres medium	CM0031B
------	------------------------	---------

See next page to find out more about our
Thermo Scientific™ *Brilliance*™ Chromogenic Media

Brilliance Chromogenic Media Plates

The chromogens used in our Thermo Scientific *Brilliance* chromogenic media ensure that colony colours are vivid, allowing rapid, easy differentiation and presumptive identification of the organism in question. Depending on media type, these are available in dehydrated and ready-to-use formats.


For Clinical Use



Brilliance CRE Agar

A chromogenic screening plate for the detection of carbapenem resistant Enterobacteriaceae from clinical samples.

Brilliance CRE Agar (ready poured plates)
10x90mm PO1226A



Brilliance ESBL Agar

Use this medium to identify ESBL-producing *E. coli* and the *Klebsiella*, *Enterobacter*, *Serratia* and *Citrobacter* group (KESC), direct from clinical samples in 24 hours.

Brilliance ESBL Agar (ready poured plates)
10x90mm PO5302A



Brilliance MRSA 2 Agar

Offering excellent sensitivity with the simplicity of a chromogenic plate, use this medium to detect MRSA in 24 hours.

Brilliance MRSA Agar (ready poured plates)
10x90mm PO1210A
10x90mm PO5310A




Brilliance UTI Clarity Agar

Use this clear medium for the differentiation and presumptive identification of common urinary tract infection organisms.

Brilliance UTI Clarity Agar (ready poured plates)
10x90mm PO1110A
10x90mm PO5159A

Brilliance UTI Clarity Agar
(dehydrated culture medium)
500g for 13.5 litres medium CM1106B




Brilliance UTI Agar

A chromogenic medium for the presumptive identification and differentiation of all the main microorganisms that cause urinary tract infections.

Brilliance UTI Agar (ready poured plates)
10x90mm PO0794A
10x90mm PO5120A

Brilliance UTI Agar (dehydrated culture medium)
400g for 8.9 litres medium CM0949C



Brilliance Candida Agar

A selective differential medium for the rapid isolation and identification of clinically important *Candida* spp.

Brilliance Candida Agar (ready poured plates)
10x90mm PO1034A
10x90mm PO5170A

Brilliance Candida Agar Base
(dehydrated culture medium)
500g for 16.0 litres medium CM1002B

Brilliance Candida Agar Selective Supplement
10 vials for 500mL medium SR0231E



Brilliance Listeria Agar

A medium for the selective growth and differentiation of *Listeria monocytogenes* and *Listeria* spp.

Brilliance Listeria Agar (ready poured plates)
10x90mm PO1102A
10x90mm PO5165A

Brilliance Listeria Agar Base
(dehydrated culture medium)
500g for 7.4 litres medium CM1080B
5kg for 74 litres medium CM1080T

Brilliance Listeria Selective Supplement
10 vials for 500mL medium SR0227E

Brilliance Listeria Differential Supplement
10 vials for 500mL medium SR0228E



Brilliance Salmonella Agar

A selective medium for the identification and differentiation of *Salmonella*.

Brilliance Salmonella Agar (ready poured plates)
10x90mm PO5098A

Brilliance Salmonella Agar Base
(dehydrated culture medium)
500g for 9.3 litres medium CM1092B

Salmonella Selective Supplement
10 vials for 500mL medium SR0194E



Brilliance VRE Agar

For the detection of vancomycin-resistant enterococci, use this medium to identify *E. faecium* and *E. faecalis*, direct from clinical samples in 24 hours.

Brilliance VRE Agar (ready poured plates)
10x90mm PO1175A




Brilliance GBS Agar

A clear, selective chromogenic medium for the presumptive identification of Group B Streptococci (*Streptococcus agalactiae*; GBS) from clinical samples.

Brilliance GBS Agar (ready poured plates)
10x90mm PO5320A

Chromogenic Media


For Food Laboratories



Brilliance CampyCount Agar

A chromogenic selective medium for the enumeration of *C. jejuni* and *C. coli* from poultry and related samples.


Brilliance CampyCount Agar (ready poured plates)
10x90mm PO1185A



Brilliance Enterobacter sakazakii Agar

For the differentiation and enumeration of *Cronobacter* spp. from infant formula and other food samples.


Brilliance Enterobacter sakazakii Agar [DFI]
(dehydrated culture medium)
100g for 2.3 litres medium CM1055A
500g for 11.6 litres medium CM1055B



Cronobacter sakazakii Isolation Agar

For the isolation and enumeration of *Cronobacter* (formerly *Enterobacter sakazakii*) from milk and milk products. This formulation complies with the formulation described in ISO/TS 22964:2006.

Cronobacter sakazakii Isolation Agar
(dehydrated culture medium)
500g for 18.0 litres medium CM1134B




Brilliance E. coli/coliform Selective Agar

For the detection and enumeration of *E. coli* and other coliforms from food and water samples.

Brilliance E. coli/coliform Selective Agar
(dehydrated culture medium)
500g for 17.8 litres medium CM1046B
2.5kg for 89.0 litres medium CM1046R


Brilliance E. coli/coliform Selective Agar
(ready poured plates)
10x90mm PO5176A



Brilliance Staph 24 Agar

A selective, chromogenic medium for the isolation and enumeration of coagulase-positive staphylococci in foods within 24 hours.


Brilliance Staph 24 Agar (ready poured plates)
10x90mm PO1186A



Brilliance Bacillus cereus Agar

Use for the isolation and differentiation of *B. cereus* from food and other samples.


Brilliance Bacillus cereus Agar Base
(dehydrated culture medium)
500g for 12.2 litres medium CM1036B
Brilliance Bacillus cereus Selective Supplement
10 vials for 500mL medium SR0230E



Chromogenic Coliform Agar

For the detection and enumeration of coliforms in food and water samples.

Chromogenic Coliform Agar (ready poured plates)
10x90mm PO5318A
10x55mm PO5428J



Brilliance Listeria Agar

For isolation, enumeration and identification of *Listeria* spp. and *Listeria monocytogenes* from food samples, and used as part of the Oxoid Listeria Precip method.

Brilliance Listeria Agar (ready poured plates)
10x90mm PO1102A
10x90mm PO5165A

Brilliance Listeria Agar Base
(dehydrated culture medium)
500g for 7.4 litres medium CM1080B
5kg for 7.4 litres medium CM1080T

Brilliance Listeria Selective Supplement
10 vials for 500mL medium SR0227E

Brilliance Listeria Differential Supplement
10 vials for 500mL medium SR0228E



Chromogenic Listeria Agar (ISO)

For isolation, enumeration and isolation *Listeria* species and *Listeria monocytogenes*, compliant with ISO 11290-1:1997 recommended formulation.

Chromogenic Listeria Agar (ready poured plates)
10x90mm PO5183A
10x90mm PO1196A

Chromogenic Listeria Agar Base
(dehydrated culture medium)
500g for 7.2 litres medium CM1084B

OCLA (ISO) Selective Supplement
10 vials for 500mL medium SR0226E

OCLA (ISO) Differential Supplement
(L- α -phosphatidylinositol solution)
10 vials for 500mL medium SR0244E



Brilliance Salmonella Agar

A selective medium for the identification and differentiation of *Salmonella*.

Brilliance Salmonella Agar (ready poured plates)
10x90mm PO5098A*

Brilliance Salmonella Agar Base
(dehydrated culture medium)
500g for 9.3 litres medium CM1092B

Salmonella Selective Supplement
10 vials for 500mL medium SR0194E



TBX Agar

Tryptone Bile Agar with X-Glucuronide, a chromogenic medium for the detection and differentiation of *E. coli*.

TBX Agar (ready poured plates)
10x90mm PO0727A
10x90mm PO5109A

TBX Agar (dehydrated culture medium)
500g for 13.7 litres medium CM0945B
2.5Kg for 68.5 litres medium CM0945R
5Kg for 137 litres medium CM0945T

Brucella Selective Medium

a selective medium for the cultivation and isolation of *Brucella* spp.

Blood Agar Base No. 2		
500g	for 12.5 litres medium	CM0271B
OR		
Brucella Medium Base		
500g	for 11.1 litres medium	CM0169B
OR		
Columbia Blood Agar Base		
500g	for 12.8 litres medium	CM0331B
Brucella Selective Supplement		
10 vials	for 500mL medium	SR0083A
OR		
Modified Brucella Selective Supplement alternative supplement containing natamycin as a replacement for cycloheximide		
10 vials	for 500mL medium	SR0209E
Horse Serum		
100mL		SR0035C

Buffered Glucose Brilliant Green Bile Broth – see EE Broth

Buffered Listeria Enrichment Broth (BLEB)

selective enrichment medium for the detection of *Listeria monocytogenes*, intended for use with samples of fermented products and as an alternative method for enrichment of environmental samples

Buffered Listeria Enrichment Broth		
500g	for 10.6 litres medium	CM0897B
2.5kg	for 53 litres medium	CM0897R
5kg	for 106.0 litres medium	CM0897T
Listeria Selective Enrichment Supplement		
10 vials	for 500mL medium	SR0141E
OR		
Modified Listeria Selective Enrichment Supplement alternative supplement containing amphotericin B as a replacement for cycloheximide		
10 vials	for 500mL medium	SR0213E

Buffered Charcoal Yeast Extract Agar – see Legionella BCYE Agar

Buffered Charcoal Yeast Extract Agar without L-Cysteine – see Legionella BCYE Agar without L-Cysteine

Buffered Peptone Water

500g	for 25.0 litres medium	CM0509B
2.5kg	for 125.0 litres medium	CM0509R
5kg	for 250.0 litres medium	CM0509T

Buffered Peptone Water (ISO)

a nonselective pre-enrichment broth for the isolation of *Salmonella* from food samples; the formulation complies with ISO 6579:2002

500g	for 25.0 litres medium	CM1049B
2.5kg	for 125.0 litres medium	CM1049R
5kg	for 250.0 litres medium	CM1049T

Buffered Sodium Chloride Peptone (EP/USP/JP/BP)

an isotonic diluent made to the harmonized Microbial Limit Testing chapters from the EP/USP/JP specification

500g	for 34.2 litres medium	CM0982B
------	------------------------	---------

Burkholderia cepacia Medium

a selective medium for the isolation of *Burkholderia cepacia*

Burkholderia cepacia Agar Base		
500g	for 13.7 litres medium	CM0995B
Burkholderia cepacia Selective Supplement		
10 vials	for 500mL medium	SR0189E

* Butzler Selective Medium

a selective medium for the isolation of *Campylobacter* spp. at 37°C; can be supplemented with Campylobacter Growth Supplement

Columbia Blood Agar Base		
500g	for 12.8 litres medium	CM0331B
Campylobacter Selective Supplement (Butzler)		
10 vials	for 500mL medium	SR0085E
Campylobacter Growth Supplement		
10 vials	for 500mL medium	SR0232E

Campylobacter Blood-Free Selective Medium
- see CCDA Medium

Campylobacter Growth Supplement

a liquid supplement for the enhanced growth and aerotolerance of *Campylobacter* spp.

10 vials	for 500mL medium	SR0232E
----------	------------------	---------

Candida Chromogenic Agar – see Brilliance Candida Agar

Cary-Blair Transport Medium

a transport medium for Gram-negative and anaerobic organisms

500g	for 37.6 litres medium	CM0519B
------	------------------------	---------

Casein Soya Bean Digest Agar - see Tryptone Soya Agar

Casein Soya Bean Digest Broth - see Tryptone Soya Broth

CAT Medium

an agar medium for the isolation of *Campylobacter* spp., especially *C. upsaliensis*

Campylobacter Blood-Free Selective Agar Base		
500g	for 11.0 litres medium	CM0739B
CAT Supplement		
10 vials	for 500mL medium	SR0174E
10 vials	for 2.0 litres medium	SR0174H

CCDA Medium

a medium for the isolation of *Campylobacter jejuni*, *C. coli*, *C. upsaliensis* and *C. lari*

Campylobacter Blood-Free Selective Agar Base		
500g	for 11.0 litres medium	CM0739B
2.5kg	for 55 litres medium	CM0739R
CCDA Selective Supplement		
10 vials	for 500mL medium	SR0155E
10 vials	for 2.0 litres medium	SR0155H

CDMN – see Clostridium difficile Moxalactam Norfloxacin Agar

Cefixime Rhamnose Sorbitol MacConkey Agar

selective and differential medium containing rhamnose and cefixime for the isolation of *E. coli* O157

Cefixime Rhamnose Sorbitol MacConkey Agar Base		
500g	for 8.8 litres medium	CM1005B
Cefixime Supplement		
10 vials	for 500mL medium	SR0191E

Cefixime Supplement

10 vials	0.025mg/vial	SR0191E
----------	--------------	---------

Cefixime Tellurite Sorbitol MacConkey Agar

for the isolation of *E. coli* O157

Sorbitol MacConkey Agar		
500g	for 9.7 litres medium	CM0813B
2.5kg	for 48.5 litres medium	CM0813R
C-T Supplement store below 0°C		
10 vials	for 500mL medium	SR0172E
10 vials	for 2.0 litres medium	SR0172H
10 vials	for 9.0 litres medium	SR0172N

Cephalexin – see Bordetella Selective Medium, Bordetella Selective Supplement

Cetrimide Agar – see Pseudomonas Cetrimide Agar Base

Cetrimide Selective Agar (EP/USP/JP/BP)

a solid medium recommended for use in quantitative procedures for selective and differential isolation of *Pseudomonas aeruginosa*

500g	11.0 litres medium	CM0579B
------	--------------------	---------

CFC Agar – see Pseudomonas CFC Agar

Chapman Medium – see Mannitol Salt Agar

Charcoal Agar Base – see Bordetella Selective Medium

China Blue Lactose Agar

for the differentiation and enumeration of bacteria in milk

500g	for 14.3 litres medium	CM0209B
------	------------------------	---------

Chloramphenicol Supplement

10 vials	50mg/vial	SR0078E
10 vials	200mg/vial	SR0078H

Cholera Medium TCBS

a selective isolation medium for pathogenic vibrios

500g	for 5.7 litres medium	CM0333B
------	-----------------------	---------

Christensen Agar – see Urea Agar

Christensen Broth – see Urea Broth

Chromogenic Cronobacter Isolation Agar (CCI)

a chromogenic medium for the isolation and enumeration of Cronobacter (formerly *E. sakazakii*) from foods

500g	for 16.34 litres medium	CM1122B
------	-------------------------	---------

Chromogenic Listeria Agar (ISO)

a selective, chromogenic agar for *Listeria* spp. and *L. monocytogenes*, compliant with the ISO 11290-1:1997 recommended formulation

Chromogenic Listeria Agar Base (ISO)		
500g	for 7.2 litres medium	CM1084B
OCLA (ISO) Selective Supplement		
10 vials	for 500mL medium	SR0226E
OCLA (ISO) Differential Supplement L- α -phosphatidylinositol solution		
10 vials	for 500mL medium	SR0244E

Chromogenic Salmonella Medium – see OSCM

Chromogenic TBX - see TBX Agar

Chromogenic (other) media - see *Brilliance* range, pages Chromogenic Media

CIN Agar – see Yersinia Selective Agar

Clarks & Lubbs Medium – see MRVP Medium

Clausen Medium – Dithionite Thioglycollate (HS-T) Medium

originally taken from the Nordic Pharmacopoeia, for sterility testing

500g	for 12.5 litres medium	CM0353B
------	------------------------	---------

CLED Medium

Cystine-Lactose-Electrolyte Deficient medium for urinary bacteriology

500g	for 13.8 litres medium	CM0301B
2.5kg	for 69.0 litres medium	CM0301R

CLED Medium (with Andrade's indicator)

for urinary bacteriology

500g	for 13.8 litres medium	CM0423B
2.5kg	for 69.0 litres medium	CM0423R

* Clostridium difficile Selective Medium

for the isolation of *Clostridium difficile*

Clostridium difficile Agar Base		
500g	for 7.2 litres medium	CM0601B
Clostridium difficile Selective Supplement		
10 vials	for 500mL medium	SR0096E

* Clostridium difficile Moxalactam Norfloxacin Selective Agar

a more selective alternative to Clostridium difficile Selective Medium

Clostridium difficile Agar Base		
500g	for 7.2 litres medium	CM0601B
CDMN Supplement		
10 vials	for 500mL medium	SR0173E

* CNA Staph/Strep Selective Agar

a selective medium for staphylococci and streptococci

Columbia Blood Agar Base		
500g	for 12.8 litres medium	CM0331B
OR		
Blood Agar Base No. 2		
500g	for 12.5 litres medium	CM0271B
Staph/Strep Selective Supplement		
10 vials	for 500mL medium	SR0070E
10 vials	for 9.0 litres medium	SR0070N

CN Agar – see Pseudomonas CN Agar

★ COBA Selective Medium

for the selective isolation of streptococci

Columbia Blood Agar Base		
500g	for 12.8 litres medium	CM0331B
Streptococcus Selective Supplement (COBA)		
10 vials	for 500mL medium	SR0126E

Cocoa Salmonella Recovery (CSR) Broth

A single-step enrichment medium for the recovery and growth of low levels of stressed Salmonella species from chocolate raw materials and interim products, such as cocoa powder, liquor, mass and butter, as well as finished chocolate products

500g		CM1155B
------	--	---------

Cold Filterable Tryptone Soya Broth

a gamma-irradiated, cold filterable Tryptone Soya Broth suitable for microbiological media fill trials in the pharmaceutical industry

500g	for 16.7 litres medium	CM1065B
2.5kg	for 83.3 litres medium	CM1065R
5kg	for 167.0 litres medium	CM1065T
25kg	for 833.0 litres medium	CM1065K

Cold Filterable Vegetable Peptone Broth

a gamma-irradiated, cold filterable Vegetable Peptone Broth suitable for microbiological media fill trials in the pharmaceutical industry, made using raw materials free from animal-derived material

500g	for 16.1 litres medium	VG0704B
5kg	for 161.0 litres medium	VG0704T

Columbia Agar MLT (EP/USP/JP/BP)

recommended for use in qualitative procedures as a general purpose medium for cultivation of fastidious and nonfastidious organisms

500g	12.0 litres medium	R455412
------	--------------------	---------

★ Columbia Blood Agar Base

a multipurpose medium suitable for the cultivation of fastidious organisms

500g	for 12.8 litres medium	CM0331B
2.5kg	for 64.0 litres medium	CM0331R
5kg	for 128.0 litres medium	CM0331T

Columbia CNA Agar – see CNA Staph/Strep Selective Agar

Cooked Meat Medium

an excellent medium for the primary growth and maintenance of aerobic and anaerobic organisms

500g		CM0081B
------	--	---------

Corn Meal Agar

a mycological medium for the identification of *C. albicans* and the maintenance of fungal stock cultures

500g	for 29.4 litres medium	CM0103B
------	------------------------	---------

Cronobacter Screening Broth (CSB)

a differential enrichment broth, for use in conjunction with chromogenic agars, for the isolation of *Cronobacter* (formerly *E. sakazakii*) from foods

Cronobacter Screening Broth Base		
500g	for 17.8 litres medium	CM1121B
Vancomycin Supplement (5mg)		
10 vials	for 500mL medium	SR0247E

Crossley Milk Medium

a diagnostic medium for examination of food samples for anaerobes

500g	for 4.5 litres medium	CM0213B
------	-----------------------	---------

CR-SMAC – see Cefixime Rhamnose Sorbitol MacConkey Agar

CT-SMAC – see Cefixime Tellurite Sorbitol MacConkey Agar

Cycloheximide 0.1% Solution

a solution for addition to a range of Thermo Scientific media for the suppression of fungi

100mL		SR0222C
-------	--	---------

Czapek Dox Agar Modified

for the cultivation of fungi

500g	for 11.0 litres medium	CM0097B
------	------------------------	---------

Czapek Dox Liquid Medium

for the cultivation of fungi

500g	for 15.0 litres medium	CM0095B
------	------------------------	---------

DCA – see Desoxycholate Citrate Agar

DCA (Hynes) – see Desoxycholate Citrate Agar (Hynes)

DCLS Agar (Desoxycholate Lactose Sucrose)

a modified DCA, containing sucrose to improve the accuracy of recognition of pathogenic Enterobacteriaceae

500g	for 10.0 litres medium	CM0393B
------	------------------------	---------

D-Cycloserine – see Tryptose Cycloserine Agar, Perfringens (TSC) Supplement

De Man Rogosa Sharpe Agar – see MRS Agar

Dermasel Medium

for the primary isolation and identification of dermatophyte fungi

Dermasel Agar Base		
500g	for 11.2 litres medium	CM0539B
Dermasel Selective Supplement		
10 vials	for 500mL medium	SR0075E
10 vials	for 2.0 litres medium	SR0075H

Desoxycholate Agar

a differential medium for the enumeration of coliforms and for the isolation of enteric pathogens

500g	for 11.1 litres medium	CM0163B
------	------------------------	---------



Desoxycholate Citrate Agar

a differential medium for the enumeration of coliforms and for the isolation of enteric pathogens

500g	for 10.3 litres medium	CM0035B
------	------------------------	---------

Desoxycholate Citrate Agar (Hynes)

for isolation of enteric pathogens

500g	for 9.6 litres medium	CM0227B
------	-----------------------	---------

Dextrose Tryptone Agar

for the detection and enumeration of thermophilic and mesophilic organisms in food

500g	for 18.5 litres medium	CM0075B
------	------------------------	---------

Dextrose Tryptone Broth

for the detection and enumeration of thermophilic and mesophilic organisms in food

500g	for 33.3 litres medium	CM0073B
------	------------------------	---------

Diagnostic Sensitivity Test Agar (DSTA)

a susceptibility test agar for antimicrobial susceptibility testing

500g	for 12.5 litres medium	CM0261B
2.5kg	for 62.5 litres medium	CM0261R
5kg	for 125.0 litres medium	CM0261T

Dichloran-Glycerol (DG18) (ISO) Agar

a selective, low water activity medium for xerophilic moulds from dried and semidried foods according to ISO 21527-2

Dichloran-Glycerol (DG18) (ISO) Agar base Requires chloramphenicol supplement SR0078E/SR0078H		
500g	for 15.8 litres medium	CM1150B
Chloramphenicol Supplement		
10 vials	for 500mL medium	SR0078E
10 vials	for 2.0 litres medium	SR0078H
Pre-Supplemented Dichloran-Glycerol (DG18) (ISO) Agar		
500g	for 15.8 litres medium	CM1151B

Dichloran-Glycerol (DG18) Selective Medium

a selective, low water activity medium for xerophilic moulds from dried and semidried foods

Dichloran-Glycerol (DG18) Agar Base		
500g	for 15.9 litres medium	CM0729B
Chloramphenicol Supplement		
10 vials	for 500mL medium	SR0078E
10 vials	for 2.0 litres medium	SR0078H

Dichloran Rose-Bengal Chloramphenicol Agar (DRBC)

a selective medium for yeasts and moulds associated with food spoilage

DRBC Agar Base		
500g	for 15.9 litres medium	CM0727B
Chloramphenicol Supplement		
10 vials	for 500 mL medium	SR0078E
10 vials	for 2.0 litres medium	SR0078H

Dichloran Rose-Bengal Chloramphenicol (ISO) Agar

a selective medium for yeasts and moulds associated with food spoilage, formulated according to ISO 21527-1

DRBC (ISO) Agar Base Requires chloramphenicol supplement SR0078E/SR0078H		
500g	for 15.8 litres medium	CM1148B
Chloramphenicol Supplement		
10 vials	for 500mL medium	SR0078E
10 vials	for 2.0 litres medium	SR0078H
Pre-supplemented DRBC (ISO) Agar		
500g	for 15.8 litres medium	CM1149B
2.5kg	for 79 litres medium	CM1149R

DNase Agar

for the detection of microbial deoxyribonuclease enzymes, particularly from staphylococci

500g	for 12.8 litres medium	CM0321B
------	------------------------	---------

DRBC – see Dichloran Rose-Bengal Chloramphenicol Agar

DSTA – see Diagnostic Sensitivity Test Agar

EC Broth

for the isolation of *E. coli* and other coliforms

500g	for 13.5 litres medium	CM0853B
------	------------------------	---------

EC Broth (Reduced Bile Salts)

a selective enrichment broth for the growth of *E. coli* O157 from food and environmental samples. Enhanced selectivity can be achieved with the addition of novobiocin.

500g	for 13.7 litres medium	CM0990B
2.5kg	for 68.5 litres medium	CM0990R
Novobiocin Supplement		
10 vials	for 500mL medium	SR0181E

EC Broth with MUG

for the isolation and presumptive identification of *E. coli*

100g	for 2.7 litres medium	CM0979A
500g	for 13.5 litres medium	CM0979B

Edel-Kampelmacher Medium – see Brilliant Green Agar (Modified)

★ Edwards Medium Modified

for the selective isolation of *Streptococcus agalactiae*

500g	for 12.2 litres medium	CM0027B
------	------------------------	---------

EE Broth (Buffered Glucose Brilliant Green Bile Broth)

an enrichment medium for bile-tolerant, Gram-negative bacteria

500g	for 11.5 litres medium	CM0317B
5kg	for 114.9 litres medium	CM0317T

Endo Agar

an APHA recommended medium for the confirmation of coliforms; also see Membrane Endo Agar LES

Endo Agar Base		
500g	for 13.9 litres medium	CM0479B
Basic Fuchsin Indicator		
10g		BR0050A

Enterobacteria Enrichment Broth-Mossel (EP/USP/JP/BP)

an enrichment medium for bile-tolerant, Gram-negative bacteria in the microbiological examination of pharmaceutical products

500g	for 19.2 litres medium	CM1115B
------	------------------------	---------

Cronobacter sakazakii Isolation Agar

A chromogenic medium for the isolation and enumeration of *Cronobacter* (formerly *E. sakazakii*) from milk and milk products, with formulation that complies with ISO/TS 22964:2006

500g	for 18.0 litres medium	CM1134B
------	------------------------	---------

Enterococcus Agar – see Slanetz and Bartley Medium

Eosin Methylene Blue Agar (Levine)

an isolation medium for the differentiation of Enterobacteriaceae

500g	for 13.3 litres medium	CM0069B
------	------------------------	---------

FBP Supplement – see Campylobacter Growth Supplement

Fraser Broth

a secondary, selective, diagnostic enrichment medium for the isolation of *Listeria* spp. from food and environmental samples

Fraser Broth Base		
500g	for 8.7 litres medium	CM0895B
2.5kg	for 43.6 litres medium	CM0895R
5kg	for 87.1 litres medium	CM0895T
Fraser Supplement		
10 vials	for 500mL medium	SR0156E

GBS Agar – see Group B Streptococci Medium

GC Media – see New York City Media and Thayer Martin Media

* Gardnerella vaginalis Selective Medium

for the isolation of *Gardnerella vaginalis*

Columbia Blood Agar Base		
500g	for 12.8 litres medium	CM0331B
Gardnerella vaginalis Selective Supplement		
10 vials	for 500mL medium	SR0119E

Gentamicin Supplement

10 vials	256mg/vial	SR0185E
----------	------------	---------

Giolitti-Cantoni Broth

for the enrichment of *Staphylococcus aureus* from foods

Giolitti-Cantoni Broth Base		
500g	for 9.2 litres medium	CM0523B
Potassium Tellurite 3.5%		
10x2mL		SR0030J

Group B Streptococci Medium

for the isolation and detection of Group B streptococci in clinical specimens

GBS Agar Base (Islam)		
500g	for 11.1 litres medium	CM0755B
Horse Serum		
100mL		SR0035C

GVPC Selective Medium – see Legionella GVPC Selective Medium

GVPN Selective Medium – see Legionella GVPN Selective Medium

Haemophilus Test Medium (HTM)

for the susceptibility testing of *Haemophilus influenzae*

HTM Base		
500g	for 11.6 litres medium	CM0898B
HTM Supplement store at -20°C to 0°C		
10 vials	for 500mL medium	SR0158E

Half Fraser Broth

a primary selective enrichment broth for *Listeria* spp.

Fraser Broth Base		
500g	for 8.7 litres medium	CM0895B
2.5kg	for 43.6 litres medium	CM0895R
5kg	for 87.1 litres medium	CM0895T
Half Fraser Supplement		
10 vials	for 225mL medium	SR0166E
10 vials	for 2.25 litres medium	SR0166G

Half Fraser Broth (without Ferric Ammonium Citrate)

a pre-supplemented primary selective enrichment broth for *Listeria* spp., no further additions, except Ferric Ammonium Citrate, if required

500g	for 8.7 litres medium	CM1053B
2.5kg	for 43.6 litres medium	CM1053R
10kg	for 174.4 litres medium	CM1053G

Heart Infusion Broth

a highly nutritious, liquid medium recommended for the cultivation of fastidious organisms. Also see Brain Heart Infusion Broth

500g	for 13.5 litres medium	CM1032B
------	------------------------	---------

Hektoen Enteric Agar

a differential selective medium for the isolation of *Shigella* and *Salmonella* spp.

500g	for 6.6 litres medium	CM0419B
------	-----------------------	---------

* Helicobacter pylori Selective Medium

for the selective isolation of *Helicobacter pylori*

Columbia Blood Agar Base		
500g	for 12.8 litres medium	CM0331B
Helicobacter pylori Selective Supplement (Dent)		
10 vials	for 500mL medium	SR0147E

* Hoyle Medium

for the isolation and differentiation of all types of *Corynebacterium diphtheriae*

Hoyle Medium Base		
500g	for 12.5 litres medium	CM0083B
Potassium Tellurite 3.5%		
10x2mL		SR0030J

Iron Sulphite Agar

for the detection of thermophilic anaerobic organisms

500g	for 21.7 litres medium	CM0079B
------	------------------------	---------

Irradiated Tryptone Soya Broth – see Cold Filterable Tryptone Soya Broth

Irradiated Vegetable Peptone Broth – see Cold Filterable Vegetable Peptone Broth

Islam Medium – see Group B Streptococci Medium

Iso-Sensitest Agar

a semidefined medium designed for antimicrobial susceptibility testing

500g	for 15.9 litres medium	CM0471B
2.5kg	for 79.5 litres medium	CM0471R
5kg	for 159.0 litres medium	CM0471T

Iso-Sensitest Broth

a semidefined medium designed for antimicrobial susceptibility testing

500g	for 21.4 litres medium	CM0473B
------	------------------------	---------

Kanamycin Aesculin Azide Agar

a selective medium for the isolation of enterococci in foodstuffs

Kanamycin Aesculin Azide Agar Base		
500g	for 11.7 litres medium	CM0591B
2.5kg	for 58.5 litres medium	CM0591R
Kanamycin Sulphate Supplement 10mg/vial		
10 vials	for 500mL medium	SR0092E

Karmali Selective Medium

a blood-free, selective medium for the isolation of *Campylobacter jejuni* and *Campylobacter coli* when incubated at 42°C

Campylobacter Agar Base (Karmali)		
500g	for 11.6 litres medium	CM0935B
Campylobacter Selective Supplement (Karmali)		
10 vials	for 500mL medium	SR0167E
OR		
Modified Karmali Selective Supplement alternative supplement containing amphotericin B as a replacement for cycloheximide		
10 vials	for 500mL medium	SR0205E

Kauffmann Medium – see Brilliant Green Agar

K-F Streptococcus Agar

a selective medium for the isolation and enumeration of Group D streptococci

K-F Streptococcus Agar Base		
500g	for 6.5 litres medium	CM0701B
TTC Solution (1%)		
10x5mL	for 500mL medium	SR0229K
OR		
TTC Solution (5%)		
100mL		SR0211C

Kligler Iron Agar

a differential medium for the identification of the Enterobacteriaceae

500g	for 9.1 litres medium	CM0033B
------	-----------------------	---------

Lab-Lemco Agar

nutrient medium equivalent to USA nutrient extract formulations

500g	for 21.7 litres medium	CM0017B
------	------------------------	---------

Lab-Lemco Broth

nutrient medium equivalent to USA nutrient extract formulations

500g	for 62.5 litres medium	CM0015B
------	------------------------	---------

Lactic Acid 10%

10x1mL		SR0021H
10x5mL		SR0021K

Lactose Broth

for the presumptive identification of *E. coli* and coliforms

500g	for 38.5 litres medium	CM0137B
2.5kg	for 192.5 litres medium	CM0137R

Lauryl Tryptose Broth (Lauryl Sulphate Broth)

an APHA standard methods medium for the examination of water and waste water

500g	for 14.0 litres medium	CM0451B
2.5kg	for 70.0 litres medium	CM0451R

Lauryl Tryptose Broth with MUG

for the presumptive identification of *E. coli*

500g	for 14.0 litres medium	CM0980B
------	------------------------	---------

Lauryl Tryptose Broth (Modified) with MUG

for the enumeration of presumptive *E. coli* and other coliforms from water, dairy and food samples; with additional tryptophan

500g	for 13.6 litres medium	CM0967B
------	------------------------	---------

Legionella BCYE Agar

a growth medium for the isolation of Legionellaceae

Legionella CYE Agar Base		
500g	for 20.0 litres medium	CM0655B
2.5kg	for 100 litres medium	CM0655R
Legionella BCYE Growth Supplement		
10 vials	for 100mL medium	SR0110A
10 vials	for 500mL medium	SR0110C

Legionella BCYE Agar without L-Cysteine

to be used in conjunction with BCYE medium as a confirmatory test for *Legionella* spp.

Legionella CYE Agar Base		
500g	for 20.0 litres medium	CM0655B
2.5kg	for 100 litres medium	CM0655R
Legionella BCYE Growth Supplement without L-cysteine		
10 vials	for 100mL medium	SR0175A

Legionella BMPA- α Selective Medium

a semiselective medium for the isolation of *Legionella pneumophila* from clinical and environmental specimens

Legionella CYE Agar Base		
500g	for 20.0 litres medium	CM0655B
2.5kg	for 100 litres medium	CM0655R
Legionella BMPA- α Growth Selective Supplement		
10 vials	for 100mL medium	SR0111E
10 vials	for 500mL medium	SR0111B
Legionella BCYE Growth Supplement		
10 vials	for 100mL medium	SR0110A
10 vials	for 500mL medium	SR0110C

Legionella GVPC Selective Medium

a selective medium for the isolation of *Legionella* spp. from environmental samples

Legionella CYE Agar Base		
500g	for 20.0 litres medium	CM0655B
2.5kg	for 100 litres medium	CM0655R
Legionella GVPC Selective Supplement		
10 vials	for 500mL medium	SR0152E
Legionella BCYE Growth Supplement		
10 vials	for 100mL medium	SR0110A
10 vials	for 500mL medium	SR0110C

Legionella GVPN Selective Medium

an alternative to GVPC selective medium, containing natamycin as a replacement for cycloheximide, for the selective isolation of *Legionella* spp.

Legionella CYE Agar Base		
500g	for 20.0 litres medium	CM0655B
2.5kg	for 100 litres medium	CM0655R
Legionella GVPN Selective Supplement		
10 vials	for 500mL medium	SR0215E
Legionella BCYE Growth Supplement		
10 vials	for 100mL medium	SR0110A
10 vials	for 500mL medium	SR0110C

Legionella MWY Selective Medium

a selective medium for the isolation of *Legionella pneumophila* from potable water samples

Legionella CYE Agar Base		
500g	for 20.0 litres medium	CM0655B
2.5kg	for 100 litres medium	CM0655R
Legionella MWY Selective Supplement		
10 vials	for 100mL medium	SR0118E
10 vials	for 500mL medium	SR0118B
Legionella BCYE Supplement		
10 vials	for 100mL medium	SR0110A
10 vials	for 500mL medium	SR0110C

Levine Medium – see Eosin Methylene Blue Agar

Liquid Yeast Autolysate

a source of nitrogen and vitamins suitable for addition to a range of culture media

10 vials	for 2.0 litres medium	SR0182H
----------	-----------------------	---------

Listeria Chromogenic Agar – see *Brilliance* Listeria Agar
Listeria Chromogenic Agar (ISO) – see Chromogenic Listeria Agar (ISO)
Listeria Novel Enrichment Broth – see ONE Broth-Listeria

Listeria Precis Method

Thermo Scientific™ Oxoid™ Listeria Precis™ Method is a rapid culture method for the enrichment, isolation and identification of *Listeria monocytogenes* and *Listeria* spp. in just 2 days. Validated by AFNOR to ISO 16140 for all food, animal feed and environmental samples; use Thermo Scientific Oxoid O.B.I.S Mono ID System (ID0600M) for rapid confirmation.

ONE Broth-Listeria

for selective enrichment of *Listeria* spp. from food samples in 24 hours

ONE Broth-Listeria Base		
500g	for 11.4 litres medium	CM1066B
2.5kg	for 57.0 litres medium	CM1066R
5kg	for 114 litres medium	CM1066T
ONE Broth-Listeria Selective Supplement		
10 vials	for 500mL medium	SR0234E
10 vials	for 2.25 litres medium	SR0234B
10 vials	for 2 litres medium	SR0234H
ONE Broth-Listeria Complete		
5x20L Dry-Bag	for 20 litres of medium	DB1066V

Brilliance Listeria Agar

a chromogenic agar for the selective growth and differentiation of *Listeria monocytogenes* and *Listeria* spp.

Brilliance Listeria Agar Base		
500g	for 7.4 litres medium	CM1080B
5kg	for 74 litres medium	CM1080T
Brilliance Listeria Selective Supplement		
10 vials	for 500mL medium	SR0227E
Brilliance Listeria Differential Supplement		
10 vials	for 500mL medium	SR0228E
10 vials	for 2.5 L medium	SR0228C

Listeria Selective Agar (Oxford Formulation)

a selective and diagnostic medium for the detection of *Listeria* spp.

Listeria Selective Agar Base (Oxford)		
500g	for 9.0 litres medium	CM0856B
2.5kg	for 45.0 litres medium	CM0856R
5kg	for 90 litres medium	CM0856T
Listeria Selective Supplement		
10 vials	for 500mL medium	SR0140E
10 vials	for 2.5 litres medium	SR0140B
OR		
Modified Listeria Selective Supplement (Oxford) alternative supplement containing amphotericin B as a replacement for cycloheximide		
10 vials	for 500mL medium	SR0206E
10 vials	for 2.5 litres medium	SR0206B

Listeria Selective Enrichment Medium

this medium complies with the formulation often referred to as the FDA method

Listeria Enrichment Broth		
500g	for 13.9 litres medium	CM0862B
Listeria Selective Enrichment Supplement		
10 vials	for 500mL medium	SR0141E
OR		
Listeria Selective Enrichment Supplement (Modified) includes cycloheximide; modified with 10mg/litre acriflavine		
10 vials	for 2.25 litres medium	SR0149A

Listeria Selective Enrichment Media (UVM Formulations)

a two step, selective enrichment (USDA-FSIS) method

Listeria Enrichment Broth - UVM 1 Broth (complete)		
500g		CM1054B
OR		
Listeria Enrichment Broth Base (UVM formulation) use with appropriate supplement for UVM I broth or UVM II broth		
500g	for 9.2 litres medium	CM0863B
Listeria Primary Selective Enrichment Supplement (UVM I)		
10 vials	for 500mL medium	SR0142E
10 vials	for 2.25 litres medium	SR0142B
OR		
Listeria Secondary Selective Enrichment Supplement (UVM II)		
10 vials	for 500mL medium	SR0143E

Listeria Selective Agar PALCAM – see PALCAM Agar

Liver Broth

for the examination of foods for saccharolytic or putrefactive mesophilic and thermophilic anaerobes

500g	for 7.8 litres medium	CM0077B
------	-----------------------	---------

Lysine Decarboxylase Broth Tablets (Taylor modification)

for the differentiation of the Enterobacteriaceae. Each tablet makes 5mL of medium; store at 2–8°C

100 tabs	for 0.5 litres medium	CM0308S
----------	-----------------------	---------

Lysine Iron Agar

a diagnostic medium for salmonellae, including *Salmonella arizonae*

500g	for 14.7 litres medium	CM0381B
------	------------------------	---------

Lysine Medium

for the isolation and enumeration of wild yeasts in pitching yeast

Lysine Medium Base		
500g	for 7.6 litres medium	CM0191B
Potassium Lactate		
10x1mL		SR0037H
Lactic Acid 10%		
10x1mL		SR0021H
10x5mL		SR0021K

M17 Agar

for improved growth of lactic streptococci and their bacteriophages, and the selective enumeration of *Streptococcus thermophilus* from yoghurt

500g	for 10.4 litres medium	CM0785B
Lactose Bacteriological		
500g		LP0070B

M17 Broth

for improved growth of lactic streptococci and their bacteriophages, and the selective enumeration of *Streptococcus thermophilus* from yoghurt

500g	for 13.4 litres medium	CM0817B
Lactose Bacteriological		
500g		LP0070B

MacConkey Agar

a differential medium for the isolation of coliforms and intestinal pathogens in water, dairy products and biological specimens

500g	for 9.6 litres medium	CM0007B
2.5kg	for 48.0 litres medium	CM0007R

MacConkey Agar (EP/USP/JP/BP)

a solid medium recommended for use in qualitative procedures for selection and differential isolation of Gram-negative bacilli

500g	for 10.0 litres medium	R453802
------	------------------------	---------

MacConkey Agar (without salt)

a differential medium for urine examination on which swarming of *Proteus* spp. is suppressed

500g	for 10.6 litres medium	CM0507B
2.5kg	for 53.2 litres medium	CM0507R
5kg	for 106 litres medium	CM0507T

MacConkey Agar No. 2

a modification of MacConkey Agar No. 3 containing Bile Salts No. 2 for the recognition of enterococci

500g	for 9.7 litres medium	CM0109B
------	-----------------------	---------

MacConkey Agar No. 3

a selective medium giving excellent differentiation between lactose and non-lactose fermenters with the inhibition of Gram-positive cocci

500g	for 9.7 litres medium	CM0115B
5kg	for 97.0 litres medium	CM0115T

MacConkey Broth

a differential medium containing neutral red for the detection of coliform organisms in water and milk

500g	for 12.5 litres medium	CM0005B
------	------------------------	---------

MacConkey Broth (EP/USP/JP/BP)

a liquid medium recommended for use in qualitative procedures for selective isolation and presumptive identification of Enterobacteriaceae

500g	for 14.2 litres medium	R453822
------	------------------------	---------

MacConkey Broth Purple (US Formulation)

a differential medium containing bromocresol purple for the detection of coliforms

500g	for 12.5 litres medium	CM0505B
2.5kg	for 62.5 litres medium	CM0505R

Malt Extract Agar

for the detection, isolation and enumeration of yeasts and moulds

500g	for 10.0 litres medium	CM0059B
5kg	for 100 litres medium	CM0059T

Malt Extract Broth

for the cultivation of moulds and yeasts

500g	for 25.0 litres medium	CM0057B
------	------------------------	---------

Mannitol Egg Yolk Polymyxin Agar – see MYP Agar

Mannitol Salt Agar (Chapman Medium)

for the isolation of pathogenic staphylococci

500g	for 4.5 litres medium	CM0085B
2.5kg	for 22.5 litres medium	CM0085R
5kg	for 45.0 litres medium	CM0085T

Mannitol Salt Agar (EP/USP/JP/BP)

a solid medium recommended for use in qualitative procedures for selective and differential isolation of staphylococci

500g	for 4.5 litres medium	R453902
------	-----------------------	---------

Mannitol Selenite Broth

a modification of Selenite F Broth for the selective enrichment of salmonellae

Mannitol Selenite Broth Base		
500g	for 26.3 litres medium	CM0399B
Sodium Biselenite		
100g	for 25.0 litres medium	LP0121A

Maximum Recovery Diluent (MRD)

a protective isotonic diluent

500g	for 52.6 litres medium	CM0733B
2.5kg	for 263.0 litres medium	CM0733R
5kg	for 526.0 litres medium	CM0733T

Membrane Clostridium perfringens Medium (m-CP)

a chromogenic medium for rapid isolation and presumptive identification of *Clostridium perfringens* from water samples

Membrane Clostridium perfringens (m-CP) Agar Base		
500g	for 7.0 litres medium	CM0992B
Membrane Clostridium perfringens (m-CP) Selective Supplement		
10 vials	for 500mL medium	SR0188E

Membrane Endo Agar LES

an APHA recommended medium for the confirmation of coliforms

Membrane Endo Agar LES		
500g	for 11.1 litres medium	MM0551B
Basic Fuchsin Indicator		
10g		BR0050A

Membrane Lactose Glucuronide Agar (MLGA)

a selective chromogenic medium for differentiation of *E. coli* and other coliforms

500g	for 5.7 litres medium	CM1031B
------	-----------------------	---------

Membrane Lauryl Sulphate Broth

a direct replacement for Membrane Enriched Teepol Broth requiring no supplementation

500g	for 6.6 litres medium	MM0615B
------	-----------------------	---------

Meropenem Supplement

10 vials	1mg/vial	SR0184E
----------	----------	---------

Milk Agar

plate count agar for use in the dairy industry

500g	for 20.8 litres medium	CM0021B
2.5kg	for 104.0 litres medium	CM0021R

Milk Plate Count Agar

with antibiotic free skimmed milk

500g	for 25.6 litres medium	CM0681B
2.5kg	for 128.0 litres medium	CM0681R

Minerals Modified Glutamate Medium

for the enumeration of the coliform group in water (supplied as a twin pack)

Minerals Modified Glutamate Medium Base plus Sodium Glutamate		
500g plus 300g	for 22.0 litres medium	CM0607G

MKTT-n – see Muller-Kauffmann Tetrathionate Novobiocin Broth

MLCB Agar

for the selective isolation of salmonellae (not *Salmonella* Typhi or *Salmonella* Paratyphi A)

500g	for 10.2 litres medium	CM0783B
------	------------------------	---------

MLGA – see Membrane Lactose Glucuronide Agar

MSB – see Membrane Lauryl Sulphate Broth

Modified Brilliant Green Agar – see Brilliant Green Agar (Modified)

Modified Lauryl Sulphate Tryptose Vancomycin Broth

for selective enrichment of *Cronobacter* (formerly *E. sakazakii*) from milk and milk products, with formulation that complies with ISO/TS 22964:2006

Modified Lauryl Sulphate Tryptose Broth Base (mLST)		
500g	for 7.74 litres medium	CM1133B
Vancomycin Supplement (5mg)		
10 vials	for 500mL medium	SR0247E

Modified Lauryl Tryptose Broth with MUG – see Lauryl Tryptose Broth (Modified) with MUG

Modified New York City Media – see New York City Media

Modified Semi-solid Rappaport Vassiliadis (MSRV) Medium

for detection of motile *Salmonella* spp. from food and environmental samples

MSRV Medium Base N.B. maximum shelf life of 12 months		
500g	for 15.8 litres medium	CM0910B
MSRV Selective Supplement		
10 vials	for 500mL medium	SR0161E

Modified Semi-solid Rappaport Vassiliadis (MSRV) Medium (ISO)

a semisolid medium for the detection of *Salmonella* spp. from food and environmental samples, with formulation conforming to that recommended in ISO 6579 Annex D

MSRV (ISO) Medium Base N.B. maximum shelf life of 12 months		
500g	for 15.8 litres medium	CM1112B
MSRV Selective Supplement		
10 vials	for 1 litre medium	SR0161E

Modified Tryptone Soya Broth – see Tryptone Soya Broth (Modified)

Modified Wadowsky Yee Medium – see Legionella MWY Selective Medium

MRD – see Maximum Recovery Diluent

MRS Agar (De Man, Rogosa, Sharpe)

for the growth of lactobacilli

500g	for 8.1 litres medium	CM0361B
2.5kg	for 40.5 litres medium	CM0361R
5kg	for 81 litres medium	CM0361T

MRS (De Man, Rogosa, Sharpe) (ISO) Agar

for the growth of lactobacilli. Complies with ISO 15214:1998

500g	for 7.5 litres medium	CM1153B
2.5kg	for 37.5 litres medium	CM1153R

MRS Broth (De Man, Rogosa, Sharpe)

for the growth of lactobacilli

500g	for 9.6 litres medium	CM0359B
------	-----------------------	---------

MRVP Medium (Clarks and Lubs Medium)

a glucose-phosphate medium for the differentiation of the coli-aerogenes group by the Methyl Red and Voges-Proskauer tests

500g	for 29.4 litres medium	CM0043B
------	------------------------	---------

MSRV – see Modified Semi-solid Rappaport Vassiliadis Medium

MSRV (ISO) – see Modified Semi-solid Rappaport Vassiliadis Medium (ISO)

m-TSB – see Tryptone Soya Broth (Modified)

Mueller-Hinton Agar

conforming to the M6-A2 CLSI international standard for antimicrobial susceptibility testing media

500g	for 13.2 litres medium	CM0337B
2.5kg	for 66.0 litres medium	CM0337R
5kg	for 132.0 litres medium	CM0337T

Mueller-Hinton Broth

conforming to international standards for antimicrobial susceptibility testing media

500g	for 23.8 litres medium	CM0405B
------	------------------------	---------

Muller-Kauffmann Tetrathionate Broth Base

a selective enrichment medium for the isolation of *Salmonella* and the suppression of *Proteus* spp.

500g	for 6.1 litres medium	CM0343B
2.5kg	for 30.5 litres medium	CM0343R

Muller-Kauffmann Tetrathionate Novobiocin Broth (MKTTn)

a selective enrichment broth for the isolation of *Salmonella* from food samples; complies with ISO 6579:2002

MKTTn Broth Base		
500g	for 5.8 litres medium	CM1048B
Novobiocin Supplement		
10 vials	for 250mL medium	SR0181E

MUG Supplement

4-methylumbeliferyl-β-D-glucuronide; a fluorogenic substrate for the detection of *Escherichia coli*

10 vials	50mg/vial	BR0071E
----------	-----------	---------

Mycoplasma Agar

for the selective isolation of *Mycoplasma* spp.

Mycoplasma Agar Base		
500g	for 14.1 litres medium	CM0401B
Mycoplasma Supplement G		
10 vials	for 80mL medium	SR0059C

Mycoplasma Broth

for the selective isolation of *Mycoplasma* spp.

Mycoplasma Broth Base		
500g	for 19.6 litres medium	CM0403B
Mycoplasma Supplement G		
10 vials	for 80mL medium	SR0059C

Mycoplasma pneumoniae Bi-phasic Medium

for the selective isolation of *Mycoplasma pneumoniae*

Mycoplasma Agar Base		
500g	for 14.1 litres medium	CM0401B
Mycoplasma Supplement P each vial reconstitutes to 20mL and produces 10 quantities of bi-phasic medium		
10 vials	for 10 overlays each	SR0060C

MYP Agar

a medium for the enumeration of *Bacillus cereus* in food samples

MYP Agar Base		
500g	for 11.6 litres medium	CM0929B
Polymyxin B Supplement		
10 vials	for 500mL medium	SR0099E
Egg Yolk Emulsion		
100mL		SR0047C

Neomycin Supplement

allows the growth of clostridia and other anaerobes

10 vials	150mg/vial	SR0163H
----------	------------	---------

★ New York City Media

for the isolation of *Neisseria gonorrhoeae*. Different combinations of growth and selective supplements give media with alternative growth and inhibitory properties; for more information visit www.thermoscientific.com/microbiology.

GC Agar Base		
500g	for 13.9 litres medium	CM0367B
2.5kg	for 69.5 litres medium	CM0367R
5kg	for 139 litres medium	CM0367T
Yeast Autolysate Supplement		
10 vials	for 500mL medium	SR0105B
Liquid Yeast Autolysate		
10 vials	for 2.0 litres medium	SR0182H
Laked Horse Blood		
100mL		SR0048C
LCAT Selective Supplement		
10 vials	for 500mL medium	SR0095B
10 vials	for 2.0 litres medium	SR0095H
VCAT Selective Supplement		
10 vials	for 500mL medium	SR0104B
10 vials	for 2.0 litres medium	SR0104H

Nickerson Medium – see BiGGY Agar

Nitrocefin

for the detection of β-lactamase activity

Nitrocefin with Rehydration Fluid		
2x5 vials	1mg/vial pair	SR0112C

Novobiocin Supplement

can be used to increase the selectivity of some *Salmonella* and *E. coli* O157 media

10 vials	10mg/vial	SR0181E
----------	-----------	---------

Novel Enrichment Broth-Listeria – see ONE Broth-Listeria

Nutrient Agar

500g	for 17.9 litres medium	CM0003B
2.5kg	for 89.5 litres medium	CM0003R
5kg	for 179 litres medium	CM0003T

Nutrient Agar Tablets

each tablet makes 5mL of medium

100 tabs	for 500mL medium	CM0004C
500 tabs	for 2.5 litres medium	CM0004D

Nutrient Broth

500g	for 38.5 litres medium	CM0001B
------	------------------------	---------

Nutrient Broth No. 2

500g	for 20.0 litres medium	CM0067B
------	------------------------	---------

Nutrient Gelatin

for the differentiation of microorganisms by their proteolytic effects

500g	for 3.9 litres medium	CM0635B
------	-----------------------	---------

OGYE Selective Supplement

Oxytetracycline 0.05g/vial		
10 vials	for 500mL medium	SR0073A

OGYE Medium – see Oxytetracycline Glucose Yeast Extract Medium

ONE Broth-Listeria

for selective enrichment of *Listeria* spp. from food and environmental samples in 24 hours. This medium is the selective enrichment broth used in **Listeria Precis**

ONE Broth-Listeria Base		
500g	for 11.4 litres medium	CM1066B
2.5kg	for 57.0 litres medium	CM1066R
5kg	for 114 litres medium	CM1066T

ONE Broth-Listeria Selective Supplement		
10 vials	for 500mL medium	SR0234E
10 vials	for 2 litres medium	SR0234H
10 vials	for 2.25 litres medium	SR0234B

ONE Broth-Salmonella

for selective enrichment of *Salmonella* spp. from food samples in 18 hours when used in **Salmonella Precis**; not recommended for enrichment in other methods

ONE Broth-Salmonella Base		
500g	for 20.0 litres medium	CM1091B
5kg	for 200 litres medium	CM1091T

ONE Broth-Salmonella Selective Supplement		
10 vials	for 225mL medium	SR0242E
10 vials	for 2.25 litres medium	SR0242B

OPSP Agar

for the enumeration of *Clostridium perfringens* in foods

Perfringens Agar Base (OPSP)		
500g	for 11.0 litres medium	CM0543B
2.5kg	for 55.0 litres medium	CM0543R

Perfringens (OPSP) Supplement A		
10 vials	for 500mL medium	SR0076E

Perfringens (OPSP) Supplement B		
10 vials	for 500mL medium	SR0077E

Orange Serum Agar

for the isolation and enumeration of spoilage organisms in citrus products

500g	for 13.5 litres medium	CM0657B
------	------------------------	---------

ORSAB – see Oxacillin Resistance Screening Agar (ORSA)

OSCM

a selective chromogenic medium for the presumptive identification of *Salmonella* spp. from clinical and food samples

Salmonella Chromogenic Agar Base		
500g	for 10.0 litres medium	CM1007B

Salmonella Selective Supplement		
10 vials	for 500mL medium	SR0194E

Oxacillin Resistance Screening Agar (ORSA)

a medium for screening clinical samples for the presence of oxacillin-resistant staphylococci

Oxacillin Resistance Screening Agar Base		
500g	for 4.8 litres medium	CM1008B

ORSAB Selective Supplement		
10 vials	for 500mL medium	SR0195E

Oxford Listeria Agar – see Listeria Selective Agar (Oxford)

Oxytetracycline Glucose Yeast Extract Selective Medium (OGYE)

a medium for the selective enumeration of yeasts and moulds

Oxytetracycline Glucose Yeast Extract Agar Base		
500g	for 13.5 litres medium	CM0545B
2.5kg	for 67.5 litres medium	CM0545R

OGYE Selective Supplement		
10 vials	for 500mL medium	SR0073A

PALCAM Agar

a selective and diagnostic medium for the detection of *Listeria* spp.

PALCAM Agar Base		
500g	for 7.2 litres medium	CM0877B
2.5kg	for 36.0 litres medium	CM0877R

PALCAM Selective Supplement		
10 vials	for 500mL medium	SR0150E
10 vials	for 2.5 litres medium	SR0150B

Peptone Salt Broth – see Maximum Recovery Diluent

Peptone Water (Alkaline) – see Alkaline Peptone Water

Peptone Water (Andrade)

a nutrient base containing Andrade's indicator to which carbohydrates and other diagnostic reagents may be added for use in fermentation studies

500g	for 33.3 litres medium	CM0061B
------	------------------------	---------

Peptone Water

a basal medium to which carbohydrates and an indicator may be added for fermentation studies

500g	for 33.3 litres medium	CM0009B
------	------------------------	---------

PEMBA – see Bacillus cereus Selective Agar

Perfringens (OPSP) Agar – see OPSP Agar

Perfringens TSC Agar – see Tryptose Sulphite Cycloserine Agar

Perfringens SFP Agar – see Shahadi Ferguson Agar

Phosphatidylinositol Supplement – see Chromogenic Listeria Differential Supplement (ISO)

Plate Count Agar

500g	for 28.6 litres medium	CM0325B
2.5kg	for 143.0 litres medium	CM0325R
5kg	for 286 litres medium	CM0325T

Plate Count Agar APHA – see Standard Plate Count Agar

Plate Count Agar with antibiotic free skimmed milk – see Milk Plate Count Agar

Plate Count Agar (Water) – see Water Plate Count Agar

Polymyxin B Supplement

10 vials	50,000 IU/vial	SR0099E
----------	----------------	---------

Potato Dextrose Agar

acid pH medium for the isolation of dermatophytes and other moulds and yeasts. Also, used during the preparation of *Aspergillus niger* for use in the Harmonized Microbial Limits Tests, according to EP/USP/JP.

500g	for 12.8 litres medium	CM0139B
------	------------------------	---------

★ Preston Selective Agar

for the isolation of *Campylobacter jejuni* and *Campylobacter coli* from human, animal, avian and environmental specimens, this medium can be supplemented with *Campylobacter Growth Supplement*

Campylobacter Agar Base		
500g	for 13.5 litres medium	CM0689B
Preston Campylobacter Selective Supplement		
10 vials	for 500mL medium	SR0117E
OR		
Modified Preston Campylobacter Selective Supplement alternative supplement containing amphotericin B as a replacement for cycloheximide		
10 vials	for 500mL medium	SR0204E
Campylobacter Growth Supplement		
10 vials	for 500mL medium	SR0232E

Preston Selective Broth

for the recovery of *Campylobacter jejuni* and *Campylobacter coli*, this medium can be supplemented with *Campylobacter Growth Supplement*

Nutrient Broth No. 2		
500g	for 20.0 litres medium	CM0067B
Preston Campylobacter Selective Supplement		
10 vials	for 500mL medium	SR0117E
OR		
Modified Preston Campylobacter Selective Supplement alternative supplement containing amphotericin B as a replacement for cycloheximide		
10 vials	for 500mL medium	SR0204E
Campylobacter Growth Supplement		
10 vials	for 500mL medium	SR0232E

Pseudomonas CFC Agar

a selective medium for the isolation of *Pseudomonas* spp.

Pseudomonas Agar Base		
500g	for 10.3 litres medium	CM0559B
2.5kg	for 51.5 litres medium	CM0559R
Pseudomonas CFC Supplement		
10 vials	for 500mL medium	SR0103E

Pseudomonas CN Agar

a selective medium for the isolation of *Pseudomonas aeruginosa*

Pseudomonas Agar Base		
500g	for 10.3 litres medium	CM0559B
2.5kg	for 51.5 litres medium	CM0559R
Pseudomonas CN Supplement		
10 vials	for 500mL medium	SR0102E

Pseudomonas cepacia Medium – see *Burkholderia cepacia* Medium

Pseudomonas Cetrinide Agar Base

a medium for the selective isolation and identification of *Pseudomonas aeruginosa* from a range of samples. This formulation is recommended in USP, EP, and AOAC methods.

500g	for 11.0 litres medium	CM0579B
------	------------------------	---------

R2A Agar

a medium for the enumeration of heterotrophic bacteria from water samples

500g	for 27.6 litres medium	CM0906B
------	------------------------	---------

Raka-Ray Medium

a medium for the isolation of lactic acid bacteria in beer and brewing processes. Requires supplementation with phenylethanol, sorbitan mono-oleate and cycloheximide.

Raka-Ray Agar Base		
500g	for 6.5 litres medium	CM0777B

Rappaport-Vassiliadis (RV) Enrichment Broth

a selective enrichment broth for the isolation of salmonellae. **N.B.** maximum shelf life of 12 months

500g	for 16.7 litres medium	CM0669B
------	------------------------	---------

Rappaport Vassiliadis Soya Enrichment Broth MLT (EP/USP/JP/BP) - see RVS Enrichment Broth MLT (EP/USP/JP/BP)

Rappaport-Vassiliadis Soya (RVS) Peptone Broth

a buffered, selective enrichment broth for the isolation of salmonellae. **N.B.** maximum shelf life of 12 months

500g	for 18.7 litres medium	CM0866B
2.5kg	for 93.5 litres medium	CM0866R

Reinforced Clostridial Agar

for the cultivation and enumeration of clostridia and other anaerobes

500g	for 9.5 litres medium	CM0151B
2.5kg	for 47.5 litres medium	CM0151R

Reinforced Clostridial Medium MLT (EP/USP/JP/BP)

a medium intended for use in qualitative procedures for cultivation of clostridia and other anaerobic organisms

500g	for 13.1 litres medium	R455402
------	------------------------	---------

Reinforced Clostridial Medium (RCM)

use this semisolid medium for the enumeration of clostridia and other anaerobes

500g	for 13.2 litres medium	CM0149B
5kg	for 132 litres medium	CM0149T

Rogosa Agar

for the selective isolation and enumeration of lactobacilli

500g	for 6.1 litres medium	CM0627B
------	-----------------------	---------

Rothe Broth – see Azide Dextrose Broth

Rose-Bengal Chloramphenicol Selective Medium

a medium for the selective enumeration of moulds and yeasts from food

Rose-Bengal Chloramphenicol Agar Base		
500g	for 15.6 litres medium	CM0549B
2.5kg	for 78.0 litres medium	CM0549R
5kg	for 156 litres medium	CM0549T
Chloramphenicol Supplement		
10 vials	for 500mL medium	SR0078E
10 vials	for 2.0 litres medium	SR0078H

RPF – see Baird Parker RPF

RV – see Rappaport-Vassiliadis Enrichment Broth

RVS – see Rappaport-Vassiliadis Soya Peptone Broth

★ requires the addition of blood
Not all products available in all territories. Please contact your local office for details.

RVS Enrichment Broth MLT (EP/USP/JP/BP)

a liquid medium recommended for use in qualitative procedures for selective enrichment and isolation of *Salmonella* spp.

500g	for 18.4 litres medium	R455432
------	------------------------	---------

Ryan's Agar – see *Aeromonas* Medium

Sabouraud Dextrose Agar

acid pH medium for the isolation of dermatophytes and other moulds and yeasts. This can be made selective using chloramphenicol.

500g	for 7.7 litres medium	CM0041B
2.5kg	for 38.5 litres medium	CM0041R
5kg	for 77.0 litres medium	CM0041T

Sabouraud Dextrose Liquid Medium (Sabouraud Dextrose Broth) EP/USP/JP

a liquid medium recommended for testing sterile and nonsterile products and for the determination of fungistatic activity in pharmaceutical products

500g	for 16.7 litres medium	CM0147B
------	------------------------	---------

Sabouraud Maltose Agar

for use where maltose is the preferred carbohydrate source, it may be made selective by the addition of chloramphenicol

500g	for 7.7 litres medium	CM0541B
------	-----------------------	---------

Salmonella Chromogenic Medium - see OSCM

Salmonella Chromogenic Medium II - see *Brilliance*

Salmonella Agar

Salmonella Novel Enrichment Broth – see ONE

Broth-*Salmonella*

Salmonella Precis Method

a rapid culture method for the enrichment, isolation and identification of *Salmonella* spp. in just 2 days. Validated by AFNOR to ISO 16140 for all food, animal feed and environmental samples; use *Salmonella* Latex (DR1108A) for rapid confirmation.

ONE Broth-Salmonella

an 18 hour selective enrichment broth for *Salmonella* spp., only for use with the *Salmonella* Precis Method

ONE Broth-Salmonella Base		
500g	for 20.0 litres medium	CM1091B
5kg	for 200 litres medium	CM1091T
ONE Broth-Salmonella Selective Supplement		
10 vials	for 225mL medium	SR0242E
10 vials	for 2.25 litres medium	SR0242B
ONE Broth-Salmonella Complete		
10x20L Dry-Bag	for 20 litres of medium	DB1091W
10x20L Dry-Bag with filter	for 20 litres of medium	DB1091M

Brilliance Salmonella Agar

a novel chromogenic medium incorporating Inhibigen technology, for the presumptive identification of *Salmonella* spp. from food and clinical samples

Brilliance Salmonella Agar Base		
500g	for 9.3 litres medium	CM1092B
5kg	for 93 litres medium	CM1092T
Salmonella Selective Supplement		
10 vials	for 500mL medium	SR0194E

Salmonella/*Shigella* Agar – see SS Agar

Salmonella/*Shigella* Modified – see SS Agar Modified

Salt Meat Broth Tablets

an enrichment medium for the isolation of staphylococci from grossly contaminated specimens

100 tabs	each tab for 5mL	CM0094C
----------	------------------	---------

Schaedler Anaerobe Agar

a thioglycollate-free medium for the growth of aerobic and anaerobic organisms

500g	for 12.5 litres medium	CM0437B
------	------------------------	---------

Schaedler Anaerobe Broth

a thioglycollate-free broth for the general growth of anaerobes, for use in blood culture and antibiotic MIC studies

500g	for 18.9 litres medium	CM0497B
------	------------------------	---------

Schiemann's CIN Agar – see *Yersinia* Selective (CIN) Agar

Selenite Broth (10mL)

a ready-to-use broth for the isolation of salmonellae from faeces and food products

24x10mL		LR0039A
---------	--	---------

Selenite Broth with Mannitol – see Mannitol Selenite Broth

Selenite Cystine Broth

a selective enrichment medium for the isolation of specimens undergoing examination for salmonellae

Selenite Cystine Broth Base		
500g	for 26.3 litres medium	CM0699B
Sodium Biselenite		
100g	for 25.0 litres medium	LP0121A

Selenite F Broth

for the selective enrichment of salmonellae from faeces and food products

Selenite Broth Base (Lactose)		
500g	for 26.3 litres medium	CM0395B
Sodium Biselenite		
100g	for 25.0 litres medium	LP0121A

Sensitest Agar

for antimicrobial susceptibility testing

500g	for 15.6 litres medium	CM0409B
2.5kg	for 78.0 litres medium	CM0409R

Shahadi Ferguson Agar

a selective medium for the presumptive identification and enumeration of *Clostridium perfringens*

Perfringens Agar Base (TSC/SFP)		
500g	for 10.9 litres medium	CM0587B
Perfringens (SFP) Supplement		
10 vials	for 500mL medium	SR0093E
Egg Yolk Emulsion		
100mL		SR0047C

* Sheep Blood Agar Base

formulated to give improved haemolytic reactions with sheep blood

500g	for 12.5 litres medium	CM0854B
------	------------------------	---------

SIM Medium

a motility-indole medium for the differentiation of the Enterobacteriaceae

500g	for 16.7 litres medium	CM0435B
------	------------------------	---------

Simmons Citrate Agar

for the differentiation of the Enterobacteriaceae

500g	for 21.7 litres medium	CM0155B
------	------------------------	---------

* Skirrow Selective Medium

a selective medium for the isolation of *Campylobacter* spp. at 42°C, this medium can be used with the Campylobacter Growth Supplement

Campylobacter Selective Supplement (Skirrow)

10 vials	for 500mL medium	SR0069E
10 vials	for 2.0 litres medium	SR0069H

Blood Agar Base No. 2

500g	for 12.5 litres medium	CM0271B
------	------------------------	---------

OR

Columbia Blood Agar Base

500g	for 12.8 litres medium	CM0331B
------	------------------------	---------

Campylobacter Growth Supplement

10 vials	for 500mL medium	SR0232E
----------	------------------	---------

Slanetz and Bartley Medium (Enterococcus Agar)

for the detection and enumeration of faecal enterococci

500g	for 11.9 litres medium	CM0377B
2.5kg	for 59.5 litres medium	CM0377R

Sorbitol MacConkey Agar

a selective differential medium for detection of *E. coli* O157. This medium can be made more selective by the addition of C-T Supplement – see Cefixime Tellurite Sorbitol MacConkey Agar

500g	for 9.7 litres medium	CM0813B
2.5kg	for 48.5 litres medium	CM0813R

Sorbitol MacConkey Agar with BCIG

selective and differential medium for the detection of *Escherichia coli* O157, incorporating BCIG which can be used with C-T Supplement

500g	for 9.7 litres medium	CM0981B
------	-----------------------	---------

C-T Supplement store below 0°C

10 vials	for 500mL medium	SR0172E
10 vials	for 2.0 litres medium	SR0172H

Sorbitol MacConkey Agar with Cefixime and Rhamnose – see Cefixime Rhamnose Sorbitol MacConkey Agar

Soybean Casein Digest Agar – see Tryptone Soya Agar

Soybean Casein Digest Medium – see Tryptone Soya Broth

SS Agar

for the isolation of salmonellae and shigellae

500g	for 7.9 litres medium	CM0099B
------	-----------------------	---------

SS Agar Modified

for the isolation of salmonellae and shigellae

500g	for 8.8 litres medium	CM0533B
------	-----------------------	---------

STAA Agar

medium for the isolation of *Brochothrix thermosphacta* from food samples

STAA Agar Base

500g	for 13.5 litres medium	CM0881B
------	------------------------	---------

STAA Selective Supplement with actidione

10 vials	for 500mL medium	SR0151E
----------	------------------	---------

STA Agar

medium for the isolation of *Brochothrix thermosphacta* from food samples

STAA Agar Base

500g	for 13.5 litres medium	CM0881B
------	------------------------	---------

STA Selective Supplement without actidione

10 vials	for 500mL medium	SR0162E
----------	------------------	---------

Standard Plate Count Agar APHA

a standard medium complying with the APHA formulation for use in testing milk, water, food and dairy products

500g	for 21.3 litres medium	CM0463B
5kg	for 213 litres medium	CM0463T

Staphylococcus Medium No. 110

a selective medium for the isolation and differentiation of pathogenic staphylococci

500g	for 3.3 litres medium	CM0145B
------	-----------------------	---------

Staph/Strep Selective Medium - see CNA Staph/Strep Selective Agar

Streptococcus Selective Agar - see COBA Selective Medium

Stuart Transport (Modified) Medium

a semisolid, non-nutritional transport medium

500g	for 31.3 litres medium	CM0111B
------	------------------------	---------

TBX Agar

Tryptone Bile Agar with X-Glucuronide, a chromogenic medium for the detection of *E. coli*

500g	for 13.7 litres medium	CM0945B
2.5kg	for 68.5 litres medium	CM0945R
5kg	for 137 litres medium	CM0945T

TCBS Cholera Medium – see Cholera Medium TCBS

Tergitol-7 Agar

for the enumeration and identification of coliform bacteria

Tergitol-7 Agar Base		
500g	for 9.2 litres medium	CM0793B
TTC Solution (0.125%)		
10x2mL	for 100mL medium	SR0148A
TTC Solution (5%)		
100mL		SR0211C

Tetrathionate Broth Base

for the selective enrichment of salmonellae

500g	for 6.5 litres medium	CM0029B
------	-----------------------	---------

Tetrathionate Broth Base (USA)

an American formulation for the selective enrichment of specimens undergoing examination for salmonellae

500g	for 10.9 litres medium	CM0671B
------	------------------------	---------

Tetrathionate (Muller-Kauffmann) – see Muller-Kauffmann Tetrathionate Broth Base

★ Thayer Martin Media

for the isolation of *Neisseria gonorrhoeae*. Different combinations of growth and selective supplements give media with alternative growth and inhibitory properties; for more information visit www.thermoscientific.com/microbiology.

GC Agar Base		
500g	for 13.9 litres medium	CM0367B
2.5kg	for 69.5 litres medium	CM0367R
5 kg	for 139 litres medium	CM0367T
Haemoglobin Powder Soluble		
500g	for 50.0 litres medium	LP0053B
Vitox plus diluent		
2x5 vials	for 500mL medium	SR0090A
2x5 vials	for 2.0 litres medium	SR0090H
VCN Selective Supplement		
10 vials	for 500mL medium	SR0101E
VCNT Selective Supplement		
10 vials	for 500mL medium	SR0091E
10 vials	for 2.0 litres medium	SR0091H

Thioglycollate Broth, alternative

a non-viscous thioglycollate medium for the culture of turbid biological products

500g	for 17.2 litres medium	CM0391B
------	------------------------	---------

Thioglycollate Medium

500g	for 16.9 litres medium	CM0173B
5kg	for 169.0 litres medium	CM0173T

Thioglycollate Medium (Brewer)

500g	for 25.0 litres medium	CM0023B
------	------------------------	---------

Tinsdale Medium

a medium for the isolation and identification of *Corynebacterium diphtheria*

Tinsdale Agar Base		
500g	for 11.1 litres medium	CM0487B
Tinsdale Supplement		
10 vials	for 200mL medium	SR0065A

Todd-Hewitt Broth

for the cultivation of streptococci prior to serological grouping

500g	for 13.7 litres medium	CM0189B
------	------------------------	---------

Tomato Juice Agar

for the cultivation of lactobacilli

500g	for 9.6 litres medium	CM0113B
------	-----------------------	---------

Tributyrin Agar

100mL		PM0004C
-------	--	---------

Trichomonas Medium

a medium for the cultivation of *Trichomonas vaginalis*

Trichomonas Medium Base		
500g	for 13.3 litres medium	CM0161B
Horse Serum		
100mL		SR0035C

Trichomonas Medium No. 2

a ready-to-use medium for the cultivation of *Trichomonas vaginalis*

24x5mL		LR0027A
--------	--	---------

Triple Sugar Iron Agar

a composite medium for the differentiation of the Enterobacteriaceae

500g	for 7.7 litres medium	CM0277B
------	-----------------------	---------

Tryptone Bile Agar

for the detection and enumeration of *E. coli* biotype 1 in foods

500g	for 13.7 litres medium	CM0595B
------	------------------------	---------

Tryptone Bile X-Glucuronide Medium – see TBX Agar

Tryptone Glucose Extract Agar

a USA standard methods medium for the examination of dairy products

500g	for 20.8 litres medium	CM0127B
------	------------------------	---------

Tryptone Soya Agar (Soybean Casein Digest Agar)

a general purpose medium for the growth of a wide variety of organisms; complies with formula specified in EP/USP/JP/BP

500g	for 12.5 litres medium	CM0131B
2.5kg	for 62.5 litres medium	CM0131R
5kg	for 125.0 litres medium	CM0131T

Tryptone Soya Broth (Soybean Casein Digest Medium)

recommended for general laboratory use and sterility testing; complies with formula specified in EP/USP/JP/BP

500g	for 16.7 litres medium	CM0129B
2.5kg	for 83.3 litres medium	CM0129R
5kg	for 167.0 litres medium	CM0129T

Tryptone Soya Broth (Cold Filterable) – see Cold Filterable Tryptone Soya Broth

Tryptone Soya Broth Modified (mTSB)

a broth for *E. coli* O157, the selectivity of which can be increased by using supplements as appropriate

Modified Tryptone Soya Broth (mTSB) Base		
500g	for 15.2 litres medium	CM0989B
Novobiocin Supplement		
10 vials	for 500mL medium	SR0181E

TSC Agar – see Tryptose Sulphite Cycloserine Agar

Tryptone Water

for the detection of indole production

500g	for 33.3 litres medium	CM0087B
------	------------------------	---------

★ Tryptose Blood Agar Base

a nutritious medium specially formulated for use as a blood agar for the growth of fastidious organisms

500g	for 16.7 litres medium	CM0233B
------	------------------------	---------

Tryptose Phosphate Broth

for the cultivation of fastidious bacteria and as an adjuvant to tissue culture media

500g	for 16.9 litres medium	CM0283B
5kg	for 169.0 litres medium	CM0283T

Tryptose Sulphite Cycloserine (TSC) Agar

a selective medium for the presumptive identification and enumeration of *Clostridium perfringens*

Perfringens Agar Base (TSC/SFP)		
500g	for 10.9 litres medium	CM0587B
Perfringens (TSC) Supplement		
10 vials	for 500mL medium	SR0088E
Egg Yolk Emulsion		
100mL		SR0047C

Universal Beer Agar

for the isolation of beer spoilage organisms; can be supplemented with cycloheximide

Universal Beer Agar		
500g	for 8.1 litres medium	CM0651B
Cycloheximide 0.1% Solution		
100mL		SR0222C

Urea Agar

for the detection of urease production

Urea Agar Base (Christensen Agar Base)		
500g	for 20.8 litres medium	CM0053B
Urea 40% Solution		
10x5mL	for 100mL medium	SR0020K

Urea Broth

for the detection of urease production

Urea Broth Base (Christensen Broth Base)		
500g	for 55.6 litres medium	CM0071B
Urea 40% Solution		
10x5mL	for 100mL medium	SR0020K

UVM Broths – see Listeria Selective Enrichment Media (UVM Formulations)

Vancomycin Supplement

10 vials	3mg/vial	SR0186E
10 vials	5mg/vial	SR0247E

Violet Red Bile (Lactose) Agar

for the presumptive enumeration of coliforms in foods

500g	for 13.0 litres medium	CM0107B
2.5kg	for 65.0 litres medium	CM0107R
5kg	for 130.0 litres medium	CM0107T

Violet Red Bile (Lactose) Agar (ISO)

for the presumptive enumeration of coliforms in foods; meeting the ISO recovery specification

500g	for 13.0 litres medium	CM0968B
------	------------------------	---------

Violet Red Bile Agar (VRBA) with MUG

allows presumptive identification of *E. coli*

500g	for 13.0 litres medium	CM0978B
------	------------------------	---------

Violet Red Bile Glucose Agar

for the presumptive enumeration of Enterobacteriaceae

500g	for 13.0 litres medium	CM0485B
2.5kg	for 65.0 litres medium	CM0485R
5kg	for 130.0 litres medium	CM0485T

Violet Red Bile Glucose Agar (EP/USP/JP/BP)

a solid medium recommended for use in quantitative and qualitative procedures for selective and differential isolation of bile-tolerant, Gram-negative bacteria

500g	for 11.7 litres medium	R455302
------	------------------------	---------

Violet Red Bile Glucose Agar (ISO)

for the presumptive enumeration of Enterobacteriaceae in foods; meeting the ISO recovery specification

500g	for 13.0 litres medium	CM1082B
2.5kg	for 65 litres medium	CM1082R

Vogel-Johnson Agar

a USP formulation for the detection of *Staphylococcus aureus*

500g	for 8.2 litres medium	CM0641B
Potassium Tellurite 3.5%		
10x2mL		SR0030J

VRBA – see Violet Red Bile (Lactose) Agar

VRBA (ISO) – see Violet Red Bile (Lactose) Agar (ISO)

VRBGA – see Violet Red Bile Glucose Agar

VRBGA (ISO) – see Violet Red Bile Glucose Agar (ISO)

VRE Media

selective media for the isolation of vancomycin-resistant enterococci (VRE) and high-level aminoglycoside-resistant enterococci (HLARE) from clinical samples; for more information visit www.thermoscientific.com/microbiology.

VRE Agar Base		
500g	for 11.7 litres medium	CM0985B
Meropenem Supplement		
10 vials	for 500mL medium	SR0184E
Gentamicin Supplement		
10 vials	for 500mL medium	SR0185E
Vancomycin Supplement		
10 vials	for 500mL medium	SR0186E

Brilliance VRE Agar – see page Chromogenic Media

Wadowsky Yee Medium – see MWY Medium

Water Plate Count Agar (ISO)

a medium for the enumeration of culturable microorganisms from water

500g	for 20.8 litres medium	CM1012B
------	------------------------	---------

* Wilkins-Chalgren Anaerobe Agar

for the growth of anaerobes, and recommended for antimicrobial susceptibility testing of anaerobes; can be made more selective with the addition of Neomycin Supplement, N-S Anaerobe Supplement or G-N Anaerobe Supplement

Wilkins-Chalgren Anaerobe Agar		
500g	for 11.6 litres medium	CM0619B
Neomycin Supplement 150mg/vial		
10 vials	for 2.0 litres medium	SR0163H
N-S Anaerobe Selective Supplement for the selective isolation of nonsporing anaerobes		
10 vials	for 500mL medium	SR0107B

Wilkins-Chalgren Anaerobe Broth

for the growth of anaerobes, and recommended for antimicrobial susceptibility testing of anaerobes

500g	for 15.2 litres medium	CM0643B
------	------------------------	---------

Wilson and Blair – see Bismuth Sulphite Agar

WL Nutrient Agar

for the determination of the microbiological flora in brewing and fermentation processes; can be supplemented with cycloheximide

500g	for 6.7 litres medium	CM0309B
Cycloheximide 0.1% Solution		
100mL		SR0222C

WL Nutrient Broth

for the determination of the microbiological flora in brewing and fermentation processes; can be supplemented with cycloheximide

500g	for 8.3 litres medium	CM0501B
Cycloheximide 0.1% Solution		
100mL		SR0222C

Wort Agar

for the cultivation and enumeration of yeasts

500g	for 10.0 litres medium	CM0247B
------	------------------------	---------

XLD Agar (EP/USP/JP/BP)

a solid medium recommended for use in qualitative procedures for selective and differential isolation of enteric pathogens

500g	for 9.0 litres medium	R459902
------	-----------------------	---------

XLD Medium

a selective medium for the isolation of salmonellae and shigellae

500g	for 9.4 litres medium	CM0469B
2.5kg	for 47.0 litres medium	CM0469R
5kg	for 94 litres medium	CM0469T

XLT-4 Agar

a highly selective medium for the isolation of salmonellae and shigellae from clinical, environmental and food samples, according to Miller and Tate

XLT-4 Agar Base		
500g	for 8.5 litres medium	CM1061B
XLT-4 Selective Supplement		
100mL	for 21.7 litres medium	SR0237C

Yeast Autolysate – see Liquid Yeast Autolysate

Yeast and Mould Agar

for the isolation and maintenance of yeasts and moulds

500g	for 12.2 litres medium	CM0920B
------	------------------------	---------

Yeast Extract Agar

a plate count agar for microorganisms in water

500g	for 21.7 litres medium	CM0019B
2.5kg	for 108.5 litres medium	CM0019R
5kg	for 217 litres medium	CM0019T

Yersinia Selective Medium (CIN) Agar

for the isolation and enumeration of *Yersinia enterocolitica* from clinical specimens and foods

Yersinia Selective (CIN) Agar Base		
500g	for 8.6 litres medium	CM0653B
Yersinia Selective Supplement		
10 vials	for 500mL medium	SR0109E

For a wide range of bagged media formats in a variety of sizes, please see pages 54 - 55.

Pre-Weigh Dehydrated Culture Media

Pre-Weigh Buffered Peptone Water (BPW)		
Pack of 20	20 x 1L pouches	PW0509-1L
Pre-Weigh Buffered Peptone Water (BPW)		
Pack of 10	10 x 8L pouches	PW0509-8L
Pre-Weigh Buffered Peptone Water (BPW) (ISO)		
Pack of 20	20 x 1L pouches	PW1049-1L
Pre-Weigh Buffered Peptone Water (BPW) (ISO)		
Pack of 10	10 x 8L pouches	PW1049-8L
Pre-Weigh Violet Red Bile Glucose Agar (VRBGA)		
Pack of 20	20 x 1L pouches	PW0485-1L
Pre-Weigh Violet Red Bile Glucose Agar (VRBGA)		
Pack of 10	10 x 8L pouches	PW0485-8L
Pre-Weigh Plate count Agar (PCA)		
Pack of 20	20 x 1L pouches	PW0325-1L
Pre-Weigh Plate count Agar (PCA)		
Pack of 10	10 x 8L pouches	PW0325-8L
Pre-Weigh Mannitol Egg Yolk Polymyxin Agar (MYP)		
Pack of 20	20 x 1L pouches	PW0929-1L
Pre-Weigh Mannitol Egg Yolk Polymyxin Agar (MYP)		
Pack of 10	10 x 8L pouches	PW0929-8L

Supplementary Reagents

Supplementary Reagent products should be stored at 2–8°C unless otherwise stated on the packaging.

Broad Spectrum Beta-Lactamase Mixture

a mixture of beta-lactamase enzymes from *Bacillus cereus*, for the inactivation of beta-lactam antibiotic preparations present in specimens

10 vials		SR0113E
----------	--	---------

Egg Yolk Emulsion

a stabilized emulsion of egg yolk for use in culture media

100mL		SR0047C
-------	--	---------

Egg Yolk Tellurite Emulsion

an emulsion of egg yolk containing potassium tellurite for use in Baird-Parker Egg Yolk Media

100mL		SR0054C
-------	--	---------

Fildes Extract

a peptic digest of blood which is a rich source of growth factors, recommended for the preparation of many culture media

100mL		SR0046C
-------	--	---------

Horse Serum

an extract of horse blood for addition to culture media

100mL		SR0035C
-------	--	---------

Lactic Acid (10%)

10x1mL		SR0021H
10x5mL		SR0021K

Laked Horse Blood

haemolysed blood for addition to culture media

100mL		SR0048C
-------	--	---------

Nitrocefin

for the rapid, chromogenic detection of beta-lactamase activity; contains nitrocefin, with rehydrating fluid

2x5 vials		SR0112C
-----------	--	---------

Penase

a lactamase enzyme from *Bacillus cereus* for the inactivation of a range of penicillins. Store at –20°C.

1 vial		SR0129B
--------	--	---------

Potassium Lactate

10x1mL		SR0037H
--------	--	---------

Potassium Tellurite 3.5%

10x2mL		SR0030J
--------	--	---------

RPF Supplement

only to be used with Baird-Parker Agar Base (RPF) (CM0961). Store below 0°C.

10 vials		SR0122A
----------	--	---------

Sputasol

a liquifying agent for sputum samples

10 vials		SR0233A
----------	--	---------

TTC Solution (0.125%)

for addition to Tergitol-7 Agar

10x2mL		SR0148A
--------	--	---------

TTC Solution (1%)

for addition to K-F Streptococcus Agar

10x5mL		SR0229K
--------	--	---------

TTC Solution (5%)

for addition to Tergitol-7 Agar or K-F Streptococcus Agar

100mL		SR0211C
-------	--	---------

Urea 40% Solution

a solution for addition to culture media

10x5mL		SR0020K
--------	--	---------

Antibiotic Single Supplements

Antibiotic Single Supplements should be stored at 2–8°C unless otherwise stated on the packaging.

Ampicillin Supplement

10 vials	2.5mg/vial	SR0136E
----------	------------	---------

Cefixime Supplement

10 vials	0.025mg/vial	SR0191E
----------	--------------	---------

Cephalexin (Bordetella Selective) Supplement

10 vials	20mg/vial	SR0082E
----------	-----------	---------

Chloramphenicol Supplement

10 vials	50mg/vial	SR0078E
10 vials	200mg/vial	SR0078H

Cycloheximide Solution

100mL	0.1% solution	SR0222C
-------	---------------	---------

D-Cycloserine (Perfringens TSC) Supplement

10 vials	200mg/vial	SR0088E
----------	------------	---------

Gentamicin Supplement

10 vials	256mg/vial	SR0185E
----------	------------	---------

Kanamycin Sulphate Supplement

10 vials	10mg/vial	SR0092E
----------	-----------	---------

Meropenem Supplement

10 vials	1mg/vial	SR0184E
----------	----------	---------

Neomycin Supplement

10 vials	150mg/vial	SR0163H
----------	------------	---------

Nitrocefin with Rehydration Fluid

2x5 vials	1mg/vial pair	SR0112C
-----------	---------------	---------

Novobiocin Supplement

10 vials	10mg/vial	SR0181E
----------	-----------	---------

Oxytetracycline (OGYE Selective) Supplement

10 vials	0.05g/vial	SR0073A
----------	------------	---------

Polymyxin B Supplement

10 vials	50,000 IU/vial	SR0099E
----------	----------------	---------

Vancomycin Supplement

10 vials	3mg/vial	SR0186E
10 vials	5mg/vial	SR0247E

Blood Products

Horse Blood Oxalated

100mL		SR0049C
-------	--	---------

Horse Blood, Defibrinated *

25mL		SR0050B
100mL		SR0050C

Laked Horse Blood

sterile haemolysed blood for addition to culture media; recommended for *Corynebacterium* media

100mL		SR0048C
-------	--	---------

Sheep Blood in Alsevers Solution

25mL		SR0053B
------	--	---------

Sheep Blood Defibrinated

25mL		SR0051B
------	--	---------

Thermo Scientific™ Ready Poured Plates

Blood Agars

Blood Agar Base No. 2 with Horse Blood	
10x90mm	PB0114A
Blood Agar Base No. 2 with Sheep Blood	
10x90mm	PB0115A
Columbia Agar with Horse Blood	
10x90mm	PB0122A
Columbia Agar with Sheep Blood	
10x90mm	PB5008A
Columbia Agar with Sheep Blood PLUS	
10x90mm	PB0123A*^
10x90mm	PB5039A
Columbia Agar with Chocolate Horse Blood	
10x90mm	PB0124A

General Anaerobe Agars

Thermo Scientific™ A.R.I.A.™ with Horse Blood		
10x90mm	PB1243A	
A.R.I.A. with Horse Blood+Neomycin		
10x90mm	PB1244A	
Actinomyces Selective Agar		
10x90mm	PB1220A	
Anaerobe Blood Agar (Wilkins Chalgren) with Neomycin		
10x90mm	PB0112A	
Anaerobe Blood Agar (Wilkins Chalgren) with Nalidixic Acid and Tween™		
10x90mm	PB0113A	
Anaerobe Blood Agar with Neomycin and Aztreonam		
10x90mm	PB1092A	
Columbia Blood Agar with Neomycin		
10x90mm	PB0219A	
Schaedler Anaerobe Agar with Sheep Blood with Haemin and Vitamin K1		
10x90mm	PB5034A	
Schaedler Anaerobe KV Selective Agar with Lysed Horse Blood		
10x90mm	PO5020A	
Yeast Extract Cysteine Medium with Sheep Blood (Beerens Formulation)		
10x90mm	PB5101A	
Schaedler Anaerobe Agar / Schaedler Anaerobe KV Selective Agar		
10x90mm	Biplate	PB5204E

Susceptibility Test Agars

Columbia Agar Base with 2% Salt (4mm depth)	
10x90mm	PO0879A
Diagnostic Sensitivity Test Agar	
10x90mm	PO0130A^
10x90mm	PO5068A+
Diagnostic Sensitivity Test (D.S.T.) with chocolate Horse Blood	
10x90mm	PB0133A
Haemophilus Test Medium	
10x90mm	PO0228A^
10x90mm	PO5138A+
Iso-Sensitest Agar	
10x90mm	PO0144A^
Iso-Sensitest Agar (Deep Fill)	
10x90mm	PO0779A

Iso-Sensitest Agar with Chocolate Horse Blood	
10x90mm	PB0147A
Iso-Sensitest Agar with Horse Blood	
10x90mm	PB0146A
Iso-Sensitest Agar with Horse Blood and 20mg/L NAD	
10x90mm	PB0378A
Iso-Sensitest Agar with Lysed Horse Blood	
10x90mm	PB0145A
Iso-Sensitest Agar with Sheep Blood	
10x90mm	PB5003A
Mueller-Hinton Agar	
10x90mm	PO5007A+
10x90mm	PO0152A^
10x140mm	PO1191S*
Mueller-Hinton Agar with 2% Salt	
10x90mm	PO0799A^
10x90mm	PO5139A+
Mueller-Hinton Agar with 5% Sheep Blood	
10x90mm	PB0431A^
10x90mm	PB5007A+
Mueller-Hinton Agar with Horse Blood	
10x90mm	PB1229A^
10x90mm	PB5303A+

Thermo Scientific™ Brilliance™ and other Chromogenic Media

Brilliance Candida Agar		
10x90mm	PO1034A^	
10x90mm	PO5170A+	
Brilliance Candida Agar / Sabouraud Dextrose Agar		
10x90mm	Biplate	PO5258E
Brilliance CampyCount Agar		
10x90mm	PO1185A^	
Brilliance CRE Agar		
10x90mm	PO1226A	
Brilliance E. coli/coliform Agar		
10x90mm	PO5176A*	
Brilliance ESBL Agar		
10x90mm	PO5302A	
Brilliance ESBL/Brilliance CRE		
10x90mm	Biplate	PO1265E
Brilliance GBS Agar		
10x90mm	PO5320A	
Brilliance GBS Agar/ Columbia CNA Agar		
10x90mm	Biplate	PB5260E
Brilliance Listeria Agar		
10x90mm	PO1102A^	
10x90mm	PO5165A+	
Brilliance MRSA 2 Agar		
10x90mm	PO1210A^	
10x90mm	PO5310A+	
Brilliance MRSA 2/Brilliance Staph24 Agar		
10x90mm	Biplate	PO1258E
Brilliance MRSA 2 Agar/Columbia Agar with Sheep Blood PLUS		
10x90mm	Biplate	PB5253E
Brilliance Salmonella Agar		
10x90mm	PO5098A	

* Limited supply. Contact your local representative. + Mainland Europe ^ Only available in UK
Not all products available in all territories. Please contact your local office for details.

Brilliance Staph 24 Agar		
10x90mm		PO1186A^
Brilliance UTI Agar		
10x90mm		PO0794A^
10x90mm		PO5120A+
Brilliance UTI Clarity Agar		
10x90mm		PO1110A^
10x90mm		PO5159A+
Brilliance VRE Agar		
10x90mm		PO1175A
MLGA		
10x90mm		PO1016A
Oxoid Chromogenic Listeria Agar (ISO) – OCLA (ISO)		
10x90mm		PO5183A
10x90mm		PO1196A
Chromogenic Coliform Agar		
10x90mm		PO5318A
10x55mm		PO5428J
Oxoid Salmonella Chromogenic Medium (OSCM)		
10x90mm		PO0958A
TBX Agar		
10x90mm		PO0727A^
10x90mm		PO5109A+

Biplates

Brilliance Candida Agar/Sabouraud GC Selective Agar		
10x90mm	Biplate	PO5258E
Brilliance ESBL/Brilliance CRE Agar		
10x90mm	Biplate	PO1265E
Brilliance GBS/ Columbia CNA Agar		
10x90mm	Biplate	PB5260E
Brilliance MRSA 2 Agar/Brilliance Staph 24		
10x90mm	Biplate	PO1258E
Brilliance MRSA 2 Agar/Columbia-Agar with sheep blood		
10x90mm	Biplate	PB5253E
Karmali Selective Medium/Karmali Selective Medium		
10x90mm	Biplate	PO5219E
Karmali Selective Medium/Karmali Selective Medium MacConkey-Nährboden Nr. 3		
10x90mm	Biplate	PO5217E
C.L.E.D. Medium/MacConkey Agar No. 3		
10x90mm	Biplate	PB5200E
Columbia Agar with Blood/MacConkey Agar No. 3		
10x90mm	Biplate	PB5207E
Columbia Agar with Blood/Chocolate Agar		
10x90mm	Biplate	PB5202E
Columbia CNA Aesculin Selective Agar/Brilliance UTI		
10x90mm	Biplate	PB5220E
Columbia CNA Aesculin Selective Agar/MacConkey Agar No. 3, mod.		
10x90mm	Biplate	PB5224E
Hektoen-Enteric-Agar/DCA-Agar (Leifson)		
10x90mm	Biplate	PO5257E
Salmonella-Shigella-Agar/XLD-Agar		
10x90mm	Biplate	PO5210E

Schaedler Anaerobe Agar/Schaedler Anaerobe KV Selective Agar		
10x90mm	Biplate	PB5204E
XLD-Agar/Brilliance Salmonella Agar		
10x90mm	Biplate	PO5248E

General Purpose Media

Columbia Agar Base		
10x90mm		PO0537A*
Columbia Agar with Sheep Blood		
(2x10) 55 mm	Contact Plate	PB5084C
Columbia Agar Base with 2% Salt		
10x90mm		PO0879A
LB Agar with Ampicillin		
10x90mm		PO5311A
LB Agar with Kanamycin		
10x90mm		PO5309A
Nutrient Agar		
10x90mm		PO0155A
Plate Count Agar		
10x90mm		PO0158A
Plate Count Agar (APHA)		
10x90mm		PO5013A
R2A Agar		
10x90mm		PO0659A^
10x90mm		PO5149A+
Sabouraud Dextrose Agar		
10x90mm		PO0160A
(2x10) 55 mm	Contact Plate	PO0192C
Tryptone Soya Agar		
10x90mm		PO0163A^
10x90mm		PO5012A+
(2x10) 55 mm	Contact Plate	PO0193C
Tryptone Soya Agar with Sheep Blood		
10x90mm		PB5012A
Tryptone Soya Agar (25mL)		
10x90mm	Deep Fill	PO5073A
Tryptone Soya Agar with Disinhibitor		
10x90mm		PO5095A
(2x10) 55 mm		PO5024C
Tryptone Soya Agar with Yeast Extract (TSYE)		
10x90mm		PO5050A
TSA with Disinhibitor PLUS		
(2x10) 55 mm	Contact Plate	PO5172C

* Product is not a stock item, minimum order and lead times may apply.
+ Mainland Europe ^ Only available in UK

Urinary Tract Infection Media

Brilliance UTI Agar		
10x90mm		PO0794A^
10x90mm		PO5120A+
Brilliance UTI Clarity Agar		
10x90mm		PO1110A^
10x90mm		PO5159A+
CLED Agar		
10x90mm		PO0120A^
10x90mm		PO5009A+
CLED Agar with Andrade's		
10x90mm		PO0121A
MacConkey Agar without Salt		
10x90mm		PO0148A^
10x90mm		PO5131A+
MacConkey Agar		
10x90mm		PO0149A
10x90mm		PO5146A
MacConkey No. 2		
10x90mm		PO5110A
MacConkey No. 3		
10x90mm		PO5002A+
10x90mm		PO0495A^

Specialist Pharmaceutical Irradiated Media (Triple Wrapped)

D/E Neutralising Agar		
10x90mm	Deep Fill	PO1128B
Sabouraud Dextrose Agar		
10x55mm	Contact Plate	PO0394D
10x90mm		PO0410B
Tryptone Soya Agar		
10x55mm	Contact Plate	PO0262D
10x90mm	Deep Fill	PO0821B^
10x90mm	Deep Fill	PO5012B+
10x90mm	VHP impermeable packaging	PO0777B
5x140mm		PO0306T
Tryptone Soya Agar with Tween™ and Lecithin		
10x55mm	Contact Plate	PO0479D
Tryptone Soya Agar (Pro-Tect)		
10x55mm	Contact Plate	PO0678D
TSA with Disinhibitor PLUS		
(2x10) 55 mm	Contact Plate	PO5171D
10x90mm		PO0746B

NEW triple wrapped, irradiated plates. Sterile pack with VHP indicator, according to EP/USP

Tryptone Soya Agar		
(10x10) 90mm		PO5500B
(10x10) 55mm	Contact Plate	PO5510D
Tryptone Soya Agar w/Lecithin, Polysorbate 80, Sodium Thiosuphate, L-Histidine		
(10x10) 90mm		PO5501B
(10x10) 55mm	Contact Plate	PO5511D
Sabouraud Dextrose Agar		
(10x10) 90mm		PO5502B*
(10x10) 55mm	Contact Plate	PO5512D*

Sabouraud Dextrose Agar w/Lecithin, Polysorbate 80, Sodium Thiosuphate, L-Histidine		
(10x10) 90mm		PO5503B*
(10x10) 55mm	Contact Plate	PO5513D*

Specialist Pharmaceutical Formulations

Cetrimid-Agar EP, USP		
10x90mm		PO1168A
MacConkey Agar (EP/USP/JP/BP)		
10x90mm		PO1142A
Mannitol Salt Agar (EP/USP/JP/BP)		
10x90mm		PO1169A
Potato Dextrose Agar (EP/USP/JP/BP)		
10x90mm		PO0186A
Sabouraud Dextrose Agar pH 5.6 (EP/USP/JP/BP)		
10x90mm		PO1166A
Violet Red Bile Glucose Agar (EP/USP/JP/BP)		
10x90mm		PO1167A
XLD (EP/USP/JP/BP)		
10x90mm		PO1132A

Water Testing Media

Chromogenic Coliform Agar		
10x90mm		PO5318A
10x55mm		PO5428J
Cetrimide Agar		
10x90mm		PO0185A^
10x90mm		PO5076A+
Enterococcus Selective Agar (BAA) Bile Aesculin Azide Agar		
10x90mm		PO5062A
Kanamycin Aesculin Azide Selective Medium		
10x90mm		PO5059A
Legionella BCYE-Agar W/O L-Cysteine		
10x90mm		PO0255A^
10x90mm		PO5028A+
Legionella BCYE α Medium		
10x90mm		PO5072A
Legionella GVPC Selective Medium		
10x90mm		PO0245A^
10x90mm		PO5074A+
Membrane Clostridium Perfringens (mCP) Agar		
10x90mm		PO5163A
Nutrient Agar		
10x90mm		PO0155A^
10x90mm		PO5025A+
R2A Agar		
10x90mm		PO0659A^
10x90mm		PO5149A+
Slanetz and Bartley Medium (Enterococcus Agar)		
10x90mm		PO5018A
Tergitol 7-Lactose-TTC Agar		
10x90mm		PO5164A
Tergitol 7-Lactose-TTC Medium		
(2x10) 55 mm		PO5411J
TSC Agar		
10x90mm		PO0520A
10x90mm		PO5315A

Disposable Analytical Filters

Single-use Nalgene Filterware

Our analytical filterware provides precise, easy-to-use systems for the recovery of microorganisms and particulates for QC testing and research.

- Disposable, easy-to-use, pre-assembled sterile filter units and funnels
- Easy-to-remove triton-free cellulose nitrate (CN) membranes
- The 0.2 µm CN membrane is suitable for sterility testing
- 0.45 µm CN gridded membrane is certified for water quality testing

Analytical Test Filter Funnels, Sterile

Cat. No.	Upper cap., ml	Membrane dia., mm	Pore size, µm	Membrane color/grid
NG145-0020	100mL	47mm	0.2µm	White/none
NG145-0045	100mL	47mm	0.45µm	White/black
NG145-2020	250mL	47mm	0.2µm	White/none
NG145-2045	250mL	47mm	0.45µm	White/black

Cat. No.	Funnel cap., ml	Membrane dia., mm	Pore size, µm	Membrane color/grid
NG147-0045	100mL	47mm	0.45µm	Gray/black
NG147-2045	250mL	47mm	0.45µm	Gray/black

Analytical Filter, Units, Sterile

Cat. No. 130	Pore size, µm	Membrane color/grid
NG130-4020	0.2µm	White/none
NG130-4045	0.45µm	White/black

Analytical Filter, Funnels, Sterile

Cat. No.	Funnel cap., ml	Membrane dia., mm	Pore size, µm	Membrane color/grid
NG140-4045	150mL	47mm	0.45µm	White/black



Filter Membranes and Accessories

Membrane Filters, Non-sterile, Nylon

Cat. No.	Membrane dia., mm	Pore size, µm	Membrane color
NG02154020	47mm	0.2µm	White
NG02154045	47mm	0.45µm	White

Water Quality Membrane Filter, Sterile CN

Cat. No.	Membrane dia., mm	Pore size, µm	Membrane color
NG02054045		0.45µm	white/black grid
NG02056045		0.45µm	gray/black grid

Membrane Filters, Non-sterile, CA

Cat. No.	Membrane dia., mm	Pore size, µm	Membrane color
NG02104020	47mm	0.2µm	white
NG02104045	47mm	0.45µm	white

Filter Forceps

Cat. No.	Style
NG03990001	Bent Tip
NG03990002	Straight Tip

Filter Funnel Adapter, Non-sterile

Cat. No.
NG03970010

Vacuum Gasket, Non-sterile thermoplastic elastomer

Cat. No.
NG03950708

Filter Stopper, Non-sterile, No. 8 rubber

Cat. No.
NG03960080

Vacuum Manifold

Cat. No.
NG03450001

Reusable Filter Holders and Funnels

Nalgene Reusable Filterware

- reduces waste!

Our analytical filterware provides precise, easy-to-use systems for the recovery of microorganisms and particulates for QC testing and research.

- Replaces metal or breakable glass – major components are made of durable, transparent, autoclavable polysulfone
- Use for sterilization or analysis – filter holders come with sterilization and analytical membrane support plates
- Versatile – designed to be used with vacuum or pressure
- Graduated chambers for easy reference

Filter Holders with Receiver

Cat. No.	Capacity, upper chamber, ml	Capacity, receiver, ml
NG300-4000	250mL	250mL
NG300-4050	500mL	500mL
NG300-4100	500mL	1000mL

Filter Holders with Funnel

Cat. No.	Capacity, upper chamber, ml	Height, incl. cover ports mm
NG03104000	250mL	190mm
NG03104050	500mL	214mm

Reusable Bottle Top Filters

Cat. No.	Capacity, upper chamber, ml	Fits bottle neck size, mm
NG03202545	250mL	45mm
NG03205033	500mL	33mm
NG03205045	500mL	45mm

Filter Funnels with Clamp

Cat. No.	Capacity, ml
NG03150047	250mL



Aeromonas

Aeromonas Medium	
10x90mm	PO0325A*

Bacillus cereus

Bacillus cereus Selective Agar	
10x90mm	PO5048A
Bacillus cereus MYP Agar	
10x90mm	PO5133A

Bordetella

Bordetella Selective Medium	
10x90mm	PB5065A

Burkholderia

Burkholderia cepacia Agar	
10x90mm	PO0938A

Candida

Brilliance Candida Agar		
10x90mm	PO1034A^	
10x90mm	PO5170A+	
Brilliance Candida Agar/Sabouraud GC Agar		
10x90mm	Biplate	PO5258E
Sabouraud-Dextrose Agar		
10x90mm	PO5001A+	
10x90mm	PO0160A^	
10x55mm	Contact Plates	PO0192C
10x90mm	TWI	PO0410B
Sabouraud Dextrose Agar with Chloramphenicol (Deep fill)		
10x90mm	PO0358A	
Sabouraud Dextrose Agar with Chloramphenicol and Actidione		
10x90mm	PO0162A	
Sabouraud-Dextrose Agar with Disinhibitor		
10x90mm	PO5103A	
Sabouraud-Dextrose Agar with Chloram. & Gentamicin		
10x90mm	PO5096A	
Sabouraud-Dextrose Agar with Chloramphenicol		
10x90mm	PO5070A+	
10x90mm	PO0161A^	
10x55mm	Contact Plates	PO5094C

Campylobacter

Brilliance CampyCount Agar	
10x90mm	PO1185A
Campylobacter CAT Agar	
10x90mm	PO0839A
Campylobacter Selective Agar (Skirrow)	
10x90mm	PB0118A
Campylobacter Selective Agar (Butzler)	
10x90mm	PB5006A
Campylobacter Selective Blood Free Agar (CCDA)	
10x90mm	PO0119A^
10x90mm	PO5091A+
Karmali Selective Medium	
10x90mm	PO5041A
Skirrow Selective Medium	
10x90mm	PO5040A

Karmali Selective Medium/Karmali Selective Medium		
10x90mm	Biplate	PO5219E

Cemo

CEMO Medium with Fungizone™ and Streptomycin		
10x90mm		PO0443A*
CEMO Medium with Fungizone™		
10x90mm		PO0444A*

Clostridium difficile

Brazier's Clostridium difficile Selective Agar		
10x90mm		PB1055A
Clostridium difficile Selective Agar		
10x90mm		PB0218A*
Clostridium difficile Selective Medium		
10x90mm		PB5054A

Clostridium perfringens

Membrane Clostridium Perfringens (mCP) Agar		
10x90mm		PO5163A
TSC Agar		
10x90mm		PO0520A
10x90mm		PO5315A

Coliforms/*E. coli*

Brilliance E. coli/coliform Agar		
10x90mm		PO5176A
Brilliance ESBL Agar		
10x90mm		PO5302A
Chromogenic Coliform Agar		
10x90mm		PO5318A
10x55mm		PO5428J
China Blue Lactose Agar		
10x90mm		PO5060A
MLGA		
10x90mm		PO1016A
TBX Agar		
10x90mm		PO0727A^
10x90mm		PO5109A+

Corynebacteria

Hoyles Medium		
10x90mm		PO0143A

Dermatophytes

Dermasel Agar		
10x90mm		PO0737A^
10x90mm		PO5037A+
Dermatophyte Medium with Phenol Red		
10x90mm		PO0166A
Dermatophytes Selective Agar (Taplin)		
10x90mm		PO5087A
Dermatophyte Medium with Chloramphenicol (100mg / Litre)		
10x90mm		PO0964A

Enterobacteriaceae

Brilliance CRE Agar		
10x90mm		PO1226A
Brilliance ESBL/Brilliance CRE Agar		
10x90mm	Biplate	PO1265E
Brilliance ESBL Agar		
10x90mm		PO5302A
Endo Agar		
10x90mm		PO5005A
Eosine Methylene Blue Agar (modified) Levine		
10x90mm		PO5045A
MacConkey Agar		
10x90mm		PO5146A
MacConkey Agar (EP, USP, JP, BP)		
10x90mm		PO1142A
MacConkey Agar No. 3		
10x90mm		PO0495A^
10x90mm		PO5002A+
(2x10) 55mm	Contact Plate	PO5053C
MacConkey Agar No.2		
10x90mm		PO5110A
MacConkey Agar with Salt		
10x90mm		PO0149A
MacConkey Agar without Salt		
10x90mm		PO0148A^
10x90mm		PO5131A+
Violet Red Bile Agar		
10x90mm		PO5075A
Violet Red Bile Agar with MUG		
10x90mm		PO5031A
Violet Red Bile Glucose Agar		
10x90mm		PO5322A
(2x10) 55mm	Contact Plate	PO5323C
10x90mm		PO1167A
C.L.E.D. Medium / MacConkey Agar No. 3		
10x90mm	Biplates	PO5217E
Columbia Agar with Blood / Endo Agar		
10x90mm	Biplates	PB5200E
Columbia Agar with Blood / MacConkey Agar No. 3		
10x90mm	Biplates	PB5207E
Columbia CNA Aesculin Selective Agar/ Brilliance™ UTI		
10x90mm	Biplates	PB5220E
Columbia CNA Aesculin Selective Agar/ MacConkey Agar No. 3, mod.		
10x90mm	Biplates	PB5224E
Enterococci		
Bile Aesculin Agar		
10x90mm		PO0169A
Aesculin Blood Agar (modified)		
10x90mm		PB5023A
Brilliance VRE Agar		
10x90mm		PO1175A
Enterococcus Selective Agar (Bile Aesculin Azide Agar)		
10x90mm		PO5062A
Kanamycin Aesculin Azide Agar		
10x90mm		PO0173A^
10x90mm		PO5059A+
Slanetz and Bartley Agar		
10x90mm		PO5018A

VRE Selective Agar		
10x90mm		PO5089A

Enterohaemorrhagic *E. coli*

Sorbitol MacConkey Agar with Cefixime Tellurite		
10x90mm		PO0702A

Enterohaemolysin Agar with Blood		
10x90mm		PB5105A

Sorbit-MacConkey-Nährboden (SMAC)		
10x90mm		PO5069A
10x90mm		PO0232A

Extended Spectrum Beta-Lactamase

Brilliance ESBL Agar		
10x90mm		PO5302A

Gardnerella vaginalis

Gardnerella Selective Agar with Sheep Blood		
10x90mm		PB0134A

Gardnerella vaginalis agar with Human Blood		
10x90mm		PB5067A

Haemophilus influenzae

Columbia Agar with Chocolate Horse Blood and Bacitracin		
10x90mm		PB0220A

Chocolate Agar with Vitox		
10x90mm		PO5090A

Chocolate GC Selective Agar		
10x90mm		PB0963A

Columbia Agar with Chocolate Horse Blood		
10x90mm		PB0124A

G.C. Selective Agar (VCNT Chocolate)		
10x90mm		PB0136A

G.C. Selective Agar (VCNT) with Lysed Blood		
10x90mm		PB0135A

G.C. Selective Agar with Lysed Horse Blood and LCAT		
10x90mm		PB0226A

G.M.P. G.C. Agar (Lysed) V.C.A.T.		
10x90mm		PB0820A

Haemophilus Selective Agar (Chocolate Agar with Vitox, Bacitracin, Vancomycin, Amphotericin B)		
10x90mm		PO5097A

Lysed G.C. Selective Agar		
10x90mm		PB1205A

Columbia Agar with Blood / Chocolate Agar		
10x90mm	Biplate	PB5202E

Helicobacter pylori

Helicobacter pylori Selective Agar		
10x90mm		PB0398A

Lactobacilli

MRS Agar		
10x90mm		PO0231A^
10x90mm		PO5047A+

Legionella

Legionella Medium BCYE		
10x90mm		PO5072A

Legionella Medium BCYE (without L-cysteine)		
10x90mm		PO0255A^

Legionella BCYE α Medium		
10x90mm		PO5072A

Legionella Selective Medium BMPA		
10x90mm		PO0324A^

Legionella Selective Medium GVPC		
10x90mm		PO0245A^

Legionella Selective Medium MWY		
10x90mm		PO5074A+

Legionella Selective Medium MWY		
10x90mm		PO5071A

Listeria

Brilliance Listeria Agar		
10x90mm		PO1102A^

Listeria Selective Agar (Oxford)		
10x90mm		PO0179A^

Oxid Chromogenic Listeria Agar (ISO) – OLCA (ISO)		
10x90mm		PO1196A^

Palcam Medium		
10x90mm		PO5104A

Tryptone Soya Agar with Yeast Extract (TSYE)		
10x90mm		PO5050A

MRSA

Brilliance MRSA 2 Agar		
10x90mm		PO1210A^

Brilliance™ MRSA 2 Agar/ Columbia Agar with Sheep Blood PLUS		
10x90mm	Biplate	PB5253E

Brilliance MRSA 2 Agar/ Brilliance Staph 24 Agar		
10x90mm	Biplate	PO1258E

Mycoplasma/Ureaplasma

Mycoplasma / Ureaplasma Agar		
10x90mm		PO5081A

Neisseria gonorrhoeae

Oxid Lysed GC Selective Agar		
10x90mm		PB1205A

Oxid Chocolate GC Selective Agar		
10x90mm		PB0963A

G.C. Selective Agar (VCNT Chocolate)		
10x90mm		PB0136A

G.C. Selective Agar with Lysed Horse Blood and LCAT		
10x90mm		PB0226A

G.M.P. G.C. Agar (Lysed) V.C.A.T.		
10x90mm		PB0820A

Neisseria Selective Medium PLUS		
10x90mm		PO5004A

Pasteurella

Pasteurella Selective Medium		
10x90mm		PB5175A

Pseudomonas

Pseudomonas CFC Selective Agar		
10x90mm		PO0291A^
10x90mm		PO5132A+
Pseudomonas CN Selective Agar		
10x90mm		PO0185A^
10x90mm		PO5076A+
Cetrimide Agar (USP, EP)		
10x90mm		PO5181A
Glutamate Starch Phenol Red Selective Medium		
10x90mm		PO5128A

Salmonella/Shigella

Brilliance Salmonella Agar		
10x90mm		PO5098A
Brilliant Green Agar (Modified)		
10x90mm		PO0171A^
10x90mm		PO5033A+
DCA (Hynes)		
10x90mm		PO0126A^
10x90mm		PO5016A+
DCLS Agar		
10x90mm		PO0127A
Gassner Medium		
10x90mm		PO5021A
Hektoen Enteric Agar		
10x90mm		PO5100A
MLCB Agar		
10x90mm		PO5029A
S. S. Agar		
10x90mm		PO5022A
Salmonella Chromogenic Medium (OSCM)		
10x90mm		PO0958A
X.L.D. Agar (EP, USP, JP, BP)		
10x90mm		PO1132A
XLD Medium		
10x90mm		PO0164A^
10x90mm		PO5057A+
X.L.T4 Medium		
10x90mm		PO1158A^
10x90mm		PO5116A+
Brilliance™ Salmonella/XLD Agar		
10x90mm	Biplate	PO5248E
Hektoen Enteric Agar/DCA (Leifson)		
10x90mm	Biplate	PO5257E
S.S. Agar / X.L.D. Medium		
10x90mm	Biplate	PO5210E

Staphylococcus aureus

Baird-Parker Agar		
10x90mm		PO0168A^
10x90mm		PO5014A+
Baird Parker RPF Agar		
10x90mm		FR59951

Brilliance MRSA 2 Agar		
10x90mm		PO1210A^
10x90mm		PO5310A+
Brilliance MRSA 2 Agar/Columbia Sheep Blood Agar		
10x90mm	Biplate	PB5253E
Brilliance Staph 24 Agar		
10 plates		PO1186A^
Brilliance MRSA 2/Brilliance Staph 24		
10x90mm	Biplate	PO1258E
CNA Staph/Strep Selective Agar		
10x90mm		PB0308A
10x90mm		PB5049A
Columbia CAP with Sheep Blood		
10x90mm		PB5028A
DNase Agar		
10x90mm		PO0128A
DNase Agar with Methyl Green		
10x90mm		PO1000A
Mannitol Salt Agar		
10x90mm		PO0151A^
10x90mm		PO5027A+
Mannitol Salt Agar (EP, USP, JP, BP)		
10x90mm		PO1169A

Staphylococci and Streptococci

CNA Staph/Strep Selective Agar		
10x90mm		PB0308A
10x90mm		PB5049A
Brilliance Staph 24 Agar		
10x90mm		PO1186A
Brilliance GBS Agar		
10x90mm		PO5320A
10x90mm	Biplate	PB5260E
Columbia CAP Selective Agar with Sheep Blood		
10x90mm		PB5082A
Columbia CAP Selective Agar with Horse Blood		
10x90mm		PB1161A
Streptococcal Selective Agar COBA		
10x90mm		PB0298A

Streptococcus agalactiae

Brilliance GBS Agar		
10x90mm		PO5320A
10x90mm	Biplate	PB5260E
Edwards Medium with Sheep Blood (modified)		
10x90mm		PB5080A

Vibrio

TCBS Cholera Medium		
10x90mm		PO0194A

VRE


Brilliance VRE Agar		
10x90mm		PO1175A
VRE Selective Agar		
10x90mm		PO5089A

Yeasts and Moulds

Biggy Agar		
10x90mm		PO5011A
Brilliance™ Candida Agar		
10x90mm		PO1034A^
10x90mm		PO5170A+
Brilliance Candida / Sabouraud Dextrose Agar		
10x90mm	Biplate	PO5258E
Dichloran-Glycerol (DG18) Selective Medium		
10x90mm		PO5088A
(2x10) 55 mm	Contact Plate	PO5313C
Kimmig Medium		
10x90mm		PO5019A
DRBC (ISO) Agar		
10x90mm		PO1227A
Malt Extract Agar		
10x90mm		PO0182A^
10x90mm		PO5055A+
(2x10) 55 mm	Contact Plate	PO5079C
Malt Extract Agar with Chloramphenicol		
(2x10) 55 mm	Contact Plate	PO5314A
Oxytetracycline Glucose Yeast Extract Agar		
10x90mm		PO0183A
Potato Dextrose Agar		
10x90mm		PO0186A
Rice Agar		
10x90mm		PO5064A
Rose Bengal Chloramphenicol Agar		
10x90mm		PO0214A
Sabouraud Dextrose Agar		
10x90mm		PO0160A^
10x90mm		PO5001A+
Sabouraud Dextrose Agar with Chloramphenicol		
10x90mm		PO0161A^
10x90mm		PO5070A+
(2x10) 55 mm	Contact Plate	PO5094C
Sabouraud Dextrose Agar with Chloramphenicol (Deep Fill)		
10x90mm		PO0358A
Sabouraud Dextrose Agar with Chloramphenicol and Actidione		
10x90mm		PO0162A
Sabouraud Glucose Agar with Disinhibitor		
10x90mm		PO5103A
Sabouraud Glucose Selective Agar with Gentamicin and Chloramphenicol		
10x90mm		PO5096A
Wort Agar		
10x90mm		PO5063A
Yeast Extract Glucose Chloramphenicol Medium (YGC-Medium)		
10x90mm		PO5032A
Dichloran-Glycerol (DG18) Selective Medium		
(2x10) 55 mm	Contact Plate	PO5313C
Sabouraud Dextrose Agar		
10x55mm	Contact Plate	PO0192C
Sabouraud Dextrose Agar		
10x55mm	Contact Plate	PO0394D
Brilliance Candida Agar/Sabouraud Dextrose Agar		
10x90mm	Biplate	PO5258E

Yersinia

Yersinia Selective Agar CIN		
10x90mm		PO0287A^
10x90mm		PO5044A+



Ready-to-Use Media

Saves Time and Effort

Thermo Scientific™ Oxoid™ Prepared Media makes laboratory life easier. Our quality assured, ready-to-use media, made using Oxoid dehydrated culture media, eliminates the time-consuming processes of in-house media preparation and quality control.

Ready Poured Plates – Large range includes many products in the Thermo Scientific™ *Brilliance*™ Range of Chromogenic Media, for rapid reliable detection and differentiation of important pathogens.

Bottled Media – A wide range of general purpose media; specialist pharmaceutical media; diluents, water and peptones are available in a variety of bottle types – see page 47 for more details.

Bagged Enrichment Media – Thermo Scientific™ Bagged Enrichment Media are designed for media preparation in minutes. This scalable range of solutions eliminates the need for weighing and autoclaving of dehydrated culture media, enabling you to quickly meet testing demands. See pages 56-57 for details.

Quality control testing of Oxoid Prepared Media manufactured at our facility in Perth is validated by UKAS (UKAS testing laboratory number 2727), meeting the requirements of ISO/IEC 17025:2005–General requirements for the Competence of Testing and Calibration Laboratories.

Prepared Media QC certificates are available via thermofisher.com/microbiology for our customers' convenience.



Thermo Scientific™ Bottled Media

High Performance Media in Formats
to Suit Your Laboratory



225ml in 250ml Sirop

Broths of various volumes available in
60, 125, 250, 300, 500 and 1,000ml Sirops



100ml in 125ml Sirop

Agars of various volumes available in
60, 125, 250, 300, 500 and 1,000ml Sirops



10ml Agar Slope in 30ml Universal

Agar Slopes available in various volumes



10ml Broth in 30ml Universal with Durham's Tube

Broths available in various volumes



3ml Agar Slope in Bijou

Broths and Agars available in 3 - 5ml volumes



100ml in 100ml Vial (injectable cap)

Broths and diluents available in various size
vials 2 - 100ml



10ml in Tube

Broths and Agar Slopes available in various volumes



100ml in 100ml DIN (injectable cap)

Broths and diluents available in various size DIN
bottles 100 - 500ml

KEY TO BOTTLE TYPES

D	DIN – wide neck bottle with septum
V	Vial – narrow neck with septum
S	Sirop – screw cap bottle
Spc	Sirop – with pierceable screw cap
U	Universal – 1oz straight walled
B	Bijou – 1/4oz straight walled
T	Tube with screw cap
P	Sterile plastic container

Product/Pack Size	Bottle Volume	Bottle Type	Code
-------------------	---------------	-------------	------

Specialist Pharmaceutical Media

0.1% Peptone Water

10x100mL	100mL	D	BO0833M
10x300mL	500mL	D	BO0833Z
10x500mL	500mL	D	BO0833V

0.1% Peptone Water with 0.1% Tween™

10x100mL	125mL	Spc	BO0293M
10x300mL	500mL	Spc	BO0293Z
10x500mL	500mL	D	BO0964V

Buffered Sodium Chloride Peptone (EP/USP/JP/BP)

24x9mL		U	BO0322D
10x90mL	125mL	S	BO0322J*
10x100mL	125mL	S	BO0322M

Thioglycollate Medium (acc. to EP/USP)

10x100mL	100mL	V	BO0368M
10x100mL	100mL	D	BO0510M
10x500mL	500mL	D	BO0510V
50x9mL	9mL	T	TV5001D

Tryptone Soya Broth (acc. to EP/USP)

24x2mL	5mL	V	BO0369A
24x10mL	20mL	V	BO0369E
10x100mL	100mL	V	BO0369M
10x100mL	100mL	D	BO0509M
10x500mL	500mL	D	BO0509V
50x10mL	20mL	T	TV5002E

Specialist Pharmaceutical Selective and Differential Media

Enterobacteria Enrichment Broth-Mossel (EP/USP/JP/BP)

24x10mL		U	BO1156E
10x90mL	125mL	S	BO1156J
10x100mL	125mL	S	BO1156M

MacConkey Broth (EP/USP/JP/BP)

10x100mL	125mL	S	BO1124M
----------	-------	---	---------

RVS Enrichment Broth MLT (EP/USP/JP/BP)

24x10mL		U	BO1121E
---------	--	---	---------

Reinforced Clostridial Medium (EP/USP/JP/BP)

10x100mL	125mL	S	BO1158M
----------	-------	---	---------

Sabouraud Dextrose Agar pH 5.6 (EP/USP/JP/BP)

10x100mL	125mL	S	BO1155M
10x250mL	300mL	S	BO1155T
10x450mL	500mL	S	BO1155Z

Sabouraud Dextrose Broth (EP/USP/JP/BP)

10x100mL	125mL	S	BO0283M
----------	-------	---	---------

Product/Pack Size	Bottle Volume	Bottle Type	Code
-------------------	---------------	-------------	------

General Purpose Media

Brain Heart Infusion Broth			
24x9mL		U	BO1230D
Brain Heart Glucose Bouillon			
10x100mL		U	BO1230M
10x125mL		U	BO1230Z
50x10mL		T	TV5090E
Brilliant Green Bile Broth (Durham's Tube)			
24x10mL		U	BO0345E
50x10mL		T	TV5009E
Buffered Peptone Water			
10x90mL	125mL	S	BO0201Z
10x225mL	225mL	S	BO0201S
10x225mL	225mL	P	BO0688S
Buffered Peptone Water (ISO)			
50x9mL		T	TV5013D
Columbia Agar with Chocolated Horse Blood			
24x3mL		B	BO0341B
Cooked Meat Medium			
24x10mL		U	BO0212E
100x10mL		U	EB0212E
Nutrient Broth with Dextrose			
50x9mL		T	TV5003D
Lethen Broth			
24x10mL		U	BO0409E
Nutrient Agar (Slope)			
24x3mL		B	BO0336B
200x3mL		B	EB0336B
24x10mL		U	BO0336E
50x7,5mL		T	TV5011Z
Nutrient Broth			
24x3mL		B	BO0210B
24x10mL		U	BO0210E
10x250mL	500mL	S	BO0210T*
Plate Count Agar			
10x250mL	500mL	S	BO0195T
10x500mL	500mL	S	BO0195V
Thioglycollate Medium USP (acc. To ESP/USP)			
24x20mL		U	BO0211G
10x100mL	125mL	S	BO0211M
50x9mL		T	TV5001D
Tryptone Soya Agar			
10x100mL	125mL	S	BO0330M
10x250mL	500mL	S	BO0330T
10x500mL	500mL	S	BO0330V
Tryptone Soya Agar (Slope)			
24x15mL		U	BO0405F
Tryptone Soya Broth (acc. to EP/USP)			
24x10mL		U	BO0351E
24x20mL		U	BO0351G
10x90mL		S	BO0351U
10x100mL	125mL	S	BO0351M
10x200mL	250mL	S	BO0351R
10x500mL	500mL	S	BO0351V
50x10mL		T	TV5002E
Yeast Extract Agar			
10x100mL	125mL	S	BO0635M
10x500mL	500mL	S	BO0556V

Product/Pack Size	Bottle Volume	Bottle Type	Code
-------------------	---------------	-------------	------

Diluents, Water and Peptones

0.1% Peptone Water with 0.85% Salt			
24x9mL		U	BO0471D
0.9% Saline			
24x3mL		B	BO0334B
200x3mL		B	EB0334B
24x10mL		U	BO0334E
Alkaline Peptone Water			
24x10mL		U	BO0335E
DEV Agar			
50x18ml		T	TV5207X
Lactose Peptone Water Sugar (Durham's Tube)			
24x3mL		B	BO0435B
Maximum Recovery Diluent			
24x9mL		U	BO0348D
50x9mL		T	TV5016D
10x90mL	125mL	S	BO0348Z
10x225mL	250mL	S	BO0348S*
10x500mL	500mL	S	BO0348V
Peptone Water			
24x3mL		B	BO0208B
200x3mL		B	EB0208B
24x9mL		U	BO0208D
24x3mL		B	BO0208B
Purified Water			
24x3mL		B	BO0209B
200x3mL		B	EB0209B
300x3mL		T	BO0184B
24x5mL		B	BO0209C
200x5mL		B	EB0209C
24x10mL		U	BO0209E
10x100mL		S	BO0209M
10x200mL		S	BO0209R
Ringers Solution 1/4 strength			
24x9mL		U	BO0332D
Tryptone Water			
24x3mL		B	BO0383B
24x5mL		B	BO0383C

Water testing media

DEV-Agar			
50 x 18 mL		T	TV5207X
10 x 250 mL		S	BO1251T
Minerals Modified Glutamate (Durham Tube)			
24x10		U	BO0541E
Minerals Modified Glutamate x2 strength (Durham Tube)			
24x10		U	BO0542E
Water Plate Count Agar (ISO)			
10x100 mL		S	BO0055M
10x200 mL		S	BO0055R
TSC Agar (Clostridium perfringens Agar)			
50x20 mL		T	TV5204G
10x100 mL	125 mL	S	BO0634M

Organism Specific

Anaerobes			
Schaedler Broth (with Haemin and Vitamin K)			
50x9 mL		T	TV5008D
TSC Agar			
50x20mL		T	TV5204G

Product/Pack Size	Bottle Volume	Bottle Type	Code
Campylobacter			
Bolton Broth			
10x225mL	250mL	S	BO1070S
Clostridium perfringens			
Lactose Gelatin Medium			
24x15mL		U	BO1068F
Buffered Motility Nitrate Medium			
24x10mL		U	BO1069E
Coliforms/E. coli			
Brilliant Green Bile Broth (Durham's Tube)			
24x10mL		U	BO0345E
50x10mL		T	TV5009E
E.E. Broth			
10x100mL		S	BO0598M
10x90mL		S	BO0443Z
50x10mL		T	TV5041E
Lactose Broth			
10x90mL	125mL	S	BO0596Z
Lauryl Tryptose Broth with MUG and Tryptophan (with Durham Tube)			
50x10mL		T	TV5079E
MacConkey Broth			
10x100mL	125mL	S	BO0550M
MacConkey Broth Purple			
10x100mL	125mL	S	BO0376M
MacConkey Broth Purple (Durham's Tube)			
24x10mL		U	BO0347E
Membrane Lauryl Sulphate Broth			
24x10mL		U	BO0394E
Minerals Modified Glutamate (Durham's Tube)			
24x10mL		U	BO0541E*
Minerals Modified Glutamate x2 Strength (Durham's Tube)			
24x10mL		U	BO0542E
Violet Red Bile Agar			
10x100mL	125mL	S	BO0455M
Enterobacteriaceae			
EE Broth			
10x90mL	125mL	S	BO0443Z
10x100mL	125mL	S	BO0598M
50x10mL		T	TV5041E
Kligler Iron Slant Agar (with Urea)			
50x9ml		T	TV5004D
S.I.M. Medium			
50x10ml		T	TV5014E
Simmons Citrate Slant Agar			
50x7,5ml		T	TV5015Z
Triple Sugar Iron Agar			
50x9mL		T	TV5074D
Urea Agar (Slope)			
24x3mL		B	BO0337B
200x3mL		B	EB0337B
Urea Broth			
24x3mL		B	BO0338B
50x6mL		T	TV5007N
Violet Red Bile Agar			
10x100mL	125mL	S	BO0455M
Violet Red Bile Glucose Agar			
10x100mL	125mL	S	BO1329M

Product/Pack Size	Bottle Volume	Bottle Type	Code
Enterococci			
Bile Aesculin Agar (Slope)			
24x3mL		B	BO0339B
Aesculin Agar (Slope)			
24x3mL		B	BO0826B
Enterohaemorrhagic E. coli			
Modified Tryptone Soya Broth with Novobiocin			
100x10mL		U	BO0869E
Group B Streptococci			
Group B Strep Selective Enrichment (LIM) Broth			
24x10mL		U	BO1086E
50x5mL		T	TV5093E
Listeria			
Fraser Broth			
50x10mL		T	TV5020E
Half Fraser Broth			
10x225mL	250mL	S	BO0350S
10x225mL	250mL	P	BO0793S
ONE Broth-Listeria			
10x225mL	250mL	S	BO1066S
LEB-Listeria Enrichment Broth			
10 x 225mL	250mL	S	BO1205S
Mycoplasma			
Mycoplasma / Ureaplasma Enrichment Broth			
50x2mL		T	TV5081A
Salmonella/Shigella			
Buffered Peptone Water			
50x9mL		T	TV5013D
Lysine Decarboxylase Broth (Taylor)			
50x6mL		T	TV5028N
MKTT(n)			
50x10mL		T	TV5065E
ONE Broth-Salmonella			
10x225mL	250mL	S	BO1096S
Rappaport-Vassiliadis Broth			
24x10mL		U	BO0203E
50x10mL		T	TV5017E
Rappaport-Vassiliadis Soya Peptone Broth			
24x10mL		U	BO0499E*
50x10mL		T	TV5036E
Selenite Cystine Broth			
100x10mL		U	EB1222E
50x10ml		T	TV5018E
Selenite F Broth			
100x10mL		U	EB1221E
Selenite Lactose Broth			
50x10mL		T	TV5005I
Tetrathionate Broth (Mueller Kaufmann)			
50x10ml		T	TV5006I
Staph. aureus/MRSA			
Contrast MRSA Broth			
100x3mL in glass		U	EB1225B
100x3mL in plastic tubes			EB1226B
24x3ml in plastic tubes		U	BO1226B
Nutrient Broth with 7.5% Salt			
100x10mL		U	EB0424E

* Product is not a stock item, minimum order and lead times may apply. Please contact your Customer Services Administrator for details.

Not all products available in all territories. Please contact your local office for details.

Culture Media

Bottled Media, Bagged Media

Product/Pack Size	Bottle Volume	Bottle Type	Code
Trichomonads			
Trichomonas Medium			
200x5mL		B	EB0861C
24x5mL		B	LR0027A
Yeasts and Moulds			
Actidione Agar			
100mL			PM0118C
Rose Bengal Chloramphenicol Agar			
10x250mL	300mL	S	BO0198T*
Sabouraud Dextrose Agar			
10x100mL	125mL	S	BO0408M
10x250mL	500mL	S	BO0408T
10x500mL	500mL	S	BO0408V
Sabouraud Dextrose Agar (Slope)			
24x10mL		U	BO0342E
Sabouraud Dextrose Agar with Chloramphenicol			
10x100mL	125mL	S	BO0756M
Sabouraud Liquid Medium USP			
24x10mL		U	BO0358E
10x100mL	125mL	S	BO0283M*

Agar Slopes

Bile Aesculin Agar (Slope)			
24x3mL		B	BO0339B
Citrate Agar (Slope)			
24x3mL		B	BO0379B
Columbia Agar with Chocolate Horse Blood (Slope)			
24x3mL		B	BO0341B
Nutrient Agar (Slope)			
24x3mL		B	BO0336B
200x3mL		B	EB0336B
24x10mL		U	BO0336E
Sabouraud Dextrose (Slope)			
24x10mL		U	BO0342E
Triple Sugar Iron Agar (Slope)			
50x9mL		T	TV5074D
Tryptone Soya Agar (Slope)			
24x15mL		U	BO0405F
Urea Agar (Slope)			
24x3mL		B	BO0337B
200x3mL		B	EB0337B



Bagged Media

Thermo Scientific™ Dry-Bags™

Enrichment Media

for maximum convenience, each bag is fitted with a filling filter; each bag makes 20 litres

Buffered Peptone Water	
10 bags	DB0509M
Buffered Peptone Water (ISO)	
10 bags	DB1049M
Half Fraser Broth	
5 bags	DB0895L
Maximum Recovery Diluent	
10 bags	DB0733M
Tryptone Soya Broth	
10 bags	DB0129M
ONE Broth Salmonella	
10 bags	DB1091M
ONE Broth Listeria	
5 bags	DB1066L
ONE broth Listeria with filter	
10 Bags	DB1066M
ONE broth Salmonella with filter	
5 Bags	DB1091L

Dry-Bags Enrichment Media without Filters

Supplied without filling filters in place for a more economical option. Use with Thermo Scientific Dry-Bags Filters (DB0200A). Each bag makes 20 litres.

Buffered Peptone Water	
10 bags	DB0509W
Buffered Peptone Water (ISO)	
10 bags	DB1049W
Maximum Recovery Diluent	
10 bags	DB0733W
Half Fraser Broth	
5 bags	DB0895V
24 Listeria Enrichment Broth without filter	
5 bags	DB1107V
Lactose Broth	
10 bags	DB0137W
Buffered Listeria Enrichment Broth	
5 bags	DB0897V
Tryptone Soya Broth	
10 bags	DB0129W
ONE Broth Salmonella	
10 bags	DB1091W
ONE Broth Listeria	
5 bags	DB1066V

Thermo Scientific™ Dry-Bags Filter

a 0.2µm, autoclavable filter for use with Dry-Bags Enrichment Media. Autoclave before use.

Dry-Bags Filter	
1 filter	DB0200A

Thermo Scientific™ FitBag™ Enrichment Media

Tryptone Soya Broth	
30 x 2.7L Bags	DF0129A
Tryptone Soya Broth	
20 x 4.5L Bags	DF0129B
Tryptone Soya Broth	
10 x 9.0L Bags	DF0129C
Buffered Peptone Water	
30 x 2.7L Bags	DF0509A
Buffered Peptone Water	
20 x 4.5L Bags	DF0509B
Buffered Peptone Water	
10 x 9.0L Bags	DF0509C
UVM1 Broth	
30 x 2.7L Bags	DF0863A
UVM1 Broth	
20 x 4.5L Bags	DF0863B
UVM1 Broth	
10 x 9.0L Bags	DF0863C
Half Fraser Broth	
30 x 2.7L Bags	DF0895A
Half Fraser Broth	
20 x 4.5L Bags	DF0895B
Half Fraser Broth	
10 x 9.0L Bags	DF0895C
Buffered Listeria Enrichment Broth	
30 x 2.7L Bags	DF0897A
Buffered Listeria Enrichment Broth	
20 x 4.5L Bags	DF0897B
Buffered Listeria Enrichment Broth	
10 x 9.0L Bags	DF0897C
Buffered Peptone Water (ISO)	
30 x 2.7L Bags	DF1049A
Buffered Peptone Water (ISO)	
20 x 4.5L Bags	DF1049B
Buffered Peptone Water (ISO)	
10 x 9.0L Bags	DF1049C
24 Listeria Enrichment Broth	
30 x 2.7L Bags	DF1107A
24 Listeria Enrichment Broth	
20 x 4.5L Bags	DF1107B
24 Listeria Enrichment Broth	
10 x 9.0L Bags	DF1107C

Thermo Scientific™ QuickBag™ Enrichment Media

Tryptone Soya Broth	
3 x 2.7L Bags	DQ0129A
Tryptone Soya Broth	
2 x 4.5L Bags	DQ0129B
Buffered Peptone Water	
3 x 2.7L Bags	DQ0509A
Buffered Peptone Water	
2 x 4.5L Bags	DQ0509B
UVM1 Broth	
3 x 2.7L Bags	DQ0863A
UVM1 Broth	
2 x 4.5L Bags	DQ0863B
Half Fraser Broth	
3 x 2.7L Bags	DQ0895A
Half Fraser Broth	
2 x 4.5L Bags	DQ0895B
Buffered Listeria Enrichment Broth	
3 x 2.7L Bags	DQ0897A
Buffered Listeria Enrichment Broth	
2 x 4.5L Bags	DQ0897B
Buffered Peptone Water (ISO)	
3 x 2.7L Bags	DQ1049A

Buffered Peptone Water (ISO)	
2 x 4.5L Bags	DQ1049B
24 Listeria Enrichment Broth	
3 x 2.7L Bags	DQ1107A
24 Listeria Enrichment Broth	
2 x 4.5L Bags	DQ1107B

Peristaltic Pumps

Simple pump for filling up to 2 Dry-Bags at once		
Each	Starter pack - includes peristaltic pump one pump head and tubing	DB2000A
Complex pump for filling and dispensing		
Each		DB3000A
Additional pump head and fittings		
Each		DB1500A
 tubing for peristaltic pump		
Each		DB1750A

Thermo Scientific™ ReadyBags™ Bagged Media

each Thermo Scientific ReadyBags bag contains 3 litres of medium ready to use

Buffered Peptone Water	
4x3 litres	BM0104T
Buffered Peptone Water (ISO)	
4x3 litres	BM1104T
Maximum Recovery Diluent	
4x3 litres	BM0204T

Bag Accessories

only suitable for use with Thermo Scientific ReadyBags Bagged Media

Bag Connectors (autoclavable)	
25	BM9901A

suitable for use with ReadyBags and Dry-Bags Bagged Media; can be autoclaved for reuse



Thermo Scientific™ DiluFlux™ Pro Automated Gravimetric Dilutor

Automated Diluter with dual pump	
DiluFlux Pro-3 Automated Gravimetric Dilutor w/ simple (non-robotic) dispensing arm	DB4100A
DiluFlux Pro-5 Automated Gravimetric Dilutor w/ robotic dispensing arm	DB4150A
Consumables	
BagFilter™ 400 sterile side-filter bag filter bag	DB4011A
BagPage™ 400 sterile full-surface filter bag	DB4012A
BagLight™ 400 sterile plain bag	DB4013A
RollBag™ 1300 sterile sampling bag	DB4014A

Culture Media

Dip Slides, Biochemical Reagents

Options & accessories for Diluflux Automated Dilutor	
Pedal	DB4007A
Trocar	DB4002A
5L autoclavable bottle	DB4005A
500g calibrated COFRAC weight	DB4006A
Dispensing platform	DB4001A
Stand for broth bags	DB4003A
Mobile tube holding arm	DB4004A
Printer	DB4008A
Stainless steel nozzle	DB4009A
Tube weight	DB4010A
Ink Ribbon for Printer	DB4101A
Thermal Printer	DB4102A
BagOpen™ 400 for Diluflow™	DB4103A
Weighing platform	DB4104A
Tongs platform	DB4105A
Adjustable distribution arm system for DiluFLux Pro-3 Instrument	DB4106A
Jumbo Kit BagOpen 3500 for DiluFlux Pro-5 Instrument	DB4107A
GeckoGrip	DB4108A
Stand for One broth bag	DB4109A
Booster kit	DB4110A
Full dispensing small assembly for booster Ø 3,2 mm	DB4111A
Small silicone tubing for booster Ø 3,2 mm	DB4112A
Independent stand for broth bags	DB4113A
Air filter for GL45 cap	DB4114A
Nozzles tubing Ø 3,2 mm	DB4115A
Full dispensing assembly for booster kit Ø 6,4 mm	DB4116A
Silicone tubing for booster Ø 6,4 mm	DB4117A
Full dispensing assembly Ø 6,4 mm	DB4118A
Silicone tubing Ø 6,4 mm	DB4119A
Full dispensing assembly Ø 3,2 mm	DB4120A
Silicone tubing Ø 3,2 mm	DB4121A
Set of silicone tubing for pump Ø 3,2 mm	DB4122A
Set of silicone tubing for pump Ø 6,4 mm	DB4123A
Adjustable distribution arm system for DiluFlux Pro-5 Instrument	DB4124A



Homogenizer Laboratory Blender

Homogenizer Laboratory Blender	
Each	DB5000A
Options & accessories for Homogenizer Laboratory Blender	
BagOpen™ 400, bag opener stand for holding bags	DB5001A
BagRack™ 400, stainless steel stackable storage racks for 10 bags	DB5002A
BagClip™ 400, closing clips for lab blender bags 400 mL	DB5003A



Thermo Scientific™ Dip Slide™ Media

Store at 10–25°C

MacConkey Agar / CLED Medium	
10 slides	DS0103A
100 slides	DS0102B
10 slides	GFD01A*
MacConkey Agar No. 3 / CLED Medium	
10 slides	DS0120A
Plate Count Agar / MacConkey Agar No. 3	
10 slides	DS0166A
Plate Count Agar / MacConkey Agar No. 3 (with added germicide inhibitors)	
10 slides	DS0167A
Plate Count Agar / Violet Red Bile Glucose Agar	
10 slides	DS0168A
10 slides	GFD06A*
Plate Count Agar / Violet Red Bile Glucose Agar (with added germicide inhibitors)	
10 slides	DS0169A
Plate Count Agar with OGYE Agar	
10 slides	DS0170A
TTC (Red Spot)	
10 slides	DS0147A
TTC (Red Spot) / Malt Extract Agar	
10 slides	DS0155A
CLED/MacConkey / Malt Extract	
1x10 Slides	GFD02A*
CLED/MacConkey / Cetrimide	
1x10 Slides	GFD03A*
TSA/Malzextract	
1x10 Slides	GFD04A*
PCA/MacConkey with MUG	
1x10 Slides	GFD05A*
Plate count agar / Plate count Agar	
1x10 Slides	GFD07A*

Biochemical Reagents

Barbitone Complement Fixation Test Diluent

tablets for the simple and convenient preparation of the special diluent for use in complement fixation test

100 tabs	each tab for 100mL	BR0016G
----------	--------------------	---------

Basic Fuchsin Indicator

a dye for use with Endo Agar Base (CM0479) and Membrane Endo Agar LES (MM0551)

10g		BR0050A
-----	--	---------

MUG Supplement

4-methylumbelliferyl-β-D-glucuronide; a fluorogenic substrate for the detection of *Escherichia coli*

10 vials	50mg/vial	BR0071E
----------	-----------	---------

Phosphate Buffered Saline (Dulbecco A)

tablets for the preparation of a balanced salt solution without calcium and magnesium

100 tabs	each tab for 100mL	BR0014G
----------	--------------------	---------

Ringers Solution Tablets

tablets for the preparation of isotonic diluents for bacterial cells and bacteriological specimens; store in a cool, dry place

Ringers Solution		
100 tabs	each tab for 500mL	BR0052G

Saline Tablets

tablets for the preparation of 0.85% isotonic saline; store in a cool, dry place

100 tabs	each tab for 500mL	BR0053G
----------	--------------------	---------

Collection and Transportation Products

Thermo Scientific™ MicroTest™ Transport System

Media in 15mL tubes with glass beads

M4

contains gelatin, vancomycin, amphotericin B and colistin; for viruses, chlamydiae, ureaplasma, mycoplasma. Store at -25 to +8°C

12x3mL	R12502
72x3mL	R12500

M4 Female with 2 female plastic shafted swabs	
100x3mL	R12550

M4RT

contains gelatin, gentamicin and amphotericin B; for viruses and chlamydiae; store at 2–30°C

12x3mL	R12506
72x3mL	R12505

M4RT Female with 2 female plastic shafted swabs	
100x3mL	R12552

M4RT w/Micro-Tip Flocked Swab with 1 micro-tip nylon flocked swab	
100x3mL	R12566

M4RT Male/Female with 1 female plastic shafted swab and 1 male stainless steel shafted swab	
100x3mL	R12578

M5

contains vancomycin, amphotericin B, colistin and protein stabilizers; for viruses, chlamydiae, ureaplasma, mycoplasma. Store at -25 to +8°C

12x3mL	R12516
72x3mL	R12515

Other Transport Swabs

Store at -2 to 25°C

Amies Agar Gel	
10x50	TS0001A

Amies Agar Gel with Charcoal	
10x50	TS0002A

Amies Liquid	
10x50	TS0003A

Amies Agar Gel Mini Tip: Firm Applicator Wire	
10x50	TS0012A

Laboratory Preparations and Biological Extracts

Intended for use in microbiological culture media, only.

Agar Bacteriological (Agar No. 1)

a clear, purified, high gel-strength bacteriological agar compatible with all culture media; normal working concentration +1.5%w/v.

500g	LP0011B
5kg	LP0011T

Agar Technical (Agar No. 3)

a general purpose agar; normal working concentration 1.2%.

500g	LP0013B
5kg	LP0013T

Agar (Purified)

an exceptionally bright, clear agar, particularly recommended for electrophoresis studies

500g	LP0028B
------	---------

Bile Salts

a standardized ox-bile extract, used as a selective inhibitory agent in culture media

250g	LP0055J
------	---------

Bile Salts No. 3

a specifically refined fraction of bile acid salts

250g	LP0056J
------	---------

Casein Hydrolysate (Acid)

a purified hydrolysate which provides nitrogenous material for culture media

500g	LP0041B
------	---------

Dextrose Bacteriological

a sugar suitable for inclusion in bacteriological culture media

500g	LP0071B
------	---------

Gelatin Bacteriological

readily soluble in water giving a clear solution free from sulphite and preservatives

500g	LP0008B
------	---------

Haemoglobin Powder Soluble

a specially prepared haemoglobin for use in culture media

400g	LP0053C
------	---------

Lab-Lemco Powder (Beef Extract)

500g	LP0029B
5kg	LP0029T

Lactalbumin Hydrolysate

a pancreatic digest of lactalbumin which contains high levels of essential amino acids

500g	LP0048B
------	---------

Lactose Bacteriological

500g	LP0070B
------	---------

Culture Media

Laboratory Preparations and Biological Extracts, Veggietones

Liver Desiccated Bacteriological

prepared under controlled conditions to ensure maximum retention of nutritive properties

250g		LP0026J
------	--	---------

Liver Digest Neutralized

a biologically standardized papaic digest of ox liver for use as a nutrient source

500g		LP0027B
------	--	---------

Malt Extract

an extract of sprouted grain with conserved nitrogen and carbohydrate constituents; recommended for use in media for yeast and mould growth

500g		LP0039B
------	--	---------

Peptone Bacteriological

a nutritious, all-purpose peptone, for use whenever a high quality bacteriological peptone is required

500g		LP0037B
5kg		LP0037T

Peptone Bacteriological Neutralized

500g		LP0034B
------	--	---------

Peptone Mycological

a mycological peptone for isolation and diagnosis of fungi

500g		LP0040B
------	--	---------

Peptone P

a peptic digest complying with the USP specification of peptic digests of animal tissue; the high sulphur content makes it ideal for use in culture media

500g		LP0049B
------	--	---------

Peptone Special

designed to encourage the growth of the most demanding organisms

500g		LP0072B
------	--	---------

Proteose Peptone

a speciality peptone containing proteoses as defined in the USP; particularly suitable for use in toxin production or culture media for bacteria with wide ranging nutritional requirements

500g		LP0085B
------	--	---------

Purified Agar – see Agar (Purified)

Skim Milk Powder

thermophile free

500g		LP0031B
------	--	---------

Sodium Biselenite

for addition to selenite media; supplied separately to avoid teratogenic effects

100g		LP0121A
------	--	---------

Sodium Chloride Bacteriological

500g		LP0005B
------	--	---------

Soya Peptone Neutralized

an enzymatic digest of soya with a high carbohydrate content

400g		LP0044C
------	--	---------

Tryptone

a pancreatic digest of casein for use where high performance and uniformity of composition are of paramount importance

500g		LP0042B
5kg		LP0042T

Tryptone T

a tryptone developed with lower ion levels, ideal for the production of the toxin of *Clostridium tetani*

500g		LP0043B
------	--	---------

Tryptose

a mixed enzymatic hydrolysate used for the growth of fastidious organisms

500g		LP0047B
------	--	---------

Yeast Extract Powder

a dried yeast autolysate which is a good source of nitrogen and vitamins

500g		LP0021B
5kg		LP0021T

Veggietones

Thermo Scientific™ Veggietones™ Media is a range of peptones and culture media manufactured from raw materials guaranteed to be both meat free and genetically modified organism (GMO) free.

Cold Filterable Vegetable Peptone Broth

a gamma-irradiated, cold filterable Vegetable Peptone Broth suitable for microbiological media fill trials in the pharmaceutical industry, made using raw materials free from animal-derived material

500g		VG0704B
5kg		VG0704T

GMO-Free Soya Peptone

a nutritious, general purpose peptone, manufactured using GMO-free soya and fungal enzymes

400g		VG0300C
5kg		VG0300T

Vegetable Peptone No. 1

highly nutritious, general purpose peptone

500g	for 16.1 litres medium	VG0700B
5kg	for 161.0 litres medium	VG0700T

Vegetable Peptone Broth

highly nutritious, meat free alternative to Tryptone Soya Broth

500g	for 16.1 litres medium	VG0701B
5kg	for 161.0 litres medium	VG0701T

Vegetable Peptone Phosphate Broth

a meat-free alternative to Tryptose Phosphate Broth, suitable for use as a component of tissue culture media or the cultivation of fastidious organisms

500g	for 16.9 litres medium	VG0200B
------	------------------------	---------

Process Simulation

Dehydrated Media

Cold Filterable Tryptone Soya Broth

a gamma-irradiated, cold filterable Tryptone Soya Broth suitable for microbiological media fill trials in the pharmaceutical industry

500g	for 16.7 litres medium	CM1065B
2.5kg	for 83.3 litres medium	CM1065R
5kg	for 167.0 litres medium	CM1065T

Cold Filterable Vegetable Peptone Broth

a gamma-irradiated, cold filterable Vegetable Peptone Broth suitable for microbiological media fill trials in the pharmaceutical industry, made using raw materials free from animal-derived material

500g	for 16.1 litres medium	VG0704B
5kg	for 161.0 litres medium	VG0704T

BioProcess Containers (BPCs)

BPCs are intended for use in process simulations, either as a liquid placebo prior to sterile filtration, or as a growth medium for a solid placebo added downstream of processing.

Cold Filterable Tryptone Soya Broth in BPCs

ready-to-use cold filterable Tryptone Soya Broth in BPCs, which simply plug into media fill processing lines.

1L		BP1065A
10L		BP1065C
20L		BP1065E

Cold Filterable Vegetable Peptone Broth in BPCs

ready-to-use cold filterable Vegetable Peptone Broth, a vegetable based alternative to TSB, in BPCs

1L		BP0104A
10L		BP0104C
20L		BP0104E

BPC Tubing Connectors

ReadyMate™ to male Lynx™	BP0010A
ReadyMate to female Lynx	BP0020A
ReadyMate to male Kleenpak™	BP0030A
ReadyMate to female Kleenpak	BP0040A
ReadyMate to male Opta™	BP0050A
ReadyMate to female Opta	BP0060A
ReadyMate to Steam-Thru™	BP0070A
ReadyMate to 2x ReadyMate	BP0080A
ReadyMate to 5x ReadyMate	BP0090A

Environmental Monitoring

Air Sampler, Environmental Monitoring Swabs, Prepared Media, Dehydrated Media

Air Sampler

The Air Sampler comes in a carrying case, with a sampling head, NIST calibration certificate, manual, warranty and battery charger.

Aluminium head for 90mm plates	
each	EM0100A
Stainless Steel head for 90mm plates	
each	EM0100S
Aluminium head for 55mm plates	
each	EM0355A
Stainless Steel head for 55mm plates	
each	EM0355S

Air Sampler Accessories

Aluminium Sampling Head and Adaptor for 90mm plates	
each	EM0190A
Stainless Steel Sampling Head and Adaptor for 90mm plates	
each	EM0290A
Aluminium Sampling Head for 55mm plates	
each	EM0155A
Stainless Steel Sampling Head for 55mm plates	
each	EM0255A
Soft Carrying Case	
each	EM0104A
Aluminium Carrying Case	
each	EM0105A
Wall/Table Base Plate	
each	EM0109A
Floor Tripod	
each	EM0106A
Battery Charger (110V)	
each	EM0102A
Battery Charger (220V)	
each	EM0101A
Air Sampler PC Connecting Cable	
each	EM0111A
Air Sampler PC Connecting Software	
each	EM0110A
Calibration	
each	EM0100C
12 Month Extended Service Agreement (including single calibration)	
each	EM0900A
Politecnic Calibration Equipment	
each	EM0100D

Qualification Documentation

IQ	
each	EM0800A
OQ	
each	EM0800B
PQ	
each	EM0800C

Thermo Scientific™ Contact Plate Applicator

NEW Contact Plate Applicator meets ISO 18593 for surface sampling; stainless steel weight is 500g +/- 5 g and the timer is 10s +/- 0.1s.

each		EM2000A
CPA Replacement Timer Unit		
each		EM2001A

Environmental Monitoring Swabs

Ideal for use in the environmental monitoring of clean rooms, isolators and other critical areas; triple wrapped, irradiated and VHP proof.

Single Tube Swab Kits

200		TS0954C
-----	--	---------

Prepared Media

NEW triple wrapped, irradiated plates. Sterile pack with VHP indicator, according to EP/USP.

Tryptone Soya Agar		
100 plates	90mm	PO5500B
100 plates	Contact Plate (55mm)	PO5510D
Tryptone Soya Agar w/Lecithin, Polysorbate 80, Sodium Thiosuphate, L-Histidine		
100 plates	90mm	PO5501B
100 plates	Contact Plate (55mm)	PO5511D
Sabouraud Dextrose Agar		
100 plates	90mm	PO5502B*
100 plates	Contact Plate (55mm)	PO5512D*
Sabouraud Dextrose Agar w/Lecithin, Polysorbate 80, Sodium Thiosuphate, L-Histidine		
100 plates	90mm	PO5503B*
100 plates	Contact Plate (55mm)	PO5513D*

Triple-wrapped irradiated media

Tryptone Soya Agar		
10x90mm		PO0821B ^
10x90mm	VHP Proof	PO0777B
10x55mm	Contact Plate	PO0262D
10x90mm		PO5012B
10x90mm		PO5171D
Sabouraud Dextrose Agar		
10x90mm		PO0410B
10x55mm	Contact Plate	PO0394D

Dehydrated Media

R2A Agar

a medium for the enumeration of heterophilic bacteria from water samples

500g	for 27.6 litres medium	CM0906B
------	------------------------	---------

Sabouraud Dextrose Agar

acid pH medium for isolation of dermatophytes and other yeasts and moulds. This can be made selective using chloramphenicol.

500g	for 7.7 litres medium	CM0041B
2.5kg	for 38.5 litres medium	CM0041R
5kg	for 77.0 litres medium	CM0041T

Tryptone Soya Agar

a general purpose medium for growth of a wide variety of organisms; meets the formula specified in EP/USP/JP

500g	for 12.5 litres medium	CM0131B
2.5kg	for 62.5 litres medium	CM0131R
5kg	for 125.0 litres medium	CM0131T

Biochemical Identification

Biochemical identification products should be stored at 2–8°C, unless otherwise stated.

Thermo Scientific™ BactiCard™ Tests

BactiCard Candida Identification Kit

For presumptive identification of *Candida albicans*

25 tests	R21106
----------	--------

BactiCard E. coli Identification Kit

For presumptive identification of *Escherichia coli*

25 tests	R21108
----------	--------

BactiCard Neisseria Identification Kit

For presumptive identification of pathogenic *Neisseria* spp. and *Moraxella catarrhalis*

25 tests	R21110
----------	--------

BactiCard Strep Identification Kit

For presumptive identification of streptococci and related genera

25 tests	R21112
----------	--------

Thermo Scientific™ Oxoid™ Biochemical Identification System (O.B.I.S™)

O.B.I.S. campy Identification Kit

Comprising Gram lysis and L-alanyl aminopeptidase tests, to aid in the identification of Campylobacteraceae

60 tests	ID0800M
----------	---------

O.B.I.S. mono Identification Kit

For the differentiation of *Listeria monocytogenes* from other *Listeria* spp.

60 tests	ID0600M
----------	---------

O.B.I.S. PYR Identification Kit

For PYRase activity in streptococcal organisms and *Citrobacter* spp.

60 tests	ID0580M
----------	---------

O.B.I.S. Salmonella Identification Kit

For PYRase and NPA activity in organisms which are presumptive *Salmonella* spp.

60 tests	ID0570M
----------	---------

Touch Sticks

Beta-Lactamase (Nitrocefin) Identification Sticks

For the detection of beta-lactamases produced by Gram-positive bacteria. Store below –10°C.

100 sticks	BR0066A
------------	---------

Diagnostic Discs and Strips

Diagnostic discs and strips may also be stored in the freezer (–20°C). For information on disc dispensers, please see page 5.

A.L.A. Discs

for demonstration of porphyrin synthesis in *Haemophilus* spp.

25 discs	R21115
----------	--------

Thermo Scientific™ An-Indent™ Discs

For the presumptive identification of anaerobic Gram-negative bacteria

One cartridge each: Erythromycin 60µg; Rifampicin 15µg; Colistin 10µg; Penicillin 2 units; Kanamycin 1000µg; Vancomycin 5µg	
6x50 discs	DD0006A

Bacitracin Discs

For the differentiation of Lancefield Group A streptococci from other beta-haemolytic streptococci

50 discs	DD0002T
5x50 discs	DD0002B

Beta Lysin Discs*

For detection of CAMP factor in Group B streptococci

25 discs	R21120
----------	--------

C. albicans Test Kit

Reagent impregnated disc test for presumptive identification of *Candida albicans*

50 tests	R30851401
----------	-----------

Catarrhalis Test Discs

For rapid detection of butyrate esterase for presumptive identification of *Moraxella catarrhalis*

25 discs	R21121
----------	--------

Cefoxitin Discs

For the detection of methicillin resistance

5x50 discs	10µg	DD0026B
------------	------	---------

Cefpodoxime Combination Disc Kit

For the detection of extended spectrum beta-lactamases

250 tests	DD0029B
-----------	---------

Furazolidone Discs

For the differentiation of micrococci

5x50 discs	100µg	DD0028B
------------	-------	---------

Hippurate Discs

To detect the hydrolysis of sodium hippurate by Group B streptococci, *Gardnerella vaginalis* and *Campylobacter jejuni*; requires Ninhydrin Reagent - see page 62

25 discs	R21085
Ninhydrin Reagent	
25 ml	R21238

Kanamycin Discs

For the presumptive identification of Gram-negative anaerobic bacilli

5x50 discs	1000µg	DD0027B
------------	--------	---------

LAP disc w/Reagent Test

For rapid differentiation of *Aerococcus* and *Leuconostoc* from *Streptococcus*, *Enterococcus*, *Lactococcus* and *Pediococcus* using reagent impregnated discs

50 tests	R30168501
----------	-----------

Metronidazole Discs

For the presumptive identification of *Gardnerella vaginalis*

50 discs	50µg	DD0008T
----------	------	---------

*NOT CE marked

Not all products available in all territories. Please contact your local office for details.

Thermo Scientific™ Microbact™ Oxidase Strips

50 tests		MB0266A
----------	--	---------

Microdase Discs

A modified oxidase test for differentiation of staphylococci from micrococci

25 discs		R21132
----------	--	--------

Niacin Reagent Strips*

For differentiation and identification of mycobacteria

25 strips		R21090
-----------	--	--------

Nitrocefin Discs

For rapid detection of beta-lactamase production

25 strips		R211667
-----------	--	---------

O129 Discs

For the differentiation of *Vibrio* from other Gram-negative rods

50 discs	10µg	DD0014T
50 discs	150µg	DD0015T

ONPG Discs

For rapid detection of beta-galactosidase activity

50 discs		DD0013T
----------	--	---------

Optochin Discs

For the presumptive identification of *Streptococcus pneumoniae*

50 discs		DD0001T
5x50 discs		DD0001B

Pathotec CO Reagent

Impregnated paper strips for detection of cytochrome oxidase enzyme

4x25 strips		R38191
-------------	--	--------

PRO Discs w/Reagent*

For identification of yeast and anaerobic microorganisms primarily *Candida albicans* and *Clostridium difficile*

25 tests		R211357
----------	--	---------

PYR Discs w/Reagent

Detects PYRase produced by enterococci, allowing rapid differentiation of enterococci from Lancefield Group D streptococci

50 tests		R30854301
100 tests		R30854401

PYR/Esculin Discs*

For presumptive identification of Group A streptococci and enterococci; requires BactiDrop PYR Reagent - see page 68

25 tests		R21138
----------	--	--------

SPS Discs

For presumptive identification of *Peptostreptococcus anaerobius*

50 discs		DD0016T
----------	--	---------

Sulphonamide Discs

For the presumptive identification of *Gardnerella vaginalis*

50 discs	1000µg	DD0011T
----------	--------	---------

X&V Factor Discs

For the differentiation of the *Haemophilus* species

X Factor Discs		
50 discs		DD0003T
5x50 discs		DD0003B
V Factor Discs*		
50 discs		DD0004T
5x50 discs		DD0004B
X+V Factor Discs*		
50 discs		DD0005T
5x50 discs		DD0005B

* best stored at -10°C to -20°C

Cartridge Container

Complete with desiccant

each		HP0055A
------	--	---------

Thermo Scientific Microbact Kits

Microbact GNB 12A Kit

For the identification of Enterobacteriaceae; strip format

60 tests		MB1132A
----------	--	---------

Microbact GNB 12B Kit

For use with GNB 12A or 12E, to identify miscellaneous Gram-negative bacilli; strip format

60 tests		MB1133A [†]
----------	--	----------------------

Microbact GNB 24E Kit

For the identification of Enterobacteriaceae and other Gram-negative bacilli; microplate format

40 tests		MB1131A
----------	--	---------

Microbact Reagents and Accessories for GNB kits

Microbact Reagent TDA	
10mL	MB0180A
Microbact Reagent VP I	
10mL	MB0181A
Microbact Reagent VP II	
10mL	MB0184A
Microbact Reagent Nitrate A	
10mL	MB0186A
Microbact Reagent Nitrate B	
10mL	MB0187A
Microbact Reagent Indole Kovacs	
10mL	MB0209A
Microbact Reagent Set D contains the 6 reagents above	
6x10mL	MB1082A
Microbact Reagent Spot Indole DMACA	
10mL	MB1448A

Microbact Oxidase Strips	
50	MB0266A
Microbact Holding Tray 24E/12E	
each	MB5128A
Microbact Holding Tray 12A/B	
each	MB5370A

Microbact Listeria 12L Kit

20 tests	MB1128A
----------	---------

Microbact Listeria 12L Haemolysin Reagent

20 tests	MB1249A
----------	---------

Microbact Staphylococcal 12S Kit

20 tests	MB1561A
----------	---------

Microbact Staph Fast Blue Reagent Set

2x10mL	MB1588A
--------	---------

Microbact Accessories

Microbact Mineral Oil

Store at 10–30°C.

50mL	MB1093A
------	---------

Microbact Identification Package

Microbact programme for Microsoft® Windows®

each	MB1244A
------	---------

Thermo Scientific™ RapID™ Manual Identification System

RapID ANA II Panel

4hr ID of medically important anaerobes; requires Spot Indole and Inoculation Fluid-1mL

20 panels	R8311002
-----------	----------

RapID CB Plus Panel

4hr ID of *Corynebacterium* and other Gram-positive coryneform bacilli; requires Nitrate A and Nitrate B Reagents and Inoculation Fluid-2mL

20 panels	R8311008
-----------	----------

RapID NF Plus Panel

4hr ID of oxidase-positive, Gram-negative bacilli, including *Vibrio* spp.; requires Spot Indole, Nitrate A and Nitrate B Reagents and Inoculation Fluid-1mL

20 panels	R8311005
-----------	----------

RapID NH Panel

4hr ID of *Neisseria*, *Moraxella*, *Haemophilus* and related organisms; requires Spot Indole, Nitrate A and Nitrate B Reagents and Inoculation Fluid-1mL

20 panels	R8311001
-----------	----------

RapID ONE Panel

4hr ID of oxidase negative, Gram-negative bacilli; requires Spot Indole and Inoculation Fluid-2mL

20 panels	R8311006
-----------	----------

RapID SS/u Panel

2hr ID of urinary tract pathogens; requires Spot Indole and Inoculation Fluid-1mL

20 panels	R8311004
-----------	----------

RapID Staph Plus Panel

4hr ID of staphylococci and other catalase-positive, Gram-positive cocci; requires Nitrate A and Nitrate B Reagents and Inoculation Fluid-2mL

20 panels	R8311009
-----------	----------

RapID STR Panel

4hr ID of streptococci and related genera; requires Inoculation Fluid-1mL

20 panels	R8311003
-----------	----------

RapID Yeast Plus Panel

4hr ID of yeast; requires Inoculation Fluid-2mL

20 panels	R8311007
-----------	----------

RapID Accessories

Thermo Scientific™ RapID™ ERIC™ (Electronic RapID Code Compendium) Microsoft® Windows® based software for use with RapID panels; 3.5" CD in plastic case

each	R8323600
------	----------

RapID Nitrate A Reagent

15mL	R8309003
------	----------

RapID Nitrate B Reagent

15mL	R8309004
------	----------

RapID Spot Indole Reagent

15mL	R8309002
------	----------

RapID Inoculation Fluid-1mL for use with RapID ANA II, RapID NF Plus, RapID NH, RapID SS/u and RapID STR panels

20x1mL	R8325102
--------	----------

RapID Inoculation Fluid-2mL for use with RapID CB Plus, RapID ONE, RapID Staph Plus and RapID Yeast Plus panels

20x2mL	R8325106
--------	----------

For McFarland Equivalence Turbidity Standards, see page 70

Immunological Tests

Immunological tests should be stored at 2–8°C, unless otherwise stated.

Agglutinating Antisera

Agglutinating sera are packaged 2mL/vial

Bordetella Agglutinating Sera

Agglutinating sera for serological identification of *Bordetella* cultures by rapid slide agglutination test

<i>Bordetella parapertussis</i>	R30165601
<i>Bordetella pertussis</i>	R30165501

Brucella Agglutinating Sera

Agglutinating sera for serological identification of *Brucella* cultures by tube agglutination test. Also for use as positive control for stained *Brucella* suspensions.

<i>Brucella abortus</i>	R30164801
<i>Brucella melitensis</i>	R30164901

E. coli Agglutinating Sera

Agglutinating sera for serological identification of *E. coli* serotypes by slide and/or tube agglutination tests

<i>E. coli</i> Polyvalent 2 (types O26, O55, O111, O119, O126)	R30954901
<i>E. coli</i> Polyvalent 3 (types O86, O114, O125, O127, O128)	R30955001
<i>E. coli</i> Polyvalent 4 (types O44, O112, O124, O142)	R30955101

Haemophilus influenzae Agglutinating Sera

Agglutinating sera for identification of *H. influenzae* by rapid slide agglutination and CIE procedures

<i>Haemophilus influenzae</i> type a	R30166001
<i>Haemophilus influenzae</i> type b	R30166101
<i>Haemophilus influenzae</i> type c	R30166201
<i>Haemophilus influenzae</i> type d	R30166301
<i>Haemophilus influenzae</i> type e	R30166401
<i>Haemophilus influenzae</i> type f	R30166501

Neisseria meningitidis Agglutinating Sera

Agglutinating sera for identification of *N. meningitidis* by rapid slide agglutination and CIE procedures

<i>Neisseria meningitidis</i> Polyvalent Groups A-D	R30166601
<i>Neisseria meningitidis</i> Group A	R30166801
<i>Neisseria meningitidis</i> Group B	R30167501
<i>Neisseria meningitidis</i> Group C	R30166901
<i>Neisseria meningitidis</i> Group D	R30167001
<i>Neisseria meningitidis</i> Group W135	R30167401
<i>Neisseria meningitidis</i> Group X	R30167101
<i>Neisseria meningitidis</i> Group Y	R30167201
<i>Neisseria meningitidis</i> Group Z	R30167301
<i>Neisseria meningitidis</i> Group X-Z, W135	R30166701

Proteus Agglutinating Sera

Agglutinating sera for use as positive control for stained *Proteus* suspensions by tube agglutination test

<i>Proteus</i> OX2	R30165701
<i>Proteus</i> OX19	R30165801
<i>Proteus</i> OXK	R30165901

Salmonella H Agglutinating Sera

Agglutinating sera for identification of flagella H antigens of *Salmonella* strains by slide and/or tube agglutination tests

<i>Salmonella</i> H Polyvalent Phases 1 & 2	R30858501
<i>Salmonella</i> H Polyvalent Phase 2 (1, 2, 5, 6, 7)	R30163201
<i>Salmonella</i> H (E) Polyvalent for eh, enx, enz15	R30160601
<i>Salmonella</i> H (G) Polyvalent for gm, gp, fg, gq, gst, fgt, gms, gmt, gt, gmq, gpu	R30161001
<i>Salmonella</i> H (L) Polyvalent for lv, lw, lz13, lz28, lz38	R30161801
<i>Salmonella</i> H (1, 2)	R30163301
<i>Salmonella</i> H (1, 5)	R30163401
<i>Salmonella</i> H (1, 6)	R30163501
<i>Salmonella</i> H (1, 7)	R30163601
<i>Salmonella</i> H (a)	R30160201
<i>Salmonella</i> H (b)	R30160301
<i>Salmonella</i> H (c)	R30160401
<i>Salmonella</i> H (d)	R30160501
<i>Salmonella</i> H (eh)	R30160701
<i>Salmonella</i> H (enx)	R30160801
<i>Salmonella</i> H (enz15)	R30160901
<i>Salmonella</i> H (f, g)	R30161101
<i>Salmonella</i> H (gm)	R30161201
<i>Salmonella</i> H (gp)	R30161301
<i>Salmonella</i> H (gq)	R30161401
<i>Salmonella</i> H (gst)	R30161501
<i>Salmonella</i> H (i)	R30161601
<i>Salmonella</i> H (k)	R30161701
<i>Salmonella</i> H (lv)	R30161901
<i>Salmonella</i> H (lw)	R30162001
<i>Salmonella</i> H (mt)	R30162101
<i>Salmonella</i> H (r)	R30162201
<i>Salmonella</i> H (y)	R30162301
<i>Salmonella</i> H (z)	R30162401
<i>Salmonella</i> H (z4, z23)	R30162501
<i>Salmonella</i> H (z6)	R30162601
<i>Salmonella</i> H (z10)	R30162701
<i>Salmonella</i> H (z27)	R30162801
<i>Salmonella</i> H (z29)	R30162901
<i>Salmonella</i> H (z36)	R30163001
<i>Salmonella</i> H (z38)	R30163101
<i>Salmonella</i> H Rapid Diagnostic 1 (b, d, E, r)	R30959801
<i>Salmonella</i> H Rapid Diagnostic 2 (b, E, k, L)	R30959901
<i>Salmonella</i> H Rapid Diagnostic 3 (d, E, G, k)	R30160101

Salmonella O Agglutinating Sera

Agglutinating sera for identification of somatic O antigens of *Salmonella* strains by slide and/or tube agglutination tests

<i>Salmonella</i> O Polyvalent (Group A - G)	R30858101
<i>Salmonella</i> O Polyvalent (Group A - S)	R30858201
<i>Salmonella</i> O Factor 2 (Group A)	R30956701
<i>Salmonella</i> O Factors 3,10,15,19 (Group E)	R30956801
<i>Salmonella</i> O Factor 4 (Group B)	R30956901
<i>Salmonella</i> O Factor 5 (Group B)	R30957001
<i>Salmonella</i> O Factor 6,7 (Group C1)	R30957101
<i>Salmonella</i> O Factor 8 (Group C2)	R30957201
<i>Salmonella</i> O Factor 9 (Group D)	R30957301
<i>Salmonella</i> O Factor 10 (Group E1)	R30957501
<i>Salmonella</i> O Factor 11 (Group F)	R30957601
<i>Salmonella</i> O Factor 13, 22 (Group G)	R30957701
<i>Salmonella</i> O Factor 14 (Group H)	R30957801
<i>Salmonella</i> O Factor 15 (Group E2)	R30957901
<i>Salmonella</i> O Factor 16 (Group I)	R30958001
<i>Salmonella</i> O Factor 17 (Group J)	R30958101
<i>Salmonella</i> O Factor 18 (Group K)	R30958201



<i>Salmonella</i> O Factor 19 (Group E4)	R30958301
<i>Salmonella</i> O Factor 20 (Group C2)	R30958401
<i>Salmonella</i> O Factor 21 (Group L)	R30958501
<i>Salmonella</i> O Factor 27 (Group B)	R30958601
<i>Salmonella</i> O Factor 28 (Group M)	R30958701
<i>Salmonella</i> O Factor 30 (Group N)	R30958801
<i>Salmonella</i> O Factor 35 (Group O)	R30958901
<i>Salmonella</i> O Factor 38 (Group P)	R30959001
<i>Salmonella</i> O Factor 39 (Group Q)	R30959101
<i>Salmonella</i> O Factor 40 (Group R)	R30959201
<i>Salmonella</i> O Factor 41 (Group S)	R30959301
<i>Salmonella</i> O/Vi	R30957401

Shigella Agglutinating Sera

Agglutinating sera for identification of *Shigella* cultures by rapid slide agglutination tests

<i>Shigella boydii</i> Polyvalent 1 (1 to 6)	R30163901
<i>Shigella boydii</i> Polyvalent 2 (7 to 11)	R30164001
<i>Shigella boydii</i> Polyvalent 3 (12 to 15)	R30164101
<i>Shigella dysenteriae</i> Polyvalent (1-10)	R30163701
<i>Shigella flexneri</i> Polyvalent (1 to 6, X & Y)	R30163801
<i>Shigella sonnei</i> Phases 1 & 2	R30164201

Vibrio cholerae Agglutinating Sera

Agglutinating sera for identification of *Vibrio cholerae* by slide and/or tube agglutination tests

<i>Vibrio cholerae</i> Polyvalent	1	R30165001
<i>Vibrio cholerae</i> Inaba	1	R30165101
<i>Vibrio cholerae</i> Ogawa	1	R30165201

Stained Suspensions

Stained Suspensions are packaged 5mL/vial

Brucella

Stained suspensions of febrile antigen for detection, identification and quantification of specific antibodies to *Brucella* in sera

<i>Brucella abortus</i>	1	R30855501
<i>Brucella melitensis</i>	1	R30953501
Febrile Negative Control for use with all stained suspensions	2mL	R30165401

Proteus

Stained suspensions of febrile antigen for detection, identification and quantification of rickettsial antibodies in sera (somatic antigens of some strains of *Proteus* are shared with some *Rickettsia* spp.)

Proteus OX2	1	R30953601
Proteus OX19	1	R30953701
Proteus OXK	1	R30953801
Febrile Negative Control for use with all stained suspensions	2mL	R30165401



Salmonella

Stained suspensions of febrile antigen for detection, identification and quantification of specific antibodies to *Salmonella* in sera

<i>Salmonella Enteritidis</i> H (g,m)	1	R30952901
<i>Salmonella Paratyphi</i> A-H (a)	1	R30855101
<i>Salmonella Paratyphi</i> A-O (1,2,12)	1	R30953301
<i>Salmonella Paratyphi</i> B-H (b)	1	R30855201
<i>Salmonella Paratyphi</i> B-O (1,4,5,12)	1	R30855401
<i>Salmonella Paratyphi</i> C-H (c)	1	R30952801
<i>Salmonella Paratyphi</i> C-O (6,7)	1	R30953401
<i>Salmonella Typhi</i> H (d)	1	R30855001
<i>Salmonella Typhi</i> O (Group D) (9,12)	1	R30855301
<i>Salmonella Typhi</i> Vi	1	R30953901
<i>Salmonella Typhimurium</i> H (i)	1	R30953101
Febrile Negative Control for use with all stained suspensions	2mL	R30165401
Febrile Positive Control for use with <i>Salmonella</i> stained suspensions, groups a, b, c, d	2mL	R30165301

Blood and Serum Tests

Blood and serum tests should be stored at 2–8°C.

Coagulation and Serological Tests

Thermo Scientific™ Thrombo-Wellcotest™ Test Kit

Semiquantitative latex test for detection of fibrin and fibrinogen degradation products in serum and urine

20 tests		R30852601
----------	--	-----------

Thrombo-Wellcotest Collection Tubes

For use with Thrombo-Wellcotest kit

100 tubes		R30853001
-----------	--	-----------

Thrombo-Wellcotest Glycine Saline Buffer

For use with Thrombo-Wellcotest kit

100mL		R30853101
-------	--	-----------

Thrombo-Wellcotest Latex

For use with Thrombo-Wellcotest kit

3mL		R30852301
-----	--	-----------

Thymune-M Haemagglutination

A microwell test for semiquantitation of thyroid microsomal antibodies in serum

50 tests		R30850501
----------	--	-----------

Thymune-T Haemagglutination

A microwell test for semiquantitation of thyroglobulin antibodies in serum

50 tests		R30850601
----------	--	-----------

Other

Purified Phytohaemagglutinin

PHA Purified used to stimulate mitotic division of lymphocytes maintained in a cell culture and facilitate cytogenetic studies of chromosomes		
2mg		R30852801

PHA Reagent Grade used to stimulate mitotic division of lymphocytes maintained in a cell culture and facilitate cytogenetic studies of chromosomes		
45mg		R30852701

Thermo Scientific™ Oxoid™ Toxin Detection Kits

Toxin detection tests should be stored at 2–8°C.

Oxoid BCET-RPLA Toxin Detection Kit

For the detection of *Bacillus cereus* enterotoxin (diarrhoeal type) by reversed passive latex agglutination

20 tests		TD0950A
----------	--	---------

Oxoid PET-RPLA Toxin Detection Kit

For the detection of *Clostridium perfringens* enterotoxin by reversed passive latex agglutination

20 tests		TD0930A
----------	--	---------

Oxoid SET-RPLA Toxin Detection Kit

For the detection of staphylococcal enterotoxins by reversed passive latex agglutination

20 tests		TD0900A
----------	--	---------

Oxoid TST-RPLA Toxin Detection Kit

For the detection of staphylococcal toxic-shock syndrome toxin (TST) by reversed passive latex agglutination

20 tests		TD0940A
----------	--	---------

Oxoid VET-RPLA Toxin Detection Kit

For the detection of *Escherichia coli* heat-labile enterotoxin and *Vibrio cholerae* enterotoxin by reversed passive latex agglutination

20 tests		TD0920A
----------	--	---------

Oxoid VTEC-RPLA Toxin Detection Kit

For the detection of verocytotoxins VT1 and VT2 produced by *E. coli* cultured from food and faecal samples by reversed passive latex agglutination

20 tests		TD0960A
----------	--	---------



Thermo Scientific™ ProSpecT™ *C. difficile* Toxin A/B Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *Clostridium difficile* Toxin A/B in fresh, frozen or Cary-Blair preserved faecal specimens

96 wells		R244596
----------	--	---------

ProSpecT Shiga Toxin *E. coli* Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *E. coli* shiga Toxin 1 and 2 in fresh or preserved faecal specimens

96 wells		R2474096
----------	--	----------

Thermo Scientific™ Xpect™ *C. difficile* Toxin A/B Test

A rapid, chromatographic immunoassay for the direct, qualitative detection of *Clostridium difficile* Toxin A &/or B in faecal specimen

20 tests		R24650
----------	--	--------

Agglutination Tests

Latex and other agglutination tests for identification and differentiation of micro-organisms. Thermo Scientific™ Dryspot™ products are stored at 2–25°C; other agglutination tests should be stored at 2–8°C.

Bacterial Antigen Test

Thermo Scientific™ Wellcogen™ Bacterial Antigen Kit

A rapid latex agglutination test for direct detection of antigens in body fluids: Group B streptococci; *H. influenzae* type b; *S. pneumoniae*; *N. meningitidis* groups

30 tests		R30859602
----------	--	-----------

Wellcogen Bacterial Antigen Kit Components and Accessories

Wellcogen Haemophilus influenzae b Test Kit
for direct detection of Haemophilus influenzae type b antigens present in body fluids or blood cultures

30 tests		R30858801
----------	--	-----------

Wellcogen Neisseria meningitidis A,C,Y,W135 Test Kit
for direct detection of Neisseria meningitidis groups A, C, Y, W 135 antigens present in body fluids or blood cultures

30 tests		R30859203
----------	--	-----------

Wellcogen Neisseria meningitidis B/E. coli K1
for direct detection of Neisseria meningitidis group B antigens or E. coli K1 antigens present in body fluids, blood cultures or on culture media

30 tests		R30859502
----------	--	-----------

Wellcogen Strep B Rapid Latex Agglutination Test
for direct detection of antigens to Group B streptococci present in body fluids or blood cultures

30 tests		R30858701
----------	--	-----------

Wellcogen Streptococcus pneumoniae Test Kit
for direct detection of Streptococcus pneumoniae capsular antigens present in body fluids or blood cultures

30 tests		R30859001
----------	--	-----------

Streptex/Wellcogen Disposable Cards

50 slides		R30368601
-----------	--	-----------

Wellcogen EDTA Solution (0.1M)
for treatment of serum samples to be tested using Wellcogen kits

10mL		R30164501
------	--	-----------

Campylobacter

Campylobacter Test

For the identification of enteropathogenic *Campylobacter* spp. from selective media

50 tests		DR0155M
----------	--	---------

Clostridium difficile

Oxoid *C. difficile* Test Kit

A rapid test for the identification of *Clostridium difficile*

50 tests	DR1107A
----------	---------

Cryptococcus

Cryptococcus Test Kit

For quantitative or semiquantitative detection of polysaccharide antigens associated with *Cryptococcus neoformans* infection in serum or CSF

50 tests	R30851501
----------	-----------

Enterohaemorrhagic *E. coli*

Thermo Scientific™ DrySpot™ *E. coli* O157

For the identification of *E. coli* serogroup O157

120 tests	DR0120M
-----------	---------

DrySpot Serocheck

Individual agglutinating latexes to screen for *E. coli* serogroups O26, O91, O103, O111, O128 and O145

DrySpot Serocheck O26	
15 tests	DR0310M†
DrySpot Serocheck O91	
15 tests	DR0320M†
DrySpot Serocheck O103	
15 tests	DR0330M†
DrySpot Serocheck O111	
15 tests	DR0340M†
DrySpot Serocheck O128	
15 tests	DR0350M†
DrySpot Serocheck O145	
15 tests	DR0360M†

DrySpot Seroscreen

For identifying *E. coli* serogroups O26, O91, O103, O111, O128 and O145 using a single test reagent

60 tests	DR0300M
----------	---------

Oxoid *Escherichia coli* O157 Latex Test

For the identification of *E. coli* serogroup O157

100 tests	DR0620M
-----------	---------

Wellcolex *E. coli* O157 Kit

For the identification of *E. coli* serogroup O157

50 tests	R30959501
----------	-----------

RIM™ *E. coli* O157:H7 Kit

For the identification of *E. coli* serotypes O157 and H7

50 tests	R24250
----------	--------

Wellcolex Accessories

Wellcolex Disposable Cards

50 cards	R30368701
----------	-----------

Haemophilus influenzae

See Wellcogen Bacterial Antigen Kit

Infectious mononucleosis

Oxoid Infectious Mononucleosis Test

For detection of infectious mononucleosis heterophile antibodies in serum and in plasma

100 tests	DR0780M
-----------	---------

Legionella

DrySpot Legionella Kits

Latex agglutination kits for the differentiation of predominant Legionella serogroups

DrySpot Legionella Pneumophila Serogroup 1 Kit	
60 tests	DR0200M
DrySpot Legionella Pneumophila Serogroups 2–14 Kit	
60 tests	DR0210M
DrySpot Legionella Species Kit	
60 tests	DR0220M

Oxoid Legionella Latex Test

For the differentiation of predominant Legionella serogroups grown on selective media

50 tests	DR0800M
----------	---------

Supplementary Reagents for the Oxoid Legionella Latex Test

Legionella Pneumophila Serogroup 1	
50 tests	DR0801M
Legionella Pneumophila Serogroups 2–14	
50 tests	DR0802M
Legionella Species Test Reagent	
50 tests	DR0803M

Listeria

Oxoid Listeria Test Kit

A rapid latex agglutination test for the presumptive identification of *Listeria* spp.

100 tests	DR1126A
-----------	---------

Meningitis

See Wellcogen Bacterial Antigen Kit

Salmonella

Oxoid Salmonella Test Kit

A rapid latex agglutination test for the presumptive identification of *Salmonella* spp.

100 tests	DR1108A
-----------	---------

Wellcolex Colour Salmonella Kit

A rapid latex agglutination test for detection and presumptive serogrouping of Salmonella from Selenite Broth and/or solid culture media

50 tests	R30858301
200 tests	R30858302

Shigella

Wellcolex Colour Shigella Kit

A rapid latex agglutination test for detection and species identification of Shigella from solid culture media

50 tests	R30858401
----------	-----------

†Product is not a stock item, minimum order and lead times may apply. Please contact your Customer Services Administrator for details.

Not all products available in all territories. Please contact your local office for details.

Staphylococcus and MRSA

Thermo Scientific™ BactiStaph™ Tests

For detection of clumping factor and protein A in identification of *Staphylococcus aureus*

150 tests		R21144
450 tests		R21143

Coagulase Plasma (lyophilized)

For detection of coagulase enzyme in staphylococci

5mL		R21050
15mL		R21051
25mL		R21052
6x5mL		R21060

DrySpot Staphytest Plus

A latex kit for identification of *S. aureus*

120 tests		DR0100M
-----------	--	---------

Oxid PBP2' Test

A latex agglutination test for the detection of penicillin binding protein 2' (PBP2') in isolates of *S. aureus* as an aid in identifying MRSA

50 tests		DR0900A
----------	--	---------

Thermo Scientific™ Staphaurex™ Latex Agglutination Test

White latex particle coated with human fibrinogen, for detecting clumping factor, and specific IgG, for detection of protein A, for identification of *S. aureus*

120 tests		R30859901
400 tests		R30859902

Staphaurex Plus

Yellow latex particle coated with human fibrinogen, for detecting clumping factor, and specific IgG, for detection of protein A, and surface antigens characteristic of MRSA and MSSA

150 tests		R30950102
450 tests		R30950201

Staphaurex Accessories

Staphaurex Disposable Cards for use with Staphaurex or Staphaurex Plus		
50 cards		R30369001

Staphylase Test

A rapid test kit for the identification of *S. aureus*

100 tests		DR0595A
-----------	--	---------

Thermo Scientific™ Staphytest™ Plus Test

A rapid latex test for the identification of *S. aureus*

100 tests		DR0850M
500 tests		DR0850B

Streptococcus Group B

See Wellcogen Bacterial Antigen Kit



Streptococcus

Oxid Streptococcal Grouping Kit

A rapid latex agglutination test for the identification of streptococcal groups A, B, C, D, F and G

50 tests		DR0585A
----------	--	---------

Oxid Streptococcal Grouping Kit: Supplementary Reagents

Extraction Enzyme		
2 vials		DR0593G
Streptococcal Grouping Reagent A		
50 tests		DR0586G
Streptococcal Grouping Reagent B		
50 tests		DR0587G
Streptococcal Grouping Reagent C		
50 tests		DR0588G
Streptococcal Grouping Reagent D		
50 tests		DR0589G
Streptococcal Grouping Reagent F		
50 tests		DR0590G
Streptococcal Grouping Reagent G		
50 tests		DR0591G
Polyvalent Positive Control		
50 tests		DR0592G
Disposable Reaction Cards		
50 cards		DR0500G

Thermo Scientific™ PathoDx™ Strep Grouping Kit

A latex agglutination assay for the identification of β -haemolytic streptococci of Lancefield groups A, B, C, F, and G from primary culture plates or pure culture from broth

60 tests		R62025
----------	--	--------

PathoDx Strep Grouping Accessories

PathoDx Strep Grouping Universal Kit for use with PathoDx Strep A, B, C, F and G Group Latexes		
60 tests		R62076
PathoDx Strep Grouping Positive Control		
3mL		R24042
PathoDx Strep A Group Latex		
3mL		R62030
PathoDx Strep B Group Latex		
3mL		R62031
PathoDx Strep C Group Latex		
3mL		R62032
PathoDx Strep F Group Latex		
3mL		R62034
PathoDx Strep G Group Latex		
3mL		R62035

PathoDx Extraction Reagents

For use with PathoDx Strep Grouping Kit

PathoDx Extraction Reagent 1		
7mL		R62050
14mL		R62052
PathoDx Extraction Reagent 2		
7mL		R62055
14mL		R62057
PathoDx Extraction Reagent 3		
14mL		R62060

PathoDxtra Strep Grouping Kit

A rapid latex agglutination test for identification of B-haemolytic streptococci of Lancefield Groups A, B, C, D, F and G from primary culture plates

PathoDxtra Strep Grouping Kit	
60 tests	DR0700M
PathoDxtra Strep Group A Latex	
60 tests	DR0701G
PathoDxtra Strep Group B Latex	
60 tests	DR0702G
PathoDxtra Strep Group C Latex	
60 tests	DR0703G
PathoDxtra Strep Group D Latex	
60 tests	DR0704G
PathoDxtra Strep Group F Latex	
60 tests	DR0705G
PathoDxtra Strep Group G Latex	
60 tests	DR0706G
PathoDxtra Polyvalent Positive Control	
2.8mL	DR0707G
PathoDxtra Extraction Reagents 1 x reagent 1 (4mL), 1 x reagent 2 (4mL), 2 x reagent 3 (10mL) for use with PathoDxtra Strep Grouping Kit	
	DR0709M
PathoDxtra Strep Grouping Reagent Set 1 x reagent 1 (4mL), 1 x reagent 2 (4mL), 1 x reagent 3 (10mL), Polyvalent positive control, mixing sticks and reaction cards	
60 tests	DR0710M
PathoDxtra Reaction Cards	
50 cards	DR0720G

Thermo Scientific™ Streptex™ Latex Agglutination Test

For identification of Lancefield Groups A, B, C, D, F and G streptococci, using enzyme extraction

50 tests	R30950501
200 tests	R30164701

Streptex Rapid Latex Agglutination Test

For identification of Lancefield Groups A, B, C, F and G, using acid extraction

50 tests	R30950555
----------	-----------

Streptex Accessories

Streptex/Wellcogen Disposable Cards	
50 slides	R30368601
Streptex Extraction Enzyme	
25 tests	R30951001
Streptex Acid Extraction Kit. Extraction reagent 1, Extraction reagent 2, Extraction reagent 3, Reaction tubes, pipettes and mixing sticks	
50 tests	R30951301
Streptex Latex Group A	
50 tests	R30950601
Streptex Latex Group B	
50 tests	R30950701
Streptex Latex Group C	
50 tests	R30950801
Streptex Latex Group D	
50 tests	R30950901
Streptex Latex Group F	
50 tests	R30951101
Streptex Latex Group G	
50 tests	R30951201
Streptex Positive Control Polyvalent Control	
2mL	R30164601

Streptococcus pneumoniae

See also Wellcogen Bacterial Antigen Testing

DrySpot Pneumo Kit

A latex agglutination assay for confirmation of *Streptococcus pneumoniae* from culture plates or positive blood cultures

60 tests	DR0420M
----------	---------

Treponema pallidum

Oxoid TPHA Test

A sensitive and specific passive haemagglutination test for the detection of antibody to *Treponema pallidum*

200 tests	DR0530M
-----------	---------

Enzyme Immunoassays

Enzyme immunoassays should be stored at 2–8°C.

Gastroenteritis: Bacterial

Thermo Scientific™ IDEIA™ Amplified Hp StAR™ Assay

For detection of *Helicobacter pylori* antigen in human stool samples; may be used to diagnose infection, monitor treatment or confirm eradication

96 wells	K663011-2
----------	-----------

ProSpecT Campylobacter Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *Campylobacter* specific antigen in fresh or preserved faecal specimens

96 wells	R2476096
----------	----------

ProSpecT C. difficile Toxin A/B Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *Clostridium difficile* Toxin A/B in fresh, frozen or Cary-Blair preserved faecal specimens

96 wells	R244596
----------	---------

ProSpecT Shiga Toxin E. coli Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *E. coli* shiga Toxin 1 and 2 in fresh or preserved faecal specimens

96 wells	R2474096
----------	----------

Gastroenteritis: Viral

IDEIA Norovirus

A qualitative EIA for direct detection of Norovirus antigen in human faecal specimens

96 wells	K604411-2
----------	-----------

ProSpecT Adenovirus Microplate Assay

A direct, qualitative enzyme immunoassay for detection of Adenovirus specific antigen in faecal specimens

96 wells	R240096
----------	---------

ProSpecT Astrovirus Microplate Assay

A unique, amplified enzyme immunoassay for detection of Astrovirus in faecal specimens

96 wells	R240196
----------	---------

ProSpecT Rotavirus Microplate Assay

A direct enzyme immunoassay for detection of Rotavirus in faecal specimens

96 wells	R240396
----------	---------

Lyme Disease

IDEIA *Borrelia burgdorferi*, IgG

For detection of human serum IgG antibodies to *Borrelia burgdorferi*

96 wells	K602911-2
----------	-----------

IDEIA *Borrelia burgdorferi*, IgM

For detection of human serum IgM antibodies to *Borrelia burgdorferi*

96 wells	K603011-2
----------	-----------

IDEIA Lyme Neuroborreliosis

For detection of intrathecally produced human IgG and IgM antibodies to *Borrelia burgdorferi* in paired CSF and serum samples

96 wells	K602811-2
----------	-----------

Parasitology

ProSpecT *Cryptosporidium* Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *Cryptosporidium* specific antigens in fresh or preserved faecal specimens

96 wells	R2454096
----------	----------

Sample Panel
for use with ProSpecT *Cryptosporidium* Microplate Assay; contains 5 positive and 5 negative vials

10 vials	R244803
----------	---------

ProSpecT *Entamoeba histolytica* Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *E. histolytica* specific antigen in fresh faecal specimens

96 wells	R2456096
----------	----------

ProSpecT *Giardia*/*Cryptosporidium* Microplate Assay

A direct, qualitative enzyme immunoassay for simultaneous detection of *Giardia* and *Cryptosporidium* specific antigen in fresh or preserved faecal specimens

96 wells	R2458496
----------	----------

ProSpecT *Giardia* Controls
for use with ProSpecT *Giardia* Microplate Assays; contains 2 positive and 2 negative vials

4 vials	R24280
---------	--------

ProSpecT *Giardia* EZ Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *Giardia* specific antigen in fresh or preserved faecal specimens

96 wells	R2458596
----------	----------

ProSpecT *Giardia* Controls
for use with ProSpecT *Giardia* Microplate Assays; contains 2 positive and 2 negative vials

4 vials	R24280
---------	--------

ProSpecT *Giardia* Microplate Assay

A direct, qualitative enzyme immunoassay for detection of *Giardia* specific antigen in fresh or preserved faecal specimens

96 wells	R2458096
----------	----------

ProSpecT *Giardia* Controls
for use with ProSpecT *Giardia* Microplate Assays; contains 2 positive and 2 negative vials

4 vials	R24280
---------	--------

Sexually Transmitted Diseases

IDEIA *PCE* Chlamydia

A highly sensitive and specific dual amplification EIA for detection of genus specific *Chlamydia* lipopolysaccharides

192 wells	K603211-2
-----------	-----------

IDEIA Chlamydia Accessories

IDEIA Chlamydia Blocking Reagents
a selective blocking test for confirmation of IDEIA *PCE* Chlamydia

40 tests	S604130-2
----------	-----------

IDEIA Chlamydia Specimen Collection Kit, Male/Female

50 tests	S600730-2
----------	-----------

IDEIA Chlamydia Transport Medium concentrated (x10) transport medium for use with IDEIA *PCE* Chlamydia

25mL	S603830-2
------	-----------

Veterinary Use

Mycoplasma hyopneumoniae ELISA Kit

For detection of antibodies to *Mycoplasma hyopneumoniae* in porcine serum or colostrum

96 wells	K004311-9
960 wells	K004321-9

Pasteurella multocida Toxin ELISA Kit

A double antibody sandwich EIA for detection of toxigenic *Pasteurella multocida* in extracts of primary mixed bacterial cultures of nasal swabs taken from pigs

96 wells	K000911-9
960 wells	K000921-9

Accessory Reagents for PMT Serology reagents required to convert *Pasteurella multocida* Toxin ELISA from a bacterial assay to a serological assay

96 tests	K003811-9
----------	-----------

ELISA Accessories

IDEIA Wash Buffer

A x25 concentrated buffer solution for use with Thermo Scientific™ IDEIA™ Amplified Hp StAR™ (K663011-1), IDEIA Norovirus (K604411-2), IDEIA Lyme Neuroborreliosis (K602811-2), IDEIA Borrelia burgdorferi IgG (K602911-2), IDEIA Borrelia burgdorferi IgM (K603011-2), Pasteurella multocida ELISA Kit (K000911-9) and Mycoplasma hyopneumoniae ELISA kit (K004311-9)

50mL	S612630-2
------	-----------

ProSpecT Microplate Wash Buffer

A x10 concentrated buffer solution for use in the wash procedure of the Thermo Scientific™ ProSpecT™ Microplate Assay range except ProSpecT Shiga Toxin E.coli (STEC) (R2474096)

1L	R2450001
----	----------

ProSpecT Microplate Wash Buffer

A x10 concentrated buffer solution for use in the wash procedure of ProSpecT Shiga Toxin E.coli (STEC) Microplate Assays (R2474096)

1L	R24500011
110mL	R2450002

Food Allergen ELISA Assays

ELISA Systems™ kits are a comprehensive range of microwell ELISA assays for detecting and quantifying allergen proteins and residues in food and swab samples.

ELISA Assay Kits

Each kit contains 48 antibody-coated microwells, extraction solution, wash solution, conjugate, substrate, stop solution, controls, microwell holder, certificate of analysis and instructions.

Almond Residue Detection Kit (ELISA Systems code ESARD-48)	
48 wells	ESARD-48
Buckwheat Residue Detection Kit (ELISA Systems code ESBWPRD-48)	
48 wells	ESBWPRD
Casein Residue Detection Kit (ELISA Systems code ESCASPRD-48)	
48 wells	ESCASPRD
Crustacean Residue Detection Kit (ELISA Systems code ESCRURD-48)	
48 wells	ESCRURD
Egg Residue Detection Kit (ELISA Systems code ESERD-48)	
48 wells	ESERD-48
Enhanced Egg Residue Detection Kit (ELISA Systems code ESEGG-48)	
48 wells	ESEGG-48
Gliadin Residue Detection Kit (ELISA Systems code ESGLISS-48)	
48 wells	ESGLISS-48
Hazelnut Residue Detection Kit (ELISA Systems code ESHRD-48)	
48 wells	ESHRD-48
Lupin Residue Detection Kit (ELISA Systems code ESLFP-48)	
48 wells	ESLFP-48
Milk (β-Lactoglobulin) Residue Detection Kit (ELISA Systems code ESMRDLG-48)	
48 wells	ESMRDLG
Mustard Seed Protein Residue Detection Kit (ELISA Systems code ESMUS-48)	
48 wells	ESMUS
Peanut Residue (Sensitive) Detection Kit (ELISA Systems™ code ESPRDS-48)	
48 wells	ESPRDS-48
Sesame Seed Protein residue Kit (ELISA Systems code ESSESE-48)	
48 wells	ESSESE-48

Soy Protein Residue Detection Kit (ELISA Systems code ESSOYPRD-48)	
48 wells	ESSOYPRD
Total Milk Residue Assay (ELISA Systems code ESTMLK-48)	
48 tests	ESTMLK-48

Swab Sampling Kit

Use with ELISA Systems allergen assays for environmental sampling. Each kit contains 10 sterile swabs, 10 vials of sterile swab wetting solution and instructions.

Food Allergen Swab Kit (ELISA Systems code ESSWB-10)	
10 swabs	ESSWB-10

Allergen Testing Sundry Items

Top-up concentrates for use with ELISA Systems allergen assays

Extraction Solution (Standard) 20x concentrate (ELISA Systems code ESEXT)		
25mL	20x conc.	ESEXT
Extraction Solution (Polyphenol) 20x concentrate (ELISA Systems code ESADDSOL)		
25mL	20x conc.	ESADDSOL
Wash Solution 20x concentrate (ELISA Systems code ESWASH)		
25mL	20x conc.	ESWASH

Immunofluorescence Assays



Immunofluorescence assays should be stored at 2–8°C.

Respiratory Viral Infections

Thermo Scientific™ IMAGEN™ Adenovirus Test

A qualitative immunofluorescence test for direct detection of Adenovirus in respiratory secretions and ophthalmic specimens, as well as detection and confirmation of Adenovirus in cell culture

50 tests	K610011-2
IMAGEN Adenovirus Positive Control Slide	
each	S611330-2

IMAGEN Human Metapneumovirus (hMPV)

A qualitative immunofluorescence test for direct detection of hMPV in respiratory specimens

50 tests	K612511-2
IMAGEN hMPV Positive Control Slide	
each	S612830-2

IMAGEN Influenza Virus A and B

A qualitative immunofluorescence test for direct detection and differentiation of Influenza A virus and Influenza B virus in respiratory specimens and cell culture

50 tests	K610511-2
IMAGEN Influenza Virus A and B Positive Control Slide	
each	S611230-2

IMAGEN Parainfluenza Virus Group

A qualitative immunofluorescence test for direct detection and confirmation of Parainfluenza virus antigens in respiratory specimens and cell culture

50 tests	K610311-2
IMAGEN Parainfluenza Virus Positive Control Slide	
each	S611130-2

IMAGEN Parainfluenza Virus Typing

A qualitative immunofluorescence test for direct detection and differentiation of Parainfluenza virus types 1, 2 and 3 antigens in respiratory specimens and cell culture

50 tests	K610411-2
IMAGEN Parainfluenza Virus Positive Control Slide	
each	S611130-2

IMAGEN Respiratory Screen

A qualitative, indirect immunofluorescence test for presumptive detection of Respiratory Syncytial Virus, Influenza A and B virus, Parainfluenza virus types 1, 2 and 3 and Adenovirus

100 tests	K612011-2
IMAGEN Respiratory Screen Positive Control Slide	
each	S612130-2

IMAGEN Respiratory Syncytial Virus (RSV)

A qualitative immunofluorescence test for direct detection of RSV in respiratory specimens and cell culture

50 tests	K610211-2
IMAGEN Respiratory Syncytial Virus Positive and Negative Control Slides	
each	S610830-2

IMAGEN Teflon™-Coated Slides

100 slides	S611430-6
------------	-----------

Thermo Scientific™ PathoDx™ RSV Kit

An immunofluorescent test for detecting RSV

100 tests	R62370
PathoDx RSV Control Slides for use with PathoDx RSV Kit	
5 slides	R62375

PathoDx Respiratory Virus Panel

An immunofluorescent assay for detection of Influenza A and B, Adenovirus, Parainfluenza 1, 2, and 3 and Respiratory Syncytial Virus (RSV)

200 Tests	R62400
PathoDx Respiratory Virus Panel Control Slides	
5 slides	R62414
PathoDx Adenovirus Reagent	
2mL	R62407
PathoDx Parainfluenza 1 Reagent	
2mL	R62408
PathoDx Parainfluenza 2 Reagent	
2mL	R62409
PathoDx Parainfluenza 3 Reagent	
2mL	R62410
PathoDx RSV Reagent	
2mL	R62411
PathoDx Influenza A Reagent	
2mL	R62405
PathoDx Influenza B Reagent	
2mL	R62406
PathoDx Respiratory Virus Panel Negative Control Reagent	
10mL	R62412
PathoDx Respiratory Virus Panel Screening Reagent	
10mL	R62416
PathoDx Respiratory Virus Panel Mounting Fluid	
15mL	R62413

Sexually Transmitted Diseases

IMAGEN Chlamydia Kit

A qualitative, direct immunofluorescence test for detection of *Chlamydia* antigen

50 tests	K610111-2
IMAGEN Chlamydia Collection Pack	
20 kits	S610730-2
IMAGEN Chlamydia Positive Control Slide	
each	S610930-2

PathoDx Chlamydia Culture Confirmation Kit

An immunofluorescence test specific for the detection and identification of all *Chlamydia* spp. from culture

100 tests	R62210
-----------	--------

PathoDx Chlamydia Trachomatis FA Direct Test Kit

An immunofluorescence test specific for the detection of *Chlamydia trachomatis* from urethral and endocervical specimens. 15 minute test

100 tests	R62220
PathoDx Chlamydia Control slide for use with <i>Chlamydia</i> direct immunofluorescent assays. Each slide has 2 wells: 1 positive and 1 negative	
5 slides	R62205

IMAGEN Herpes Simplex Virus (HSV) Kit

A qualitative, direct immunofluorescence test for identification and differentiation of HSV1 and HSV2 in cell culture

50 tests	K610611-2
IMAGEN HSV Positive Control Slide	
each	S611030-2

PathoDx Herpes Typing Kit

An immunofluorescent assay for the detection and typing of Herpes Simplex Virus Types I and II

100 tests	R62250
PathoDx Herpes Typing Control Slides	
5 Slides	R62251

Lateral Flow Assays

Lateral flow assays should be stored at 2–8°C, unless otherwise stated.

RAPID Hp StAR

An amplified immunochromatographic assay for direct detection of *Helicobacter pylori* antigen in faecal samples, to diagnose infection or confirm eradication

20 tests	K671411-2
----------	-----------

Xpect C. difficile Toxin A/B Test

A rapid, chromatographic immunoassay for the direct, qualitative detection of *Clostridium difficile* Toxin A &/or B in faecal specimen

20 tests	R24650
----------	--------

Xpect Cryptosporidium Test

A rapid, chromatographic immunoassay for the direct, qualitative detection of *Cryptosporidium* antigens in fresh, frozen or preserved faecal specimens

20 tests	R2451020
Xpect Cryptosporidium Control for the Xpect Cryptosporidium or Xpect Giardia/Cryptosporidium Test Kits; contains 1 positive and 1 negative vial	
2 vials	R2451030

Xpect Flu A & B Test

A rapid, chromatographic immunoassay for the direct, qualitative detection of Influenza A and Influenza B viral antigens in respiratory specimens

20 tests	R24600
Xpect Flu A & B Control Swabs includes 20 Flu A swabs and 20 Flu B swabs	
20 tests	R246003

Xpect Giardia Test

A rapid, chromatographic immunoassay for the direct, qualitative detection of *Giardia* antigens in fresh, frozen or preserved faecal specimens

20 tests	R2450020
----------	----------

Xpect Giardia/Cryptosporidium Test

A rapid chromatographic immunoassay for the direct, qualitative detection of *Giardia* and *Cryptosporidium* antigens in fresh, frozen or preserved faecal specimens

20 tests	R2450520
Xpect Giardia/Cryptosporidium Control Kit for the Xpect Cryptosporidium, Xpect Giardia or Xpect Giardia/Cryptosporidium Test Kits; contains 1 Giardia positive, 1 Cryptosporidium positive and 1 negative vial	
3 vials	R2450530

Xpect Legionella Test

A rapid, chromatographic immunoassay for the direct detection of *L. pneumophila* antigen in urine samples; store at 2–30°C

20 tests	R24680
----------	--------

Xpect Rotavirus Test

A rapid chromatographic immunoassay for qualitative detection of rotavirus antigen in faecal samples

20 tests	R24655
Xpect Rotavirus Control Kit contains 1 negative and 1 positive control	
2 vials	R246551

Xpect RSV Test

A rapid, chromatographic immunoassay for the direct, qualitative detection of Respiratory Syncytial Virus (RSV) directly from nasopharyngeal specimens in neonatal and paediatric patients

30 tests	R24601
Xpect RSV Control Kit contains 1 negative and 1 positive control	
2 vials	R246012



Reagents and Stains

Reagents

Thermo Scientific™ BactiDrop™ Reagents

Thermo Scientific BactiDrop Reagents and Stains are presented in crushable ampoules encased in protective plastic sheaths with convenient dropper dispenser tips; store at 20–25°C.

BactiDrop Desoxycholate (10%) for bile solubility reaction	
50 amp.	R21508
BactiDrop Ferric Chloride (10%) for detection of phenylalanine deamination	
50 amp.	R21514
BactiDrop Indole Kovacs for detection of the ability to split indole from the tryptophan molecule; for use with Indole broth	
50 amp.	R21522
BactiDrop Indol, Ehrlichs	
50 amp.	R21520
BactiDrop Nitrate A to test for the ability to reduce nitrate to nitrite or free nitrogen gas, in conjunction with Nitrate Reagent B	
50 amp.	R21536
BactiDrop Nitrate B to test for the ability to reduce nitrate to nitrite or free nitrogen gas, in conjunction with Nitrate Reagent A	
50 amp.	R21538
BactiDrop Oxidase to test for the presence of the oxidase enzyme	
50 amp.	R21540
BactiDrop Potassium Hydroxide (10%) for preparation of clinical specimens for microscopic examination of fungal elements	
50 amp.	R21524
BactiDrop PYR for rapid, presumptive identification of Group B streptococci and enterococci; use with PYR/Esculin disc* - see page 57	
50 amp.	R21544
BactiDrop Spot Indole (DMAC) to test for the ability to split indole from the tryptophan molecule, by the spot test method	
50 amp.	R21550

Microbact Reagents

Store at 2–8°C.

Microbact Reagent Indole Kovacs	
10mL	MB0209A
Microbact Reagent Spot Indole DMACA	
10mL	MB1448A

Other Reagents

Store at 2–8°C; unless otherwise stated.

Ninhydrin Reagent detects sodium hippurate hydrolysis; use with Hippurate disc - see page 56	
25mL	R21238
Nitrate Reagent A for AFB* use in conjunction with Nitrate Reagent B for AFB	
25mL	R21243
Nitrate Reagent B for AFB* use in conjunction with Nitrate Reagent A for AFB	
25mL	R21244
Spot Indole Reagent* p-dimethylaminocinnamaldehyde (DMACA); for detection of the ability to split indole from the tryptophan molecule; store at 2–30°C	
25mL	R21245

Stains

BactiDrop Stains

Store ampoules at 20–25°C; unless otherwise stated.

BactiDrop Acridine Orange a fluorescent stain for detecting bacteria in body fluids	
50 amp.	R21502
BactiDrop Calcofluor White a fluorescent stain for fungi, <i>Microsporidium</i> , <i>Acanthamoeba</i> and <i>Pneumocystis</i> ; store at 15–30°C	
50 amp.	R21507
BactiDrop Dobell & O'Connor Iodine for staining protozoa	
50 amp.	R21510
BactiDrop India Ink for staining encapsulated microorganisms	
50 amp.	R21518
BactiDrop Lactophenol Aniline Blue for staining fungal elements	
50 amp.	R21526
BactiDrop Lugol's Iodine for trichrome stain procedure	
50 amp.	R21528

Other Stains

Gram stains are stored at 15–30°C; other stains at 20–25°C.

Gram Stain Kit includes 1 each of crystal violet, decolorizer, iodine, safranin and a tray	
4x250mL	R40080
Gram Crystal Violet	
250mL	R40052
5x250mL	R40053
3.78 litres	R40073
Gram Decolorizer	
250mL	R40054
5x250mL	R40055
3.78 litres	R40075
Gram Iodine	
250mL	R40056
5x250mL	R40057
3.78 litres	R40077
Gram Safranin	
250mL	R40058
5x250mL	R40059
3.78 litres	R40079
Gram Stain Kit Tray holds 4x250mL bottles	
each	R40081
TB Auramine-Rhodamine T* a fluorescent stain for microscopic detection of AFB	
250mL	R40090
5x250mL	R40190
TB Decolorizer, Truant-Moore* (5% Acid Alcohol) for TB Auramine O and TB Auramine-Rhodamine T fluorescent stains	
250mL	R40107
5x250mL	R40207
TB Kinyoun AFB Stain Kit* includes 2 Kinyoun, 2 Decolorizer and 2 Methylene Blue	
6x250mL	R40112
TB Potassium Permanganate * a counterstain for AFB Fluorescent Stain	
250mL	R40092
5x250mL	R40192
Trichrome (Wheatley) * for detection and differentiation of intestinal parasites	
250mL	R40025

Stain Controls

Store stain controls at 20–25°C.

QC-Slide AFB Stain Control* heat fixed quality control slide for acid-fast bacilli stain procedures	
10 slides	R40144
45 slides	R40146
QC-Slide Gram Stain Control* heat fixed quality control slide for Gram stain procedures	
10 slides	R40140
45 slides	R40142

Laboratory Supplies

Microbact Mineral Oil Store at 10–30°C.	
50mL	MB1093A
Polysorbate 80 Hydrolysis Substrate Concentrate Store at 2–8°C.	
5mL	R21276

Equivalence Turbidity Standards

Store at 20–25°C.

McFarland Equivalence Turbidity Standard Set*

1 each of 0.5, 1.0, 2.0, 3.0 and 4.0, with a visual comparison card to aid interpretation; tube dimensions 15x103mm	
5 tubes	R20421

Individual McFarland Equivalence Turbidity Standard Tubes*

McFarland Equivalence Turbidity Standard 0.5*	
each	R20410
McFarland Equivalence Turbidity Standard 1.0*	
each	R20411
McFarland Equivalence Turbidity Standard 2.0*	
each	R20412
McFarland Equivalence Turbidity Standard 3.0*	
each	R20413
McFarland Equivalence Turbidity Standard 4.0*	
each	R20414
McFarland Equivalence Turbidity Standard 5.0*	
each	R20415

Sample Prep for Molecular and non-molecular workflows

Prepare even the most challenging samples with options to pool, clean, purify and automate your sample preparation for molecular detection and other techniques.

Pathatrix

Applied Biosystems™ Pathatrix™ Auto System & Reagents

Item Number	Description
PTXAUTOINS	Pathatrix Auto Instrument
4480724	Lysis Buffer, FS
APE500SDP	Pathatrix 10 Pooling E. coli O157:H7 Kit
POOL1010MLN	Pathatrix 10 Pooling Kit
APS500P	Pathatrix 10 Pooling Salmonella spp. Kit
APDES250P	Pathatrix 5 Pooling Dual E. coli/Salmonella Kit
APE250SDP	Pathatrix 5 Pooling E. coli O157:H7 Kit
POOL510MLN	Pathatrix 5 Pooling Kit
APL250P	Pathatrix 5 Pooling Listeria spp. Kit
APS250P	Pathatrix 5 Pooling Salmonella spp. Kit
ACART5	Pathatrix Auto Cartridge
ACART	Pathatrix Auto Cartridge
ACARTHOLD	Pathatrix Auto Cartridge Holder
ACARTRACK	Pathatrix Auto Cartridge Rack
APCAT50	Pathatrix Cationic/General Viral Capture Kit
APSAK500P	Pathatrix Cronobacter spp. 10-Pooling Kit
APSAK250P	Pathatrix Cronobacter spp. 5-Pooling Kit
APSAK50	Pathatrix Cronobacter spp. Kit, same day
APD50	Pathatrix Dual (Listeria/Salmonella spp.) Kit
APE50SD	Pathatrix E. coli O157:H7 Kit, same day
ATUBERACK	Pathatrix Elution Vessel Holder
PFF	Pathatrix Foam Filters
APL50	Pathatrix Listeria spp. Kit
APS50	Pathatrix Salmonella spp. Kit
APS50SD	Pathatrix Salmonella spp. Kit - Same Day
ATUBEHOLD	Pathatrix Sample Vessel Holder
MAGNETICPLATE	Pathatrix Magnetic Capture Plate

Dynabeads and BeadRetriever

Applied Biosystems™ & Invitrogen™ Dynabeads™ & BeadRetriever™ Products

Item Number	Description
71003-	Dynabeads anti-E. coli O157, 1mL
71004-	Dynabeads anti-E. coli O157, 5mL
71006-	Dynabeads anti-Listeria
71002-	Dynabeads anti-Salmonella
71011-	Dynabeads EPEC/VTEC O103
71009-	Dynabeads EPEC/VTEC O111
71007-	Dynabeads EPEC/VTEC O145
71013-	Dynabeads EPEC/VTEC O26
74003-	Dynabeads L10 Tubes
A10714	Dynabeads MAX E. coli O157 Kit, 1mL
A10715	Dynabeads MAX E. coli O157 Kit, 5mL
A14632	Dynabeads MAX EPEC/VTEC O121 Kit, 100 reactions
A14684	Dynabeads MAX EPEC/VTEC O121 Kit, 200 reactions
A14631	Dynabeads MAX EPEC/VTEC O45 Kit, 100 reactions
A14683	Dynabeads MAX EPEC/VTEC O45 Kit, 200 reactions
12001D	Dynabeads MPC™-1 (Magnetic Particle Concentrator)
A13346	Dynabeads MPC-S (Magnetic Particle Concentrator)
15903-	Dynabeads MX 12-tube Mixing Wheel
15904-	Dynabeads MX 8-tube Mixing Wheel

15902-	Dynabeads MX Mixer Base
15906-	Dynabeads MX Mixer Turret
15914D	Dynabeads MX MPC-6 Semi-automated Mixer Head
15911-	Dynabeads MX MPC-S Mixing Wheel
15907-	Dynabeads MX1 Mixer
15908-	Dynabeads MX2 Mixer
15910-	Dynabeads MX4 Mixer
15915D	Dynabeads MX7 Mixer
94701-	Dynabeads Rotary Mixer
94703-	Dynabeads Rotary Mixer Mixing Platform
94702-	Dynabeads Rotary Mixer Rod
74004-	Dynabeads Spot-On™ Slides
15952-	BeadRetriever Tube Rack
15950-	BeadRetriever System
15951-	BeadRetriever Tubes & Tips
12002D	MPC™-6 (Magnetic Particle Concentrator)

Applied Biosystems™ KingFisher™ Flex 96 Deep Well & Applied Biosystems™ PrepSEQ™ NA System

INSTRUMENT	
A32681	KingFisher Flex 96 Deep Well
Contact your local MBD representative	
CONSUMABLES/ Sample Prep KITS	
A24401	Lysis Buffer 1 + RNase for Food ID
4480466	PrepSEQ Nucleic Acid Extraction Kit for Food and Environmental Testing - 100 reactions
4428176	PrepSEQ Nucleic Acid Extraction Kit for Food and Environmental Testing - 300 reactions
PLASTICS	
4388476	MagMAX Express-96 Deep Well Plates
4388487	MagMAX Express-96 Deep Well Tip Combs
4388475	MagMAX Express-96 Standard Plates
ADDITIONAL ACCESSORY	
AM10050	Magnetic-Ring Stand (96 well)
4485076	96-well Pre-Clarification Tray

PrepSEQ Rapid Spin

Applied Biosystems™ PrepSEQ™ Rapid Spin Sample Preparation Kit.

- Fast, efficient sample preparation for PCR-based assays
- Optimal performance on a wide range of sample types
- Simple protocol reduces time-to-results
- Microcentrifuge-based—no additional equipment required
- Ideal for low to medium sample throughput

Designed to be used with the Applied Biosystems™ MicroSEQ™ and Applied Biosystems™ TaqMan™ assays. Can be used with other kits.

Item Number	Description
4468304	PrepSEQ Rapid Spin Sample Prep Kit - Extra Clean & Bead Beating
4407760	PrepSEQ Rapid Spin Sample Preparation Kit
4413269	PrepSEQ Rapid Spin Sample Preparation Kit - Extra Clean
4426715	PrepSEQ Rapid Spin Sample Preparation Kit - Extra Clean with Proteinase K
4426714	PrepSEQ Rapid Spin Sample Preparation Kit with Proteinase K
4480715	Proteinase K, FS

Food Molecular detection and quantification

A recently expanded range of market-leading molecular instrumentation and PCR technology for detection of food pathogens and other targets including pathogen serotypes, meat species and GMOs.

Applied Biosystems 7500 Fast Real-Time PCR Instrument

INSTRUMENT	
Laptop: A30299 Desktop: A30304	Applied Biosystems™ 7500 Fast Real-Time PCR System with Applied Biosystems™ RapidFinder™ Express Software Version 2.0 Contact your local MBD representative.
SOFTWARE	
A28811	RapidFinder Express Software Version 2.0 - Full Installation
ACCESSORIES	
A24820	Precision Plate Holder for 0.1 mL Tube Strips (for Applied Biosystems 7500 Fast System)
4330015	Applied Biosystems™ MicroAmp™ Cap Installing Tool
4358293	MicroAmp Fast 8-Tube Strip, 0.1 mL
PT0690	Precision Plate holder for MicroSEQ and SureTect Assays on Applied Biosystems 7500 Fast Instrument (to be used ONLY with PT0695 or PT0696)
PT0696	PCR carry plate for MicroSEQ Assays on Applied Biosystems 7500 Fast Instrument (to be used ONLY with PT0690)
4323032	MicroAmp Optical 8-Cap Strips

QuantStudio™ 5 Real-Time PCR Instrument

INSTRUMENT	
A36320	Applied Biosystems™ QuantStudio™ 5 Instrument with a 0.1mL block, packaged with power cord desktop computer, monitor, RapidFinder Analyzer v1.0 software and blank USB flashdrive.
A36328	QuantStudio 5 Instrument with a 0.1mL block, packaged with power cord laptop computer, RapidFinder Analyzer v1.0 software and blank USB flashdrive.
A36329	QuantStudio 5 Instrument with a 0.1mL block, packaged with power cord, No Computer, RapidFinder Analyzer v1.0 software and blank USB flashdrive.
A24811	SimpliAmp Thermal Cycler For use with the Applied Biosystems 7500 Fast and QuantStudio 5 PCR Instruments and the RapidFinder Salmonella species, Typhimurium and Enteritidis Multiplex PCR Assay Contact your local MBD representative.
SOFTWARE	
A34777	RapidFinder Analysis Software V1.0 for use with the QuantStudio 5 Instrument 0.1 mL block Full installation Contact your local MBD representative.

Applied Biosystems™ MicroSEQ™/Applied Biosystems™ RapidFinder™ Assays

The MicroSEQ product range of assays are based on TaqMan® real-time PCR. Real-time PCR is a proven method for pathogen detection and has been applied successfully to a wide range of foodborne pathogens. The method tests for food pathogens by amplifying and detecting a DNA target sequence that is specific to the organism under investigation. In this way, real-time PCR can detect pathogens quickly, simply, and with outstanding sensitivity and specificity. Since the MicroSEQ assays uses fast PCR chemistry, the PCR run time takes only 40 minutes. The MicroSEQ range and RapidFinder STEC solutions have been validated by AOAC RI PTM, AFNOR, and DAFF.

Item Number	Description
4427409	MicroSEQ E. coli O157:H7 Detection Kit
4445656	MicroSEQ E. coli O157:H7 Detection Starter Kit with PrepSEQ Nucleic Acid Extraction Kit
4445657	MicroSEQ E. coli O157:H7 Detection Starter Kit with PrepSEQ Rapid Spin Sample Preparation Kit
4403874	MicroSEQ Listeria monocytogenes Detection Kit
4415045	MicroSEQ Listeria monocytogenes Detection Starter Kit with PrepSEQ Rapid Spin Sample Preparation Kit
4412637	MicroSEQ Listeria monocytogenes Detection Starter Kit with PrepSEQ Nucleic Acid Extraction Kit
4445659	MicroSEQ Listeria spp. Detection Starter Kit with PrepSEQ Rapid Spin Sample Preparation Kit
4427410	MicroSEQ Listeria spp. Detection Kit
4445658	MicroSEQ Listeria spp. Detection Starter Kit with PrepSEQ Nucleic Acid Extraction Kit
4403930	MicroSEQ Salmonella spp. Detection Kit
4412639	MicroSEQ Salmonella spp. Detection Starter Kit with PrepSEQ Nucleic Acid Extraction Kit
4415044	MicroSEQ Salmonella spp. Detection Starter Kit with PrepSEQ Rapid Spin Sample Preparation Kit - Extra Clean
4480467	RapidFinder STEC Screening Assay Starter Kit - with PrepSEQ Nucleic Acid Extraction Kit
4476886	RapidFinder STEC Screening Assay
4476901	RapidFinder STEC Confirmation Assay
A33227	RapidFinder Salmonella species, Typhimurium and Enteritidis Multiplex PCR Assay

Applied Biosystems™ TaqMan® Pathogen Detection Kits

TaqMan® Pathogen Detection Kits are part of a systems approach that helps user to detect food pathogens both rapidly and accurately. The kits use molecular methods (i.e., Real-Time PCR) to simplify handling and to reduce the time-to-result. The kits have undergone rigorous internal testing. All TaqMan® Pathogen Detection Kits use a uniform protocol to facilitate simplified workflows and generate consistent results. Both the PCR reaction and the PCR run protocol are standardized so that multiple pathogens can be screened on the same plate in the same PCR run.



Item Number	Description
4368736	TaqMan® Listeria monocytogenes Detection Kit with PrepMan Ultra Sample Preparation Reagent
4368723	TaqMan® Salmonella enterica Detection Kit & PrepMan Ultra Reagent (AFNOR and AOAC Validated)
4366102	TaqMan® Listeria monocytogenes Detection Kit
4366104	TaqMan® Salmonella enterica Detection Kit
4382492	TaqMan® Cronobacter sakazakii Detection Kit
4368606	TaqMan® Staphylococcus aureus Detection Kit

Meat speciation Imegen™ Real-Time PCR Tests

Reliable detection and quantification of species DNA in collaboration with Imegen, we are offering a broad portfolio of real-time PCR based species detection and quantification testing solutions, with very high sensitivity below the officially required limit of detection (LOD) of 0,1%.

For high-throughput isolation:

Item Number	Description
A24401	Lysis Buffer 1
4428176	RNase for Food ID
4480466	PrepSEQ Nucleic Acid Extraction Kit Cat. nos.
4466336	GMO Extraction Kit
A24391	RapidFinder Beef ID Kit
A24393	RapidFinder Chicken ID Kit
A15570	RapidFinder Equine ID Kit
A24392	RapidFinder Pork ID Kit
A24397	RapidFinder Poultry ID Kit
A15579	RapidFinder Quant Equine Set
A24399	RapidFinder Quant Multi-Meat Set
A24396	RapidFinder Ruminant ID Kit
A24395	RapidFinder Sheep ID Kit
A24394	RapidFinder Turkey ID Kit
A24398	RapidFinder Fish ID Kit

GMO testing portfolio

In partnership with Imegen, a GMO testing portfolio containing solutions for the complete workflow including nucleic acid extraction, GMO screening and quantification kit, instrument, technical support and GMO expertise.

Thermo Scientific GMO testing solutions provide a high level of confidence:

- applicable on processed food due to MGB proves
- high sensitivity due to high sample input
- optimal GMO coverage due to unique design set up

Item Number	Description
4466366	GMO Extraction Kit
4481972	TaqMan® GMO Maize Quantification Kit
4466334	TaqMan® GMO Screening Kit
4466335	TaqMan® Roundup Ready Soya Quantification Kit

For a complete list of Meat Species ID and GMO Real-Time PCR Assays, please contact your local Thermo Fisher Microbiology Food Safety Specialist.

Applied Biosystems™ TaqMan® Environmental Master Mixes

Optimized for analyzing environmental, food, and other challenging samples, TaqMan® Fast Environmental Master Mixes offer accurate microbial detection in the presence of high levels of problematic inhibitors.

Item Number	Description	Number of Reactions
4396838	TaqMan® Environmental Master Mix 2.0	200 rxn
4398021	TaqMan® Environmental Master Mix 2.0	400 rxn
4388044	TaqMan® Environmental Master Mix 2.0	800 rxn
4488681	TaqMan® Fast Environmental Master Mix Beads + IPC	96 rxn
4488676	TaqMan® Fast Environmental Master Mix Beads w/o IPC	96 rxn

Applied Biosystems™ RNA UltraSense™ One-Step Quantitative RT-PCR System

Specially designed for amplification and real-time detection of RNA viruses and ultra-low abundance transcripts, this optimized, ultra-concentrated system provides greater priming specificity, higher product yields, and detection over a broad dynamic range.

Item Number	Description	Number of Reactions
11732927	RNA UltraSense One-Step Quantitative RT-PCR System	100 rxn

Custom TaqMan® Real-Time PCR Assay Beads

Our food safety custom lyophilized assay beads provide a unique offering of solutions specifically designed for the food safety testing market. These assays allow you to bring the same specificity and sensitivity of PCR beyond the big 4 pathogens.

Item Number	Description
4485027	Custom TaqMan® Campylobacter Multiplex Assay Beads
4485034	Custom TaqMan® Cronobacter sakazakii Assay Beads
4485039	Custom TaqMan® Escherichia coli spp. Assay Beads
4485045	Custom TaqMan® Enterococcus spp Assay Beads
4485046	Custom TaqMan® Oncorhynchus keta Assay Beads
4485049	Custom TaqMan® Salmonella spp. Ultimate Assay Beads
4485050	Custom TaqMan® Salmonella Heidelberg Assay Beads
4485053	Custom TaqMan® Salmonella Typhimurium Assay Beads
4485059	Custom TaqMan® STX1/STX2 Assay Beads
4485063	Custom TaqMan® STEC O103 & O145 Assay Beads, MLG
4485064	Custom TaqMan® STEC O26 & O111 Assay Beads, MLG
4485065	Custom TaqMan® STEC O45 & O121 Assay Beads, MLG
4485066	Custom TaqMan® STEC STX & EAE Assay Beads, MLG
4485068	Custom TaqMan® Vibrio Multiplex Assay Beads
4485074	Custom TaqMan® Verotoxin-producing Escherichia coli VT1/VT2 Multiplex Assay Beads, HPA
4485075	Custom TaqMan® ISO STEC STX 1/2, EAE, O157 Screening Assay Beads
4485077	Custom TaqMan® Salmonella Triplex (SP/SE/ST) Assay Beads
4485082	Custom TaqMan® ISO STEC O45 & O121 Assay Beads
4485083	Custom TaqMan® ISO STEC O26/O145/O103 Assay Beads
4485084	Custom TaqMan® ISO STEC O111 & O104 Assay Beads
4485086	Custom TaqMan® Salmonella spp. Ultimate / Enteritidis Multiplex Assay Beads
A29797	Custom TaqMan® Enterobacteriaceae 8-Genus Assay Beads

For a complete list of Custom TaqMan® Real-Time PCR Assays, please contact your local Thermo Fisher Food Safety Specialist

Thermo Fisher Scientific Food Safety Custom Services

Thermo Fisher Scientific combines innovative assay design, sequencing and bioinformatics expertise with comprehensive global support to deliver personalized food safety technologies that help protect your brand and bottom line. If you do not see a solution that specifically meets your needs, we will work with you to develop a solution to streamline your microbial testing workflow.

Our Food Safety custom services include, but are not limited to:

- Novel Assay Development
- Multiplexed Assay Solutions
- Workflow Optimization
- Small-Sale Sequencing
- Validation & Verification Services



SureTect Real-Time PCR System

Speed and Performance

Now, you can go into every food safety test confidently with the Thermo Scientific™ SureTect™ Real-Time PCR System. Built on proven PCR technology and backed by world-class service and support, the SureTect System is designed to quickly and accurately detect microorganisms in a broad range of foods and associated samples. This unique solution combines speed and performance in an easy-to-use, cost-effective platform – giving you results you can be sure about.

Instruments

Applied Biosystems™ 7500 Fast Real-Time PCR instrument	
Applied Biosystems 7500 Fast Real-Time PCR System with Applied Biosystems™ RapidFinder™ Express Software Version 2.0	Laptop: A30299 Desktop: A30304

Software

RapidFinder Express Software Version 2.0 - For Applied Biosystems 7500 Fast Real-Time PCR instrument	
NOTE: RapidFinder Express Software 2.0 for Applied Biosystems 7500 and SureTect Assay Kit File for Applied Biosystems 7500 Fast are BOTH required to run the SureTect Assays on the 7500 Fast Real Time PCR instrument	A28811

Accessories

Starter Pack for SureTect Assays for SureTect PikoReal 24-well PCR Cycler (includes all required pipettes, equipment and accessories)	
each	PT0600
Starter Pack for High Throughput SureTect Assays for Applied Biosystems 7500 96 well Fast Real-Time PCR instrument (includes all required pipettes, equipment and accessories)	
each	PT0675
Precision Plate holder for MicroSEQ and SureTect Assays on Applied Biosystems 7500 Fast (to be used ONLY with PT0695 or PT0696)	
each	PT0690
PCR carry plate for SureTect Assays on Applied Biosystems 7500 Fast (to be used ONLY with PT0690)	
each	PT0695
SureTect Lysis carry plate	
each	PT0680
SureTect lysis carry plate handle	
each	PT0681

Starter packs include all consumables required

Warranty

PikoReal Standard Service Contract (year 2)	
each	PT0850
PikoReal Premium Service Contract (year 2)	
each	PT0851
AB Assurance Service Plan for Applied Biosystems 7500 Fast Real-Time PCR Instrument Instrument Service	
each	ZG11SC7500FASTAMT
AB Complete Service Plan for Applied Biosystems 7500 Fast Real-Time PCR Instrument Instrument Service	
each	ZG21SC7500FASTAMT
OQ/IPV Assurance Service for Applied Biosystems 7500 Fast Real-Time PCR instrument	
each	4365571

Assays

SureTect Salmonella Species Assay	
each	PT0100A
SureTect Listeria Species Assay	
each	PT0200A
SureTect Listeria monocytogenes Assay	
each	PT0300A
SureTect E. coli 0157:H7 Assay	
each	PT0400A
SureTect Cronobacter sakazakii Assay	
each	PT1060A

Hygiena Bax System

The Hygiena™ BAX™ System allows food laboratories to quickly and accurately detect microorganisms. Computer-generated results are clearly displayed on-screen.

The Hygiena BAX System Q7 is available in Europe, Canada, Australia and New Zealand from your local Thermo Fisher Scientific representative.

Automated BAX™ System Q7 start up package

7500 Fast PCR cycler/detector, computer/monitor/keyboard/mouse/printer, BAX software/Microsoft® Windows® based, dry block heaters (2), cooling blocks (3), capping/decapping tools, pipettes, disposables to run the system, 1 year warranty.

each		HYBASY2020
------	--	------------

BAX™ System Real-Time Campylobacter coli/jejuni/lari Assay

96 tests		HYBKIT2018
----------	--	------------

BAX™ System Real-Time Escherichia coli O157:H7 Assay

This assay can be run in parallel with the BAX System Real-Time STEC Screening Assay where STEC other than *E. coli* O157:H7 may be present.

96 tests		HYBKIT2000
----------	--	------------

BAX System Real time Salmonella Assay

96 tests		HYBKIT2006
----------	--	------------

BAX System Real Time Software V 3.6

		HYBMIS210
--	--	-----------

BAX™ System Real-Time Genus Listeria Assay

96 tests		HYBKIT2019
----------	--	------------

BAX™ System Real-Time Listeria monocytogenes Assay

96 tests		HYBKIT2005
----------	--	------------

BAX™ System Real-Time Staphylococcus aureus Assay

96 tests		HYBKIT2020
----------	--	------------

BAX™ System Real-Time STEC Screening Assay

For the detection of *stx* and *eae* genes associated with Shiga-toxin producing strains of *Escherichia coli* (STEC). Following screening, presence of genes specific for O antigens of the most common STEC can be determined using the BAX System Real-Time STEC Panel 1 and Panel 2 Assays and the BAX System Real-Time Escherichia coli O157:H7 MP Assay.

96 tests		HYBKIT2021
----------	--	------------

BAX™ System Real-Time STEC Panel 1 Assay (O26, O111, O121)

An assay for determining the presence of specific genes specific to STEC O antigens: O26, O111 and O121. For use following screening using the BAX System Real-Time STEC Screening Assay. This assay may be run in parallel with the BAX System Real-Time STEC Panel 2 Assay and the BAX System Real-Time Escherichia coli O157:H7 MP Assay.

48 tests		HYBKIT2008
----------	--	------------

BAX™ System Real-Time STEC Panel 2 Assay (O45, O103, O145)

An assay for determining the presence of specific genes specific to STEC O antigens: O45, O103 and O145. For use following screening using the BAX System Real-Time STEC Screening Assay. This assay may be run in parallel with the BAX System Real-Time STEC Panel 1 Assay and the BAX System Real-Time Escherichia coli O157:H7 MP Assay.

48 tests		HYBKIT2009
----------	--	------------

BAX™ System Real-Time Vibrio cholerae/parahaemolyticus/vulnificus Assay

96 tests		HYBKIT2010
----------	--	------------

BAX™ System Cronobacter sakazakii Assay

96 tests		HYBKIT2001
----------	--	------------

BAX™ System Escherichia coli O157:H7 MP Assay

96 tests		HYBKIT2004
----------	--	------------

BAX™ System Genus Listeria 24E Assay

96 tests		HYBKIT2003
----------	--	------------

BAX™ System Listeria monocytogenes 24E Assay

96 tests		HYBKIT2002
----------	--	------------

BAX™ System Salmonella Assay

96 tests		HYBKIT2012
----------	--	------------

BAX™ System Salmonella 2 Assay

96 tests		HYBKIT2011
----------	--	------------

BAX™ System Listeria monocytogenes Assay

96 tests		HYBKIT2017
----------	--	------------

BAX™ System Genus Listeria Assay

96 tests		HYBKIT2016
----------	--	------------

BAX™ System Accessories

Cluster Tube Caps		
120 Strips		HYBCAP2008
Cluster Tube Capping/Decapping Tool		
each		HYBMIS2011
Cluster Tubes		
10x96		HYBMIS2012
PCR Tube Holder Tray for Low Profile Tubes		
each		HYBMIS2040
Pipettor – Finn 20–200µL adjustable		
each		HYBMIS2043
Pipettor – Finn 5–50µL adjustable		
each		HYBMIS2041
Powder Free Nitrile Gloves		
100		HYBMIS2028
Tip 250µL with no Barrier		
10x96		HYBMIS2058
Tip 100µL Extended with Barrier		
10x96		HYBMIS2057
Tip 5mL for Stepper Pipette		
25		HYBMIS2055
Buffered Peptone Water Dry-Bags		
10 bags	for 10x20 litres medium	DB0509M
Buffered Peptone Water Dry-Bags without filters		
10 bags	for 10x20 litres medium	DB0509W

Microbial Characterization System

Start-Up Pack

Includes Hygiena RiboPrinter™ characterization unit; flat-screen LCD monitor; computer with CD-RW and USB ports for connecting peripherals; online help manual; RiboPrinter™ system software with data security tools for compliance with 21 CFR Part 11; USB mass storage device for data transfer; heat treatment station; sample preparation rack and hand-held mixer

instrument		HYRASY2005
------------	--	------------

Kits

EcoR I RiboPrinter Kit

48 tests	(QR0517C, QR0536C, QR3120C, QR3116C & QR3118C)	HYRASY2077
48 tests	(QR0517C, QR0536C, QR3120C & QR3116C)	HYRMCH2050*

PVU II RiboPrinter Kit

48 tests	(QR0600C, QR0536C, QR3120C, QR3116C & QR3118C)	HYRASY2078
48 tests	(QR0600C, QR0536C, QR3120C & QR3116C)	HYRASY2072*

No-Enzyme RiboPrinter Kit

48 tests	(QR0600C, QR0536C, QR3120C, QR3116C & QR3118C)	HYRASY2071
48 tests	(QR0526C, QR0536C, QR3120C & QR3116C)	HYRASY2076*

Components

DNA Preparation Reagents with ECOR I Restriction Enzyme		
48 tests	6x ECOR I enzyme vials, lysing agents A and B, lactic agent	HYRKIT2029
DNA Preparation Reagents with PVU II Restriction Enzyme		
48 tests	6x PVU II enzyme vials, lysing agents A and B	HYRKIT2032
DNA Preparation Reagents without Restriction Enzyme		
48 tests	6x lysing agents A and B, lactic agent	HYRKIT2031
DNA Prep and MP insert		
48 tests	6x probe, substrate, conjugate and DNA prep packs	HYRKIT2030
Gel Cassette & Membrane		
6+6		HYRKIT2033
MP Base		
6	Unit for holding reagents	HYRKIT2034
Purified Water		
48 tests	2 x containers, each sufficient for 3-4 batches	HYRKIT2036
Sample Set Prep Boxes		
480 tests	60 sample carriers, 1 bottle of sample buffer and 500 colony picks	HYRKIT2035

*RiboPrinter kits without water

PathoProof System

Thermo Scientific™ PathoProof™ PCR Assays

Maximising milk productivity, safeguarding milk quality and sustaining herd health are important goals of every dairy operation.

Thermo Scientific™ PathoProof™ PCR Assays identify mastitis-causing bacteria from bovine milk typically in just four hours. The assays are based on an accurate and sensitive real-time PCR. The tests can be used for analyzing fresh and preserved quarter milk samples, milk recording samples, and bulk tank milk samples.

PathoProof Complete - 16 kit

For the identification of *S. aureus*, *Staphylococcus spp.* (including all major coagulase-negative staphylococci), *S. agalactiae*, *S. dysgalactiae*, *S. uberis*, *E. coli*, *S. marcescens*, *Enterococcus spp.* (including *E. faecalis* and *E. faecium*), *K. oxytoca* and/or *K. pneumoniae*, *C. bovis*, *T. pyogenes* and *Peptoniphilus (Peptostreptococcus) indolicus*, *Staphylococcal β-lactamase* gene (penicillin resistance gene), *M. bovis*, *Mycoplasma spp.*, *Yeast*, *Prototheca spp.*

PathoProof Complete – 16 Reagent Kit for Agilent™ Mx3005P™ qPCR instrument

50 tests	PF1600S
4 x 96 tests	PF1600L
4 x 96 tests on Kingfisher Flex Magnetic Particle Processor	PF1600LKF

PathoProof Complete – 16 qPCR reagents for Agilent Mx3005P qPCR instrument

for use with PF1600S	PF1600SB
for use with PF1600L and PF1600LKF	PF1600LB

PathoProof Complete – 16 Reagent Kits for Applied Biosystems™ 7500 Fast qPCR instrument

50 tests	PF1650S
4 x 96 tests	PF1650L
4 x 96 tests on Kingfisher Flex Magnetic Particle Processor	PF1650LKF

PathoProof Complete – 16 qPCR reagents for Applied Biosystems 7500 Fast qPCR instrument

for use with PF1650S	PF1650SB
for use with PF1650L and PF1650LKF	PF1650LB

PathoProof Mycoplasma-8 kits

For the identification of *M. alkalescens*, *M. bovigenitalium*, *M. bovis*, *M. californicum*, *M. canadense*, *Mycoplasma spp.*, *S. aureus*, *S. agalactiae*.

PathoProof Mycoplasma – 8 Reagent for Agilent Mx3005P qPCR instrument Kits

50 tests	PF8500S
4 x 96 tests	PF8500L
4 x 96 tests on Kingfisher Flex Magnetic Particle Processor	PF8500LKF

PathoProof Mycoplasma – 8 qPCR reagents for Agilent Mx3005P qPCR instrument

for use with PF8500S	PF8500SB
for use with PF8500L and PF8500LKF	PF8500LB

PathoProof Mycoplasma – 8 Reagent Kits for Applied Biosystems 7500 Fast qPCR instrument

50 tests	PF8100S
4 x 96 tests	PF8100L
4 x 96 tests on Kingfisher Flex Magnetic Particle Processor	PF8100LKF

PathoProof Mycoplasma – 8 qPCR reagents for Applied Biosystems 7500 Fast qPCR instrument

for use with PF8100S	PF8100SB
for use with PF8100L and PF8100LKF	PF8100LB

PathoProof qPCR reagents must be ordered with a matched PathoProof Reagent Kit for DNA extraction.

PathoProof Major 4.2 Kit

For the identification of *M. bovis*, *S. aureus*, *S. agalactiae*, *S. uberis*.

PathoProof Major – 4.2 Reagent Kit for Agilent Mx3005P qPCR instrument

50 tests	PF0942S
4 x 96 tests	PF0942L
4 x 96 tests on Kingfisher Flex Magnetic Particle Processor	PF0942LKF

PathoProof Major - 4.2 qPCR reagents for Agilent Mx3005P qPCR instrument

for use with PF0942S	PF0942SB
for use with PF0942L	PF0942LB

PathoProof Major - 4.2 Reagent Kits for Applied Biosystems 7500 Fast qPCR instrument

50 tests	PF4250S
4 x 96 tests	PF4250L
4 x 96 tests on Kingfisher Flex Magnetic Particle Processor	PF4250LKF

PathoProof Major – 4.2 qPCR reagents for Applied Biosystems 7500 Fast qPCR instrument

for use with PF4250S	PF4250SB
for use with PF4250L and PF4250LKF	PF4250LB

PathoProof Major 3 Kit

PathoProof Major – 3 Reagent Kit for Agilent Mx3005P qPCR instrument

50 tests	PF0914S
4 x 96 tests	PF0914L
4 x 96 tests on Kingfisher Flex Magnetic Particle Processor	PF0914LKF

PathoProof Major – 3 qPCR reagents for Agilent Mx3005P qPCR instrument

for use with PF0914S	PF0914SB
for use with PF0914L	PF0914LB

PathoProof Major - 3 Reagent Kits for Applied Biosystems 7500 Fast qPCR instrument

50 tests	PF3050S
4 x 96 tests	PF3050L
4 x 96 tests on Kingfisher Flex Magnetic Particle Processor	PF3050LKF

PathoProof Major -3 qPCR Reagents for Applied Biosystems 7500 Fast qPCR instrument

for use with PF3050S	PF3050LSB
for use with PF3050L and PF3050LKF	PF3050LB

PathoProof software

PathoProof Norden™ Studio Software Licence

For up to 5 Norden Studio Users	PF0888A
---------------------------------	---------

Thermo Scientific Culti-Loops Organisms

Thermo Scientific™ Culti-Loops™ Qualitative Quality Control Organisms are available in packs of five. Culti-Loops organisms should be stored at 2–8°C, with the exception of *Campylobacter* spp. which should be kept frozen (-18°C).

Organism	Order Code (5 loops)
----------	----------------------

A	
Achromobacter xylosoxidans subsp. xylosoxidans ATCC® 27061™*	R4601012
Acinetobacter baumannii ATCC® 19606™*	R4601007
Acinetobacter baumannii ATCC® BAA-747™*	R4601009
Acinetobacter lwoffii ATCC® 17925™*	R4600503
Acinetobacter sp. ATCC® 49137™*	R4601002
Acinetobacter sp. ATCC® 49139™*	R4600501
Acinetobacter sp. ATCC® 49466™*	R4601008
Acinetobacter sp. ATCC® 9957™*	R4609204
Actinomyces odontolyticus ATCC® 17929™*	R4600500
Actinomyces viscosus ATCC® 15987™*	R4600510
Actinomyces viscosus ATCC® 43146™*	R4609459
Aerococcus viridans ATCC® 11563™*	R4601019
Aerococcus viridans ATCC® 700406™*	R4601022
Aeromonas caviae ATCC® 15468™*	R4601021
Aeromonas hydrophila ATCC® 35654™*	R4601013
Aeromonas hydrophila ATCC® 49140™*	R4601003
Aeromonas hydrophila ATCC® 7965™*	R4609350
Aeromonas hydrophila ATCC® 7966™*	R4601020
Aeromonas salmonicida ATCC® 33658™*	R4609456
Aeromonas veronii biogroup sobria ATCC® 9071™*	R4601005
Aggregatibacter aphrophilus ATCC® 33389™*	R4603020
Alcaligenes faecalis ATCC® 35655™*	R4601011
Alcaligenes faecalis subsp. faecalis ATCC® 8750™*	R4601006
Alternaria alternata TX 8025	R4601016
Aneurinibacillus aneurinolyticus ATCC® 11376™*	R4601015
Arcanobacterium haemolyticum ATCC® BAA-1784™*	R4609413
Arcanobacterium pyogenes ATCC® 19411™*	R4600504
Arcanobacterium pyogenes ATCC® 49698™*	R4600502
Aspergillus brasiliensis ATCC® 16404™*	R4601100
Aspergillus brasiliensis ATCC® 9642™*	R4601102
Aspergillus fumigatus KM 8001	R4601018
Aspergillus oryzae ATCC® 10124™*	R4601017
Atopobium vaginae ATCC® BAA-55™*	R4609450

B	
Bacillus badius ATCC® 14574™*	R4609394
Bacillus cereus ATCC® 11778™*	R4601220
Bacillus cereus ATCC® 14579™*	R4601217
Bacillus circulans ATCC® 61™*	R4601216
Bacillus licheniformis ATCC® 12759™*	R4601233
Bacillus megaterium ATCC® 14581™*	R4609395
Bacillus pumilus ATCC® BAA-1434™*	R4609379
Bacillus stearothermophilus ATCC® 10149™*	R4601202
Bacillus subtilis ATCC® 6633™*	R4601221
Bacteroides fragilis ATCC® 23745™*	R4601251
Bacteroides fragilis ATCC® 25285™*	R4601250
Bacteroides ovatus ATCC® 8483™*	R4601222
Bacteroides ovatus ATCC® BAA-1296™*	R4609374
Bacteroides ovatus ATCC® BAA-1304™*	R4609376
Bacteroides thetaiotaomicron ATCC® 29741™*	R4601260
Bacteroides uniformis ATCC® 8492™*	R4601259
Bacteroides ureolyticus ATCC® 33387™*	R4601201
Bacteroides vulgatus ATCC® 8482™*	R4601261
Bifidobacterium breve ATCC® 15700™*	R4606801
Bordetella bronchiseptica ATCC® 10580™*	R4601224
Bordetella bronchiseptica ATCC® 4617™*	R4601223
Bordetella pertussis ATCC® 12742™*	R4609213
Bordetella pertussis ATCC® 9340™*	R4601225
Brevibacillus agri ATCC® 51663™*	R4607217
Brevibacillus laterosporus ATCC® 64™*	R4607228
Brevundimonas diminuta ATCC® 11568™*	R4607219
Brevundimonas diminuta ATCC® 19146™*	R4607200
Brochothrix thermosphacta ATCC® 11509™*	R4609038
Burkholderia cepacia ATCC® 25416™*	R4605220
Burkholderia cepacia ATCC® 25608™*	R4601232

Organism	Order Code (5 loops)
C	
Campylobacter coli ATCC® 33559™*	R4609039
Campylobacter coli ATCC® 43478™*	R4609387
Campylobacter jejuni ATCC® 29428™*	R4601410
Campylobacter jejuni ATCC® BAA-1153™*	R4601440
Campylobacter jejuni subsp. jejuni ATCC® 33291™*	R4601400
Campylobacter jejuni subsp. jejuni ATCC® 33292™*	R4607070
Candida albicans ATCC® 10231™*	R4601503
Candida albicans ATCC® 14053™*	R4601501
Candida albicans ATCC® 2091™*	R4601523
Candida albicans ATCC® 36232™*	R4601239
Candida albicans ATCC® 60193™*	R4601502
Candida albicans ATCC® 66027™*	R4601495
Candida albicans ATCC® 90028™*	R4601496
Candida dubliniensis ATCC® MYA-577™*	R4609463
Candida geochares ATCC® 36852™*	R4601401
Candida glabrata ATCC® 15126™*	R4601513
Candida glabrata ATCC® 2001™*	R4601504
Candida glabrata ATCC® 66032™*	R4601511
Candida glabrata ATCC® MYA-2950™*	R4601505
Candida guilliermondii ATCC® 6260™*	R4601521
Candida kefyr ATCC® 204093™*	R4608013
Candida kefyr ATCC® 2512™*	R4601507
Candida kefyr ATCC® 66028™*	R4601512
Candida krusei ATCC® 14243™*	R4601510
Candida krusei ATCC® 34135™*	R4609464
Candida lusitanae ATCC® 34449™*	R4601243
Candida lusitanae ATCC® 42720™*	R4601506
Candida lusitanae ATCC® 66035™*	R4601516
Candida parapsilosis ATCC® 22019™*	R4601518
Candida parapsilosis ATCC® 34136™*	R4601515
Candida parapsilosis ATCC® 90018™*	R4609458
Candida tropicalis ATCC® 1369™*	R4609465
Candida tropicalis ATCC® 201380™*	R4601244
Candida tropicalis ATCC® 13803™*	R4609240
Candida tropicalis ATCC® 66029™*	R4601517
Candida tropicalis ATCC® 750™*	R4601240
Candida utilis ATCC® 9950™*	R4608012
Cellulosimicrobium cellulans ATCC® 27402™*	R4601544
Citrobacter braakii ATCC® 10625™*	R4601801
Citrobacter diversus KM 11012	R4601508
Citrobacter freundii ATCC® 8090™*	R4601800
Clostridium barati ATCC® 27638™*	R4601615
Clostridium difficile ATCC® 43255™*	R4601603
Clostridium difficile ATCC® 700057™*	R4609452
Clostridium difficile ATCC® 9689™*	R4601610
Clostridium difficile ATCC® BAA-1870™*	R4609433
Clostridium histolyticum ATCC® 19401™*	R4601601
Clostridium novyi ATCC® 7659™*	R4601605
Clostridium novyi Type A ATCC® 19402™*	R4609220
Clostridium septicum ATCC® 12464™*	R4601650
Clostridium sordellii ATCC® 9714™*	R4601602
Clostridium sporogenes ATCC® 11437™*	R4601703
Clostridium sporogenes ATCC® 19404™*	R4601700
Clostridium sporogenes ATCC® 3584™*	R4601701
Clostridium tertium ATCC® 19405™*	R4601509
Corynebacterium diphtheriae ATCC® 13812™*	R4609234
Corynebacterium jeikeium ATCC® 43734™*	R4601901
Corynebacterium minutissimum ATCC® 23348™*	R4601904
Corynebacterium pseudodiphtheriticum ATCC® 10700™*	R4601902
Corynebacterium pseudodiphtheriticum ATCC® 10701™*	R4601903
Corynebacterium renale ATCC® 19412™*	R4601906
Corynebacterium renale ATCC® BAA-1785™*	R4609415
Corynebacterium striatum ATCC® BAA-1293™*	R4609308
Corynebacterium urealyticum ATCC® 43044™*	R4609414
Corynebacterium xerosis ATCC® 373™*	R4601905
Cronobacter mytjensii ATCC® 51329™*	R4607093
Cryptococcus albidus ATCC® 66030™*	R4601710
Cryptococcus albidus var. albidus ATCC® 10666™*	R4601811
Cryptococcus albidus var. albidus ATCC® 34140™*	R4601813
Cryptococcus humicolus ATCC® 9949™*	R4601101
Cryptococcus laurentii ATCC® 18803™*	R4601241
Cryptococcus laurentii ATCC® 66036™*	R4601711
Cryptococcus laurentii ATCC® 76483™*	R4601910
Cryptococcus neoformans ATCC® 14116™*	R4601812
Cryptococcus neoformans ATCC® 204092™*	R4608016
Cryptococcus neoformans ATCC® 32045™*	R4609327
Cryptococcus neoformans ATCC® 34877™*	R4601714
Cryptococcus neoformans ATCC® 56991™*	R4609466
Cryptococcus neoformans ATCC® 66031™*	R4601712
Cryptococcus neoformans ATCC® 76484™*	R4601912
Cryptococcus uniguttulatus ATCC® 66033™*	R4601713
Curtobacterium pusillum ATCC® 19096™*	R4609418
Curvularia sp. KM 8023	R4601705

Organism	Order Code (5 loops)
E	
Edwardsiella tarda ATCC® 15947™*	R4601973
Eggerthella lenta ATCC® 43055™*	R4601951
Eikenella corrodens ATCC® 23834™*	R4601974
Eikenella corrodens ATCC® BAA-1152™*	R4601947
Elizabethkingia meningoseptica ATCC® 13253™*	R4602002
Enterobacter aerogenes ATCC® 13048™*	R4607080
Enterobacter aerogenes ATCC® 35028™*	R4607081
Enterobacter aerogenes ATCC® 35029™*	R4607082
Enterobacter aerogenes ATCC® 49701™*	R4601941
Enterobacter cloacae ATCC® 13047™*	R4607090
Enterobacter cloacae ATCC® 23355™*	R4607091
Enterobacter cloacae ATCC® 35030™*	R4601945
Enterobacter gergoviae ATCC® 33028™*	R4601940
Enterobacter hormaechei ATCC® 700323™*	R4607094
Enterococcus avium ATCC® 14025™*	R4607028
Enterococcus casseliflavus ATCC® 700327™*	R4607099
Enterococcus durans ATCC® 11576™*	R4601957
Enterococcus durans ATCC® 49479™*	R4601952
Enterococcus durans ATCC® 6056™*	R4607032
Enterococcus faecalis ATCC® 19433™*	R4601990
Enterococcus faecalis ATCC® 29212™*	R4607030
Enterococcus faecalis ATCC® 49149™*	R4607087
Enterococcus faecalis ATCC® 49452™*	R4601984
Enterococcus faecalis ATCC® 51299™*	R4601996
Enterococcus faecalis ATCC® 7080™*	R4606854
Enterococcus faecium ATCC® 35667™*	R4711956
Enterococcus faecium ATCC® 51559™*	R4609383
Enterococcus faecium ATCC® 700221™*	R4609432
Enterococcus gallinarum ATCC® 700425™*	R4601958
Enterococcus hirae ATCC® 49135™*	R4607086
Enterococcus raffinosus ATCC® 49464™*	R4601994
Enterococcus saccharolyticus ATCC® 43076™*	R4609393
Epidermophyton floccosum ATCC® 52066™*	R4601977
Erysipelothrix rhusiopathiae ATCC® 19414™*	R4601978
Escherichia coli ATCC® 10536™*	R4601968
Escherichia coli ATCC® 11229™*	R4601997
Escherichia coli ATCC® 11775™*	R4601989
Escherichia coli ATCC® 12014™*	R4609352
Escherichia coli ATCC® 13706™*	R4609469
Escherichia coli ATCC® 25922™*	R4607050
Escherichia coli ATCC® 29194™*	R4601986
Escherichia coli ATCC® 35218™*	R4601971
Escherichia coli ATCC® 35421™*	R4607084
Escherichia coli ATCC® 4157™*	R4607088
Escherichia coli ATCC® 51446™*	R4607097
Escherichia coli ATCC® 51755™*	R4601998
Escherichia coli ATCC® 8739™*	R4607085
Exiguobacterium aurantiacum ATCC® 49676™*	R4601900
Exophiala jeanselmei ATCC® 10224™*	R4601980
F	
Fingoldia magna ATCC® 29328™*	R4605002
Fluoribacter bozemanii ATCC® 33217™*	R4603951
Fluoribacter dumoffii ATCC® 33279™*	R4603952
Fonsecaea pedrosii ATCC® 28174™*	R4601981
Fusobacterium mortiferum ATCC® 25557™*	R4609207
Fusobacterium mortiferum ATCC® 9817™*	R4602009
Fusobacterium necrophorum ATCC® 25286™*	R4609375
Fusobacterium nucleatum ATCC® 10953™*	R4602012
Fusobacterium nucleatum subsp. nucleatum ATCC® 25586™*	R4602010
Fusobacterium oxysporum ATCC® 48112™*	R4601983
G	
Gardnerella vaginalis ATCC® 14018™*	R4602050
Gardnerella vaginalis ATCC® 49145™*	R4602051
Gemella morbillorum ATCC® 27824™*	R4602100
Geobacillus stearothermophilus ATCC® 10149™*	R4601202
Geobacillus stearothermophilus ATCC® 12978™*	R4603011
Geobacillus stearothermophilus ATCC® 12980™*	R4601219
Geobacillus stearothermophilus ATCC® 7953™*	R4601218
Geotrichum candidum ATCC® 34614™*	R4603001
Geotrichum capitatum ATCC® 10663™*	R4606800
Geotrichum capitatum ATCC® 28576™*	R4606805
H	
Haemophilus aphrophilus ATCC® 19415™*	R4603002
Haemophilus haemoglobinophilus ATCC® 19416™*	R4603795
Haemophilus haemolyticus ATCC® 33390™*	R4603790
Haemophilus influenzae ATCC® 33930™*	R4603903
Haemophilus influenzae ATCC® 35056™*	R4603840
Haemophilus influenzae ATCC® 35540™*	R4609253
Haemophilus influenzae ATCC® 49144™*	R4603801
Haemophilus influenzae ATCC® 49247™*	R4603830
Haemophilus influenzae ATCC® 49766™*	R4603806

Organism	Order Code (5 loops)
I	
Issatchenkia orientalis ATCC® 6258™*	R4601520
K	
Klebsiella pneumoniae C6	R4603935
Klebsiella oxytoca ATCC® 13182™*	R4603012
Klebsiella oxytoca ATCC® 43086™*	R4603006
Klebsiella oxytoca ATCC® 49131™*	R4607019
Klebsiella oxytoca ATCC® 700324™*	R4603912
Klebsiella oxytoca ATCC® 8724™*	R4603009
Klebsiella pneumoniae subsp. pneumoniae ATCC® 10031™*	R4603930
Klebsiella pneumoniae subsp. pneumoniae ATCC® 13882™*	R4607022
Klebsiella pneumoniae subsp. pneumoniae ATCC® 13883™*	R4607037
Klebsiella pneumoniae subsp. pneumoniae ATCC® 27736™*	R4607031
Klebsiella pneumoniae subsp. pneumoniae ATCC® 33495™*	R4607021
Klebsiella pneumoniae subsp. pneumoniae ATCC® 35657™*	R4603922
Klebsiella pneumoniae subsp. pneumoniae ATCC® 700603™*	R4603074
Klebsiella pneumoniae subsp. pneumoniae ATCC® 9997™*	R4603915
Klebsiella pneumoniae subsp. pneumoniae ATCC® BAA 1705™*	R4609384
Klebsiella pneumoniae subsp. pneumoniae ATCC® BAA 1706™*	R4609385
Kloecera apiculata var. apis ATCC® 32857™*	R4600397
Kloecera japonica ATCC® 58370™*	R4609398
Kocuria kristinae ATCC® BAA-752™*	R4603441
Kocuria rhizophila ATCCv 533™*	R4604070
Kocuria rhizophila ATCC® 9341™*	R4604075
Kocuria rosea ATCC® 186™*	R4603440
L	
Lactobacillus acidophilus ATCC® 314™*	R4603050
Lactobacillus acidophilus ATCC® 4356™*	R4603007
Lactobacillus brevis ATCC® 8287™*	R4603053
Lactobacillus casei ATCC® 393™*	R4603008
Lactobacillus delbrueckii subsp. lactis ATCC® 12315™*	R4603100
Lactobacillus gasserii ATCC® 19992™*	R4603051
Lactobacillus paracasei subsp. paracasei ATCC® BAA-52™*	R4609378
Lactobacillus plantarum ATCC® 8014™*	R4603052
Leclercia adecarboxylata ATCC® 23216™*	R4603926
Leclercia adecarboxylata ATCC® 700325™*	R4603927
Legionella pneumophila ATCC 33823™*	R4603949
Legionella pneumophila ATCC® 33152™*	R4603950
Leuconostoc mesenteroides ATCC® 8293™*	R4603953
Listeria grayi ATCC® 25401™*	R4603959
Listeria innocua ATCC® 33090™*	R4609005
Listeria innocua VC 32293	R4603960
Listeria monocytogenes ATCC® 7644™*	R4603970
Listeria monocytogenes ATCC® 7646™*	R4609316
Listeria monocytogenes ATCC® BAA-751™*	R4603971
Lysinibacillus sphaericus ATCC® 4525™*	R4601229
M	
Malassezia furfur ATCC® 14521™*	R4605025
Malassezia furfur ST 8036	R4604023
Microbacterium liquefaciens ATCC® BAA-1819™*	R4609417
Microbacterium paroxydans ATCC® BAA-1818™*	R4609416
Microbacterium testaceum ATCC® 15829™*	R4604040
Micrococcus luteus ATCC® 10240™*	R4604071
Micrococcus luteus ATCC® 4698™*	R4604085
Micrococcus luteus ATCC® 49732™*	R4604080
Micrococcus lylae ATCC® 27566™*	R4604010
Micrococcus sp. ATCC® 700405™*	R4604086

Quality Control Organisms

Quality Control Organism, qualitative

Organism	Order Code (5 loops)
Microsporium canis ATCC® 11621™*	R4604100
Moraxella catarrhalis ATCC® 23246™*	R4601236
Moraxella catarrhalis ATCC® 25238™*	R4601230
Moraxella catarrhalis ATCC® 25240™*	R4601231
Moraxella catarrhalis ATCC® 49143™*	R4601227
Moraxella catarrhalis ATCC® 8176™*	R4601228
Moraxella osloensis ATCC® 10973™*	R4604020
Moraxella osloensis ATCC® 19976™*	R4604028
Morganella morganii subsp. morganii ATCC® 25830™*	R4604029
Myroides odoratus ATCC® 4651™*	R4602003

N	
Neisseria gonorrhoeae ATCC® 19424™*	R4607041
Neisseria gonorrhoeae ATCC® 31426™*	R4607042
Neisseria gonorrhoeae ATCC® 35541™*	R4609235
Neisseria gonorrhoeae ATCC® 43069™*	R4607043
Neisseria gonorrhoeae ATCC® 43070™*	R4607044
Neisseria gonorrhoeae ATCC® 49226™*	R4609006
Neisseria gonorrhoeae ATCC® 49926™*	R4609019
Neisseria gonorrhoeae ATCC® 49981™*	R4607075
Neisseria lactamica ATCC® 23970™*	R4607045
Neisseria lactamica ATCC® 23971™*	R4607145
Neisseria lactamica ATCC® 49142™*	R4604027
Neisseria meningitidis serogroup C ATCC® 13102™*	R4609007
Neisseria meningitidis serogroup A ATCC® 13077™*	R4607048
Neisseria meningitidis serogroup B ATCC® 13090™*	R4607049
Neisseria meningitidis serogroup Y ATCC® 35561™*	R4607155
Neisseria mucosa ATCC® 19695™*	R4607150
Neisseria perflava ATCC® 14799™*	R4607046
Neisseria sicca ATCC® 29256™*	R4609217
Neisseria sicca ATCC® 9913™*	R4607047
Nocardia asteroides ATCC® 19247™*	R4604111
Nocardia asteroides CL 11014™*	R4609010
Nocardia brasiliensis ATCC® 19296™*	R4609011
Nocardia brasiliensis ATCC® 19297™*	R4604110
Nocardia farcinica ATCC® 3308™*	R4607156

O	
Ochrobactrum anthropi ATCC® 49187™*	R4604150
Ochrobactrum anthropi ATCC® 49687™*	R4604151
Ochrobactrum anthropi ATCC® BAA-749™*	R4604152
Oligella ureolytica ATCC® 43534™*	R4604175
Oligella urethralis ATCC® 17960™*	R4604022

P	
Paenibacillus gordonae ATCC® 29948™*	R4609381
Paenibacillus macerans ATCC® 8509™*	R4601237
Paenibacillus polymyxa ATCC® 43865™*	R4601234
Paenibacillus polymyxa ATCC® 7070™*	R4604600
Paenibacillus polymyxa ATCC® 842™*	R4601235
Parabacteroides distasonis ATCC® 8503™*	R4601226
Parabacteroides distasonis ATCC® BAA-1295™*	R4609377
Parvimonas micra ATCC® 33270™*	R4607203
Pasteurella aerogenes ATCC® 27883™*	R4607103
Pasteurella multocida subsp. multocida ATCC® 43137™*	R4607201
Pediococcus pentosaceus ATCC® 33314™*	R4609254
Penicillium chrysogenum ATCC® 10106™*	R4607202
Peptoniphilus asaccharolyticus ATCC® 29743™*	R4605001
Peptostreptococcus anaerobius ATCC® 27337™*	R4607051
Phialophora verrucosa ATCC® 28181™*	R4607204
Plesiomonas shigelloides ATCC® 14029™*	R4605035
Plesiomonas shigelloides ATCC® 51903™*	R4607218
Porphyromonas gingivalis ATCC® 33277™*	R4609008
Porphyromonas levii ATCC® 29147™*	R4601263
Prevotella melaninogenica ATCC® 25845™*	R4601255
Propionibacterium acidipropionici ATCC® 25562™*	R4607056
Propionibacterium acnes ATCC® 11827™*	R4607052
Propionibacterium acnes ATCC® 6919™*	R4607101
Proteus hauseri ATCC® 13315™*	R4607057
Proteus mirabilis ATCC® 12453™*	R4607059
Proteus mirabilis ATCC® 25933™*	R4607061
Proteus mirabilis ATCC® 29245™*	R4607054
Proteus mirabilis ATCC® 29906™*	R4605055

Organism	Order Code (5 loops)
Proteus mirabilis ATCC® 35659™*	R4607077
Proteus mirabilis ATCC® 43071™*	R4607073
Proteus mirabilis ATCC® 7002™*	R4607053
Proteus vulgaris ATCC® 49132™*	R4605151
Proteus vulgaris ATCC® 6380™*	R4607069
Proteus vulgaris ATCC® 8427™*	R4607058
Prototheca wickerhamii ATCC® 16529™*	R4608014
Providencia alcalifaciens ATCC® 51902™*	R4605169
Providencia stuartii ATCC® 33672™*	R4605170
Providencia stuartii ATCC® 49809™*	R4605171
Pseudomonas aeruginosa ATCC® 10145™*	R4607065
Pseudomonas aeruginosa ATCC® 15442™*	R4607210
Pseudomonas aeruginosa ATCC® 17934™*	R4609208
Pseudomonas aeruginosa ATCC® 27853™*	R4607060
Pseudomonas aeruginosa ATCC® 35032™*	R4607062
Pseudomonas aeruginosa ATCC® 35422™*	R4605200
Pseudomonas aeruginosa ATCC® 9027™*	R4605210
Pseudomonas aeruginosa ATCC® 9721™*	R4607068
Pseudomonas aeruginosa ATCC® BAA-1744™*	R4609392
Pseudomonas fluorescens ATCC® 13525™*	R4607216
Pseudomonas putida ATCC® 49128™*	R4607074
Pseudomonas stutzeri ATCC® 17588™*	R4607211

R	
Ralstonia pickettii ATCC® 49129™*	R4605209
Rhizopus stolonifer ATCC® 14037™*	R4607212
Rhodococcus equi ATCC® 6939™*	R4605400
Rhodotorula glutinis ATCC® 32765™*	R4607095

S	
Saccharomyces cerevisiae ATCC® 4098™*	R4606900
Saccharomyces cerevisiae ATCC® 9763™*	R4608201
Salmonella enterica serovar Typhimurium ATCC® 13311™*	R4606100
Salmonella enterica subsp. enterica serovar Anatum ATCC® 9270™*	R4606006
Salmonella enterica subsp. enterica serovar Choleraesuis ATCC® 10708™*	R4609001
Salmonella enterica subsp. enterica serovar Enteritidis ATCC® 13076™*	R4608200
Salmonella enterica subsp. enterica serovar Montevideo ATCC® 8387™*	R4606004
Salmonella enterica subsp. enterica serovar Newport ATCC® 6962™*	R4606005
Salmonella enterica subsp. enterica serovar Paratyphi A ATCC® 11511™*	R4606002
Salmonella enterica subsp. enterica serovar Paratyphi B ATCC® 8759™*	R4606003
Salmonella enterica subsp. enterica serovar Typhimurium ATCC® 14028™*	R4606000
Salmonella enterica sv Poona NCTC 4840	R4606009
Salmonella sp. not Typhi group D BF-SD	R4608202
Salmonella sp. serovar Abony NCTC 6017	R4606007
Salmonella tranorova NCTC 10252	R4606010
Scopulariopsis acremonium ATCC® 58636™*	R4608204
Serratia liquefaciens ATCC® 27592™*	R4608000
Serratia marcescens ATCC® 13880™*	R4608002
Serratia marcescens ATCC® 14756™*	R4609004
Serratia marcescens ATCC® 8100™*	R4608001
Serratia odorifera ATCC® 33077™*	R4602000
Shewanella putrefaciens ATCC® 49138™*	R4601004
Shewanella putrefaciens ATCC® 8071™*	R4607071
Shigella sonnei group D ATCC® 11060™*	R4609255
Shigella sonnei group D ATCC® 25931™*	R4608150
Shigella sonnei group D ATCC® 9290™*	R4608151
Shigella boydii serovar 1 group C ATCC® 9207™*	R4608110
Shigella flexneri serovar 2b group B ATCC® 12022™*	R4608101
Sphingobacterium multivorum ATCC® 35656™*	R4602004
Sphingobacterium spiritivorum ATCC® 33861™*	R4602001
Sporidiobolus salmonicolor ATCC® MYA-4550™*	R4609396
Sporothrix schenckii ATCC® 10212™*	R4608205
Staphylococcus aureus ATCC® 29737™*	R4607026
Staphylococcus aureus ATCC® 35548™*	R4609446
Staphylococcus aureus ATCC® 9144™*	R4609321
Staphylococcus aureus ATCC® BAA-1026™*	R4601026
Staphylococcus aureus ATCC® BAA-1708™*	R4609389
Staphylococcus aureus subsp. aureus ATCC® 12600™*	R4606507
Staphylococcus aureus subsp. aureus ATCC® 25178™*	R4607015
Staphylococcus aureus subsp. aureus ATCC® 25904™*	R4607007
Staphylococcus aureus subsp. aureus ATCC® 33592™*	R4607020
Staphylococcus aureus subsp. aureus ATCC® 49444™*	R4609023
Staphylococcus aureus subsp. aureus ATCC® 49476™*	R4606510
Staphylococcus aureus subsp. aureus ATCC® 51153™*	R4606508

Organism	Order Code (5 loops)
Staphylococcus aureus subsp. aureus ATCC® 6538P™*	R4609002
Staphylococcus aureus subsp. aureus ATCC® 6538™*	R4607016
Staphylococcus aureus subsp. aureus ATCC® BAA-976™*	R4606512
Staphylococcus aureus subsp. aureus ATCC® BAA-977™*	R4606513
Staphylococcus aureus subsp. aureus ATCC® 25923™*	R4607010
Staphylococcus aureus subsp. aureus ATCC® 29213™*	R4607011
Staphylococcus aureus subsp. aureus ATCC® 29247™*	R4607009
Staphylococcus aureus subsp. aureus ATCC® 33591™*	R4607003
Staphylococcus aureus subsp. aureus ATCC® 33862™*	R4607008
Staphylococcus aureus subsp. aureus ATCC® 43300™*	R4609022
Staphylococcus capitis subsp. capitis ATCC® 35661™*	R4608207
Staphylococcus epidermidis ATCC® 12228™*	R4606500
Staphylococcus epidermidis ATCC® 14990™*	R4606501
Staphylococcus epidermidis ATCC® 29887™*	R4608209
Staphylococcus epidermidis ATCC® 49134™*	R4607017
Staphylococcus epidermidis ATCC® 49461™*	R4608224
Staphylococcus epidermidis ATCC® 700296™*	R4609457
Staphylococcus haemolyticus ATCC® 29970™*	R4608210
Staphylococcus hominis ATCC® 27844™*	R4608211
Staphylococcus lentus ATCC® 700403™*	R4606511
Staphylococcus lugdunensis ATCC® 700328™*	R4608229
Staphylococcus saprophyticus ATCC® 15305™*	R4607014
Staphylococcus saprophyticus ATCC® 35552™*	R4608213
Staphylococcus saprophyticus ATCC® 43867™*	R4608223
Staphylococcus saprophyticus ATCC® 49453™*	R4608214
Staphylococcus saprophyticus ATCC® 49907™*	R4609020
Staphylococcus saprophyticus ATCC® BAA-750™*	R4607029
Staphylococcus sciuri subsp. sciuri ATCC® 29060™*	R4607012
Staphylococcus sciuri subsp. sciuri ATCC® 29061™*	R4608230
Staphylococcus simulans ATCC® 27851™*	R4607013
Staphylococcus xylosus ATCC® 29967™*	R4606504
Staphylococcus xylosus ATCC® 29971™*	R4606502
Staphylococcus xylosus ATCC® 35663™*	R4606505
Staphylococcus xylosus ATCC® 49148™*	R4607018
Staphylococcus xylosus ATCC® 700404™*	R4606506
Stenotrophomonas maltophilia ATCC® 13637™*	R4607067
Stenotrophomonas maltophilia ATCC® 17666™*	R4607072
Stenotrophomonas maltophilia ATCC® 49130™*	R4609009
Stenotrophomonas maltophilia ATCC® 51331™*	R4607064
Streptococcus agalactiae group B ATCC® 12386™*	R4607027
Streptococcus agalactiae group B ATCC® 13813™*	R4608250
Streptococcus agalactiae group B CL 810	R4607005
Streptococcus bovis ATCC® 33317™*	R4609454
Streptococcus criceti ATCC® 19642™*	R4606902
Streptococcus dysgalactiae subsp. equisimilis group G ATCC® 12394™*	R4608219
Streptococcus dysgalactiae subsp. equisimilis ATCC® 35666™*	R4608217
Streptococcus dysgalactiae subsp. equisimilis ATCC® 9542™*	R4608011
Streptococcus dysgalactiae subsp. equisimilis group C ATCC® 12388™*	R4607035
Streptococcus equi subsp. equi ATCC® 9528™*	R4601001
Streptococcus equi subsp. zooepidemicus ATCC® 43079™*	R4608237
Streptococcus equi subsp. zooepidemicus group C ATCC® 700400™*	R4608225
Streptococcus gallolyticus ATCC® 49147™*	R4607078
Streptococcus gallolyticus ATCC® 9809™*	R4607025
Streptococcus gallolyticus subsp. gallolyticus ATCC® 49475™*	R4608175
Streptococcus mutans ATCC® 25175™*	R4607001
Streptococcus mutans ATCC® 35668™*	R4607002
Streptococcus oralis ATCC® 35037™*	R4609461
Streptococcus oralis ATCC® 9811™*	R4608220
Streptococcus pasteurianus ATCC® 49133™*	R4607079
Streptococcus pneumoniae ATCC® 27336™*	R4607033
Streptococcus pneumoniae ATCC® 49136™*	R4607004
Streptococcus pneumoniae ATCC® 49150™*	R4606901
Streptococcus pneumoniae ATCC® 49619™*	R4609015
Streptococcus pneumoniae ATCC® 6301™*	R4606909
Streptococcus pneumoniae ATCC® 6303™*	R4608221
Streptococcus pneumoniae ATCC® 6305™*	R4607024
Streptococcus pneumoniae CL 811	R4606700

Organism	Order Code (5 loops)
Streptococcus pyogenes ATCC® 19615™*	R4607000
Streptococcus pyogenes group A ATCC® 12384™*	R4609242
Streptococcus pyogenes group A ATCC® 21547™*	R4609211
Streptococcus salivarius serotype II ATCC® 13419™*	R4608222
Streptococcus sanguinis Type 1 ATCC® 10556™*	R4607023
Streptococcus sp Gp D ATCC® 9854™*	R4606853
Streptococcus sp. group B ATCC® 12401™*	R4609453
Streptococcus sp. group D ATCC® 27284™*	R4607034
Streptococcus sp. Type 2 group F ATCC® 12392™*	R4608218
Streptococcus thermophilus ATCC® 19258™*	R4608238
Streptococcus uberis ATCC® 700407™*	R4608236
Streptococcus uberis ATCC® 9927™*	R4607036
Streptomyces albus ATCC® 17900™*	R4609024
Streptomyces griseus subsp. griseus ATCC® 10137™*	R4609013

T

Tatlockia micdadei ATCC® 33204™*	R4609205
Trichophyton equinum ATCC® 12544™*	R4608302
Trichophyton mentagrophytes ATCC® 9533™*	R4608300
Trichophyton rubrum ATCC® 28188™*	R4608301
Trichophyton tonsurans ATCC® 28942™*	R4609014
Trichophyton verrucosum ATCC® 42898™*	R4608305
Trichosporon cutaneum ATCC® 28592™*	R4608303
Trichosporon mucoides ATCC® 204094™*	R4608015

V

Veillonella parvula ATCC® 10790™*	R4609200
Veillonella parvula ATCC® 17745™*	R4609462
Vibrio parahaemolyticus ATCC® 17802™*	R4609000
Vibrio vulnificus ATCC® 27562™*	R4609017
Virgibacillus pantothenticus ATCC® 14576™*	R4609025

Y

Yarrowia lipolytica ATCC® 9773™*	R4609003
Yersinia enterocolitica ATCC® 23715™*	R4607076
Yersinia enterocolitica subsp. enterocolitica ATCC® 9610™*	R4609018
Yersinia kristensenii ATCC® 33639™*	R4609021

Z

Zygosaccharomyces bailii ATCC® MYA-4549™*	R4609397
---	----------



Quality Control Organisms

Quality Control Organism Sets, qualitative

Quality Control Organism Sets, qualitative

Thermo Scientific™ RapID™ QC sets and Thermo Scientific™ MicroScan™ sets contain five loops of each of the listed organisms. Sensititre sets contain one loop of each of the listed organisms.

Set	Order Code
MicroScan Gram-Negative QC Set/4	
Enterococcus faecalis ATCC® 29212™* Klebsiella oxytoca ATCC® 49131™* Proteus vulgaris ATCC® 49132™* Pseudomonas aeruginosa ATCC® 22853™*	R4661000
MicroScan Gram-Positive QC Set/4	
Enterococcus faecalis ATCC® 29212™* Micrococcus luteus ATCC® 49732™* Staphylococcus aureus subsp. aureus ATCC® 29213™* Streptococcus gallolyticus ATCC® 49147™*	R4661001
RapID ANA II QC Set/3	
Bacteroides uniformis ATCC® 8492™* Clostridium sordellii ATCC® 9714™* Parabacteroides distasonis ATCC® 8503™*	R4653050
RapID CB PLUS QC Set/3	
Arcanobacterium pyogenes ATCC® 19411™* Corynebacterium pseudodiphtheriticum ATCC® 10701™* Paenibacillus polymyxa ATCC® 842™*	R4653048
RapID NH QC Set/4	
Haemophilus influenzae Type a ATCC® 9006™* Haemophilus paraphrophilus ATCC® 49146™* Moraxella catarrhalis ATCC® 8176™* Oligella urethralis ATCC® 17960™*	R4653051
RapID NF PLUS QC Set/4	
Acinetobacter baumannii ATCC® 19606™* Aeromonas hydrophila ATCC® 35654™* Elizabethkingia meningoseptica ATCC® 13253™* Oligella ureolytica ATCC® 43534™*	R4653054
RapID ONE QC Set/4	
Enterobacter aerogenes ATCC® 13048™* Escherichia Coli ATCC® 25922™* Proteus vulgaris ATCC® 6380™* Pseudomonas aeruginosa ATCC® 27853™*	R4653056
RapID SS/u QC Set/5	
Enterococcus faecalis ATCC® 29212™* Escherichia coli ATCC® 25922™* Klebsiella pneumoniae subsp. pneumoniae ATCC® 13883™* Proteus mirabilis ATCC® 25933™* Serratia marcescens ATCC® 8100™*	R4653053
RapID STAPH PLUS QC Set/4	
Staphylococcus haemolyticus ATCC® 29970™* Staphylococcus saprophyticus ATCC® 35552™* Enterobacter aerogenes ATCC® 13048™* Oligella ureolytica ATCC® 43534™*	R4653047
RapID STR QC Set/4	
Enterococcus durans ATCC® 11576™* Enterococcus aecalis ATCC® 29212™* Streptococcus gallolyticus ATCC® 9809™* Streptococcus pyogenes ATCC® 19615™*	R4653052
RapID YEAST PLUS QC Set/5	
Candida albicans ATCC® 14053™* Candida glabrata ATCC® 2001™* Candida kefyr ATCC® 2512™* Cryptococcus laurentii ATCC® 66036™* Yarrowia lipolytica ATCC® 9773™*	R4653060
Sensititre Anaerobe MIC Susc QC Set/5	
Bacteroides fragilis ATCC® 25285™* Bacteroides thetaiotaomicron ATCC® 29741™* Eggerthella lenta ATCC® 43055™*	R4687021
Sensititre AQUACULTURE MIC QC Set/2	
Aeromonas salmonicida subsp. salmonicida ATCC® 33658™* Escherichia coli ATCC® 25922™*	R4687022

Set	Order Code
Sensititre MIC Breakpoint QC Set/5	
Enterococcus faecalis ATCC® 29212™* Escherichia coli ATCC® 35218™* Escherichia coli ATCC® 25922™* Pseudomonas aeruginosa ATCC® 27853™* Staphylococcus aureus subsp. aureus ATCC® 29213™*	R4687020
Sensititre Automated Gram Neg QC Set/6	
Edwardsiella tarda ATCC® 15947™* Klebsiella oxytoca ATCC® 8724™* Morganella morganii subsp. morganii ATCC® 25830™* Proteus vulgaris ATCC® 6896™* Pseudomonas aeruginosa ATCC® 10145™* Shigella sonnei ATCC® 25931™*	R4687023
Sensititre GPID ID QC Set/5	
Enterococcus faecalis ATCC® 29212™* Escherichia coli ATCC® 25922™* Kocuria rosea ATCC® 186™* Pseudomonas aeruginosa ATCC® 27853™* Staphylococcus epidermidis ATCC® 700296™*	R4687024
Sensititre H.infl/S.pneu MIC QC Set/3	
Haemophilus influenzae ATCC® 49247™* Haemophilus influenzae ATCC® 49766™* Streptococcus pneumoniae ATCC® 49619™*	R4687025
Sensititre MIC Susceptibility QC Set/9	
Enterococcus faecalis ATCC® 29212™* Enterococcus faecalis ATCC® 51299™* Escherichia coli ATCC® 35218™* Escherichia coli ATCC® 25922™* Klebsiella pneumoniae subsp. pneumoniae ATCC® 700603™* Pseudomonas aeruginosa ATCC® 27853™* Staphylococcus aureus subsp. aureus ATCC® 29213™* Staphylococcus aureus subsp. aureus ATCC® BAA-976™* Staphylococcus aureus subsp. aureus ATCC® BAA-977™*	R4687026
Sensititre YeastOne QC Set/2	
Candida parapsilosis ATCC® 22019™* Issatchenkia orientalis ATCC® 6258™*	R4687027
Vitek GNI+ QC Set/8	
Bordetella bronchiseptica ATCC® 10580™* Burkholderia cepacia ATCC® 25608™* Klebsiella pneumoniae subsp. pneumoniae ATCC® 13883™* Leclercia adecarboxylata ATCC® 23216™* Plesiomonas shigelloides ATCC® 51903™* Proteus mirabilis ATCC® 7002™* Providencia alcalifaciens ATCC® 51902™* Serratia liquefaciens ATCC® 27592™*	R4652014
Vitek* GPI QC Set/7	
Enterococcus durans ATCC® 6056™* Enterococcus faecalis ATCC® 29212™* Erysipelothrix rhusiopathiae ATCC® 19414™* Staphylococcus xylosum ATCC® 29971™* Streptococcus equi subsp. equi ATCC® 9528™* Streptococcus gallolyticus ATCC® 9809™* Streptococcus pyogenes ATCC® 19615™*	R4652010
Vitek GPS Susceptibility QC Set/4	
Enterococcus faecalis ATCC® 29212™* Enterococcus faecalis ATCC® 51299™* Escherichia coli ATCC® 35218™* Staphylococcus aureus subsp. aureus ATCC® 29213™*	R4652008
Vitek NHI QC Set/6	
Aggregatibacter aphrophilus ATCC® 7901™* Haemophilus influenzae Type a ATCC® 9006™* Moraxella catarrhalis ATCC® 25238™* Neisseria gonorrhoeae ATCC® 49981™* Neisseria lactamica ATCC® 23970™* Neisseria meningitidis serogroup C ATCC® 13102™*	R4652006

Set	Order Code
-----	------------

Vitek Urine ID Card (UID and UID-3) QC Set/5

Enterococcus faecalis ATCC® 29212™*	
Escherichia coli ATCC® 25922™*	
Klebsiella pneumoniae subsp. pneumoniae ATCC® 27736™*	R4652013
Proteus mirabilis ATCC® 7002™*	
Pseudomonas aeruginosa ATCC® 27853™*	

Vitek 2 AST-GN QC Set/4

Escherichia coli ATCC® 25922™*	
Escherichia coli ATCC® 35218™*	
Klebsiella pneumoniae subsp. pneumoniae ATCC® 700603™*	R4652018
Pseudomonas aeruginosa ATCC® 27853™*	

Vitek 2 AST-GP QC Set/7

Enterococcus faecalis ATCC® 29212™*	
Enterococcus faecalis ATCC® 51299™*	
Staphylococcus aureus subsp. aureus ATCC® 29213™*	
Staphylococcus aureus subsp. aureus ATCC® BAA-976™*	R4652019
Staphylococcus aureus subsp. aureus ATCC® BAA-977™*	
Staphylococcus aureus subsp. aureus ATCC® BAA-1026™*	
Streptococcus pneumoniae ATCC® 49619™*	

Vitek* 2 BCL 21345 QC Set/12

Aneurinibacillus aneurinilyticus ATCC® 11376™*	
Bacillus badius ATCC® 14574™*	
Bacillus circulans ATCC® 61™*	
Bacillus megaterium ATCC® 14581™*	
Bacillus pumilus ATCC® BAA-1434™*	
Brevibacillus agri ATCC® 51663™*	
Brevibacillus laterosporus ATCC® 64™*	R4652023
Enterobacter aerogenes ATCC® 13048™*	
Paenibacillus gordonae ATCC® 29948™*	
Paenibacillus macerans ATCC® 8509™*	
Paenibacillus polymyxa ATCC® 7070™*	
Staphylococcus epidermidis ATCC® 12228™*	

Vitek 2 CBC Comprehensive QC Set/11

Arcanobacterium haemolyticum ATCC® BAA-1784™*	
Cellulosimicrobium cellulans ATCC® BAA-1816™*	
Cellulosimicrobium cellulans ATCC® BAA-1817™*	
Corynebacterium renale ATCC® BAA-1785™*	
Corynebacterium urealyticum ATCC® 43044™*	
Curtobacterium pusillum ATCC® 19096™*	R4652031
Klebsiella oxytoca ATCC® 700324™*	
Microbacterium liquefaciens ATCC® BAA-1819™*	
Microbacterium paraoxydans ATCC® BAA-1818™*	
Microbacterium testaceum ATCC® 15829™*	
Ochrobactrum anthropi ATCC® BAA-749™*	

Vitek 2 GN 21341 QC Set/10

Acinetobacter baumannii ATCC® BAA-747™*	
Elizabethkingia meningoseptica ATCC® 13253™*	
Enterobacter hormaechei ATCC® 700323™*	
Klebsiella oxytoca ATCC® 700324™*	
Ochrobactrum anthropi ATCC® BAA-749™*	R4652021
Proteus vulgaris ATCC® 6380™*	
Pseudomonas aeruginosa ATCC® 9721™*	
Pseudomonas aeruginosa ATCC® BAA-1744™*	
Shigella sonnei group D ATCC® 25931™*	
Stenotrophomonas maltophilia ATCC® 17666™*	

Vitek 2 GP 21342 QC Set/9

Enterococcus casseliflavus ATCC® 700327™*	
Enterococcus saccharolyticus ATCC® 43076™*	
Kocuria kristinae ATCC® BAA-752™*	
Listeria monocytogenes ATCC® BAA-751™*	
Staphylococcus saprophyticus ATCC® BAA-750™*	R4652020
Staphylococcus sciuri subsp. sciuri ATCC® 29061™*	
Streptococcus equi subsp. zooepidemicus ATCC® 43079™*	
Streptococcus pneumoniae ATCC® 49619™*	
Streptococcus thermophilus ATCC® 19258™*	

Vitek 2 NH 21346 QC Set/9

Aggregatibacter aphrophilus ATCC® 33389™*	
Eikenella corrodens ATCC® BAA-1152™*	
Enterobacter aerogenes ATCC® 13048™*	
Haemophilus influenzae ATCC® 9007™*	
Neisseria gonorrhoeae ATCC® 19424™*	R4652024
Neisseria lactamica ATCC® 23970™*	
Oligella urethralis ATCC® 17960™*	
Paenibacillus polymyxa ATCC® 7070™*	
Staphylococcus epidermidis ATCC® 12228™*	

Vitek 2 YST 21343 QC Set/11

Candida albicans ATCC® 33389™*	
Candida glabrata ATCC® MYA-2950™*	
Candida lusitanae ATCC® 34499™*	
Candida utilis ATCC® 9950™*	
Kloeckera japonica ATCC® 58370™*	
Oligella ureolytica ATCC® 43534™*	R4652022
Prototheca wickerhamii ATCC® 16529™*	
Sporidiobolus salmonicolor ATCC® MYA-4550™*	
Staphylococcus epidermidis ATCC® 12228™*	
Trichosporon mucoides ATCC® 204094™*	
Zygosaccharomyces bailii ATCC® MYA-4549™*	

Quality Control Organisms

Quality Control Organism and Sets, quantitative, Water Quality Control Organism and Sets

Quality Control Organism, quantitative

Each pack contains ten vials. Each kit is designed to deliver <100 CFU/0.1mL. Total diluent volume per vial is 1.1mL. Thermo Scientific™ Quanti-Cult™ Plus should be stored at 2–8°C. Blue foam rack and extra diluents available separately.

Organism	Order Code
Aspergillus brasiliensis ATCC® 16404™*	R4711100
Bacillus subtilis ATCC® 6633™*	R4711221
Brevundimonas diminuta ATCC® 19146™*	R4711300
Burkholderia cepacia ATCC® 25416™*	R4715220
Candida albicans ATCC® 10231™*	R4711503
Candida albicans ATCC® 2091™*	R4711601
Clostridium sporogenes ATCC® 11437™*	R4711703
Clostridium sporogenes ATCC® 19404™*	R4711700
Enterococcus faecalis ATCC® 29212™*	R4717030
Enterococcus faecium ATCC® 35667™*	R4711956
Escherichia coli ATCC® 25922™*	R4717050
Escherichia coli ATCC® 8739™*	R4717085
Kocuria rhizophila ATCC® 9341™*	R4714075
Pseudomonas aeruginosa ATCC® 9027™*	R4715210
Salmonella enterica subsp. enterica serovar Choleraesuis ATCC® 10708™*	R4719001
Salmonella enterica subsp. enterica serovar Typhimurium ATCC® 13311™*	R4716100
Salmonella enterica subsp. enterica serovar Typhimurium ATCC® 14028™*	R4716000
Salmonella sp. serovar Abony NCTC 6017	R4716007
Serratia marcescens ATCC® 8100™*	R4718001
Shigella sonnei ATCC® 25931™*	R4716300
Staphylococcus aureus subsp. aureus ATCC® 29737™*	R4717020
Staphylococcus aureus subsp. aureus ATCCv 6538P™*	R4719002
Staphylococcus aureus subsp. aureus ATCC® 6538™*	R4717016
Staphylococcus epidermidis ATCC® 12228™*	R4716500
Streptococcus pyogenes ATCC® 19615™*	R4717000

Quality Control Organism Sets, quantitative

Quanti-Cult Plus sets contain two vials of each of the listed organisms.

Set	Order Code
Growth Promotion Testing QC Set/5	R4721000
Aspergillus brasiliensis ATCC® 16404™*	
Bacillus subtilis ATCC® 6633™*	
Candida albicans ATCC® 10231™*	
Clostridium sporogenes ATCC® 11437™*	
Pseudomonas aeruginosa ATCC® 9027™*	
Growth Promotion 2 QC Set/5	R4721001
Aspergillus brasiliensis ATCC® 16404™*	
Bacillus subtilis ATCC® 6633™*	
Candida albicans ATCC® 10231™*	
Pseudomonas aeruginosa ATCC® 9027™*	
Staphylococcus aureus subsp. aureus ATCC® 6538™*	
Nutritional Adequacy Testing QC Set/6	R4723000
Aspergillus brasiliensis ATCC® 16404™*	
Bacillus subtilis ATCC® 6633™*	
Candida albicans ATCC® 10231™*	
Escherichia coli ATCC® 8739™*	
Pseudomonas aeruginosa ATCC® 9027™*	
Staphylococcus aureus subsp. aureus ATCC® 6538™*	

Water Quality Control Organism and Sets

Water/Wastewater quality control testing microorganisms in Quanti-Cult format. WaterBugs deliver <50 CFU in a single inoculum to mimic water samples. One test/reconstituted vial.

Set	Order Code
Water QC Kit 1002 SET/3 Contents: 9 WaterBugs (3 per organism); 9 rehydration vials.	R4741002
Enterobacter cloacae subsp. cloacae ATCC® 13047™*	
Escherichia coli Strain C1	
Proteus mirabilis ATCC® 12453™*	
Water QC Kit 1003 SET/2 Contents: 10 WaterBugs (5 per organism); 10 rehydration vials.	R4741003
Escherichia coli Strain C1	
Pseudomonas aeruginosa ATCC® 27853™*	
Water QC Kit 1004 SET/2 Contents: 10 WaterBugs (5 per organism); 10 rehydration vials.	R4741004
Enterobacter cloacae subsp. cloacae ATCC® 13047™*	
Escherichia coli Strain C1	
Water QC Kit 1009 SET/1 Contents: 10 WaterBugs (10 per organism); 10 rehydration vials.	R4741009
Enterococcus faecalis ATCC® 29212™*	

WaterBugs

Water/Wastewater quality control testing microorganisms in Quanti-Cult format. Thermo Scientific™ WaterBugs™ deliver <50 CFU in a single inoculum to mimic water samples. One test/reconstituted vial. Contents: ten preserved organism vials; ten rehydration vials.

Set	Order Code
Enterobacter cloacae subsp. cloacae ATCC® 13047™*	R4757090
Escherichia coli Strain C1	R4751985
Klebsiella pneumoniae Strain C6	R4753935
Proteus mirabilis ATCC® 12453™*	R4757059
Pseudomonas aeruginosa ATCC® 27853™*	R4757060



Look for the ATCC Licensed Derivative (Registered Mark) Emblem for products derived from ATCC cultures.

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalog marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

thermoscientific.com/microbiology

Agilent is a registered trademark of Agilent Technologies

© 2017 Thermo Fisher Scientific Inc. All rights reserved. AOAC is a trademark of AOAC International. BagClip, BagFilter, BagLight, BagOpen, BagPage, BagRack and RollBag are all trademarks of INTERSCIENCE. BAX and RiboPrinter are trademarks of Hygiene LLC. ELISA Systems is a trademark of ELISA Systems Pty Ltd. ELITE MGB is a trademark of Elitech Holding B.V. Agilent is a registered trademark of Agilent Technologies. Fungizone is a trademark of E. R. Squibb and Sons. Imegen is a trademark of Instituto de Medicina Genomica. Isolator is a trademark of Alere, New York. Kleenpak is a trademark of the Pall Corporation. Lynx is a trademark of Merck KGaA. Microsoft and Windows are registered trademarks or trademarks of Microsoft Corporation. NF Validation is a trademark of Association Française de Normalisation (AFNOR). Opta is a trademark of Sartorius Stedim Biotech. ReadyMate is a trademark of GE Healthcare Companies. Steam-Thru is a trademark of Colder Products Company. TaqMan® is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. Tween is a trademark of the Croda group of companies. VERIGENE is a trademark of NANOSPHERE, Inc.

Specifications and terms are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.